




RESTORATIVE PROSTHETIC PRODUCT BOOK

13TH EDITION



LABORATORY & CLINICAL SOLUTIONS

IMPLANT PARTS | ATTACHMENTS | DIGITAL DENTISTRY

Two generations and over 40 years of service.

Founded in 1980, Preat Corporation has grown from a family-owned business into a trusted leader in dental attachments, implant components and prosthetic solutions.

Known for innovation and expert support, Preat pioneered attachment systems that enhanced function, esthetics and patient comfort. Over the years, the company expanded into custom abutments and implant parts, partnering with leading manufacturers. In 2020, Preat joined Young Innovations, further strengthening its ability to deliver innovative solutions and outstanding customer service.



**Scan to learn
more about
Preat's history!**

Thomas Bormes, Founder, with
his son, Chris Bormes, who served
as President from 2011-2020.



The Preat Promise



Clinical Choice

With a full portfolio of restorative options for single-tooth to full arch cases, we give you the freedom to choose the best solution for each patient.



Expert Support

Our restorative experts are here to help you identify implants, review your restorative options and even answer questions years after the restoration was placed – that's what we call lifecycle support and it's what we do every day.



Easy Ordering

We'll help you find the right item to fit your restorative plan. Our broad portfolio also means we can be your one-stop shop for implant restorations to make ordering even easier!



Lab Trusted

Over the decades, we've earned the trust of dental laboratories. We're here to be a true part of your implant restorative team.



Quality Tested

We know quality is essential. Preat products are tested to the highest standards and meet or exceed the fit of OEM providers.



Timely Delivery

99.4% on-time and accurate delivery means you and your patient receive exactly what you need when you need it. That's why we do what we do.

TABLE OF CONTENTS

ATTACHMENTS

OVERDENTURE BARS

Hader Bar	6
Ackermann MD & BL	11
Universal Plunger Loc	12
LOCATOR® Bar	34
Clix	26
Sagix	55
OT Equator®	52
Ceka and Ceka Revax	59

OVERDENTURE STUDS

PreatLoc™	14
O-Ring Size Chart	21
O-Rings / OSO	22
Sphero Reconstruction Kits	13
LOCATOR®	31
Ceka and Ceka Revax	59

IMPLANT ABUTMENTS / ATTACHMENTS/ PREATLOC™

O-Ring Abutment	24
Clix	26
LOCATOR® Abutment	37
LOCATOR R-Tx®	45
LOCATOR® Overdenture Implant System	50
OT Equator®	52
Sphero Flex / Sphero Block	53

EXTRACORONAL ATTACHMENTS *(for removable partial dentures)*

Sphero Retrofit Kit	55
Sagix	55
Vertex	57
American Mini (Mini Dalbo)	57
O-Ring Distal Extension Attachment	58
Mays® Unilateral	58
Universal IC Attachment	59
Ceka and Ceka Revax	59

CROWN & BRIDGE SPECIALTY

Splintastik	61
Double Mandrel BPD	61

FIBER TECHNOLOGY

TriLor® High Performance Polymer	63
--	----

GRAPHENE TECHNOLOGY

G-CAM	65
-----------------------------	----

CLINICAL SPECIALTY

Attachment Pick Up Material	68
Chairside® Block Out	68
Denture Removal Tool	68
Perma Ret	69
Pro Blue NTR Polish	69

CAD-CAM BAR SOLUTIONS

Attachment Options	70
--	----

IMPLANT PROSTHETICS

Quick Reference Chart	74
3i®	76
ACE Surgical®	83
Ankylos®	84
Astra®	86
Astra Tech EV	91
BioHorizons®	101
Bio-Lok®	100
Blue Sky Bio®	100
Branemark®	102
Camlog®	104
Glidewell Laboratories®	110
Hiossen®	111
Hi-Tec®	121
Implant Direct®	116
Keystone®	118
Lifecore®	121
MegaGen®	123
MIS®	125
Multi-Unit Components	126
Neodent®	130
Neoss®	132
Nobel Biocare®	132
Southern®	146
Sterngold®	146
Straumann®	147
Zimmer®	157

IMPLANT DRIVERS AND TOOLS

Implant Buddy	163
Implant Drivers	164
Cuff Height and Angulation Tool Set	166


RESOURCES

Alloy Chart	71
Implant Compatibility Chart	166
Implant Driver Size	163
Implant Prosthetics Quick Reference Chart	74
O-Ring Replacement Size Chart	21
Torque Value Chart	174



HADER BARS

The Hader Bar is the most popular and versatile bar attachment in the world as it is low profile, provides secure retention options, and is simple to service!



BAR PATTERNS

- » Available in standard and short heights
- » Plastic-wax is ideal for fast set investment



ALUMINUM ANALOG BAR

- » For relines, repairs, and making new prostheses



HADER BAR MANDREL

- » To assure parallelism between multiple bar segments



PROCESSING SPACER


- » The processing spacer will prevent the over extension of acrylic which inhibits flexation of the female ride



FEMALE RIDERS

Retentions of female riders


- » **Yellow**-standard
- » **White**-decreased
- » **Red**-increased
- » **Blue**-for worn bars



COMBO KIT

Everything needed for one bar case


- » 1 Plastic Bar Pattern
- » 2 Female Riders
- » 2 Metal Housings (either style)
- » 2 Processing Spacers



HADER BAR SERVICING KIT


Ideal for relines, rebases, repairs, or fabricating a new overdenture on existing bars or CAD-CAM bars

- » 1 Delrin Analog Bar
- » 2 Female Riders
- » 2 Metal Housings (either style)
- » 2 Processing Spacers



HOUSINGS



Easy and accurate alignment of riders and housings with space saving options.



METAL CLIPS

Ackermann gold and stainless and C118 gold and stainless

Why Hader Bars are the Best Option for Bar Overdentures



INDICATIONS

- » In the maxillae, the bone quality and quantity often dictates the use of a bar to protect the implants in weaker bone via splinting.
- » In the resorbed mandible, the atrophied ridge can not be engaged by the denture flange and therefore does not provide prosthesis stability.
- » Extremely divergent implants are easier to accommodate with a bar overdenture.

ADVANTAGES OF THE HADER BAR

Flexibility in Prostheses Function



- » Metal Ackermann clips provide both vertical and rotational movement, ideal for 2 abutments bar with a single straight section and for protecting abutments.
- » Plastic Hader clips provide rotational movement. Ideal for restorations utilizing multiple bar sections and/or for rigid/fully abutment supported prostheses.

Space Saving

- » Total vertical height of bar, 2.5mm.
- » Total vertical of attachments, 1.5mm.
- » Total width of attachment, 4.2mm.

Easy to Service

- » 4 retentions of plastic clips to provide patient specific retention
Replacing a Hader clip can be done in less than 30 seconds
- » 2 designs and 2 alloys of metal Ackermann clips



ATTACHMENTS

OVERDENTURE BAR ATTACHMENTS

DIRECT RETAINERS

Hader Bar Attachment

4mm† x 3mm



Dolder Bar Attachment

4mm† x 3.5mm




4 Retention clips provide rotational movement on the bar

Ackermann Bar Attachment—

5.5mm† x 4.5mm

Vertical movement



Vertical & Rotation motion

ON TOP OF BAR

LOCATOR® Attachment

4mm† x 5.5mm

8 Retention Inserts. Attachment can hinge, have vertical and rotational movement.



Clix Attachment

4mm† x 4mm

Attachment can hinge, have vertical and rotational movement.



OT Equator® Attachment

5.5mm† x 4.5mm

System can hinge, have vertical and rotational movement.



OFFSET ATTACHMENTS

Universal Plunger Loc

2.5mm† x 4mm


Can be rigid or resilient.



Sagix Attachment

3.1mm† x 3mm

Attachment can hinge and rotate.





Vertex Attachment

3.5mm† x 3mm

Allows vertical movement.



Hader Bar: the most popular and versatile bar attachment in the world. Low profile, secure retention options, and simple servicing! 3.5mm tall 4.20mm Ø housing







Hader Bar	
Item #	Product
0101701	Hader Bar Kit 2 castable bars, 6 yellow riders, 6 green processing, 1 insertion tool
0101750	Hader Bar Master Kit 6 castable bar, 6 alignment metal housings, 6 yellow riders, 6 green processing, 1 insertion tool
0101741	Hader Combo Kit 1 castable bar, 2 metal housings, 2 yellow riders, 2 green processing, and 1 insertion tool
0101742	Hader Alignment Combo Kit 1 castable bar, 2 alignment metal housings, 2 yellow riders, 2 green processing, 1 insertion tool
0101745	Hader Servicing Kit Delrin Analog, 2 metal housings, 2 yellow riders, 2 green processing, 1 insertion tool
0101746	Hader Bar Alignment Servicing Kit 1 Delrin Analog, 2 alignment housing, 2 yellow riders, 2 green processing, 1 insertion tool
0101704-6	Hader Plastic Bar Patterns (Standard)(6-Pack)
0101704-50	Hader Plastic Bar Patterns (Standard)(50-Pack)
0101709-6	Hader Plastic Bar Pattern (Short)(6-Pack)

ATTACHMENTS

OVERDENTURE BARS

10

Hader Bar




Hader Bar

Item #	Product
0101709-50	Hader Plastic Bar Pattern (Short) (50-Pack)
0101775	Hader Aluminum Analog
0101705	Hader Insertion Tool
0101703-6	Hader Green Processing Spacers (6-Pack)
0101703-24	Hader Green Processing Spacers (24-Pack)
0101703-50	Hader Green Processing Spacers (50-Pack)
0101720-6	Hader Alignment Housing (6-Pack)
0101720-24	Hader Alignment Housing (24-Pack)
0101720-50	Hader Alignment Housing (50-Pack)
0101710-6	Hader Metal Housings (6-Pack)
0101710-24	Hader Metal Housings (24-Pack)
0101710-50	Hader Metal Housings (50-Pack)
0101753-6	Hader White Riders - Reduced (6-Pack)
0101753-24	Hader White Riders - Reduced (24-Pack)
0101753-50	Hader White Riders - Reduced (50-Pack)
0101752-6	Hader Yellow Riders - Standard (6-Pack)
0101752-24	Hader Yellow Riders - Standard (24-Pack)
0101752-50	Hader Yellow Riders - Standard (50-Pack)



Hader Bar	
Item #	Product
0101754-6	Hader Red Riders - Increased (6-Pack)
0101754-24	Hader Red Riders - Increased (24-Pack)
0101754-50	Hader Red Riders - Increased (50-Pack)
0101760-6	Hader Blue Riders - Worn Bars (6-Pack)
0101760-24	Hader Blue Riders - Worn Bars (24-Pack)
0101760-50	Hader Blue Riders - Worn Bars (50-Pack)
0101755	Hader Combo Pack - (2) white riders, (2) yellow riders, (2) red riders

Ackermann Clips: Metal Riders provide both rotational and vertical movement with a choice of two styles of acrylic retentions



Metal Riders for Hader/Round Bars	
Item #	Product
0102108	Ackermann BL Clip SS
0102108-6	Ackermann BL Clip SS (6-Pack)
0102118	Ackermann MD Clip SS
0102118-6	Ackermann MD Clip SS (6-Pack)
0102119	Ackermann MD Clip GOLD
0102119-6	Ackermann MD Clip GOLD (6-Pack)
0102124	Ackermann BL Clip GOLD
0102124-6	Ackermann BL Clip GOLD (6-Pack)


UNIVERSAL PLUNGER LOC ATTACHMENT

Titanium locking pin attachment designed for resilient or rigid retention on implant overdenture bars.




- » Biocompatible Titanium
- » Simple to Service
- » Standard Techniques

One universal size may be adjusted from 7.5mm to 4mm in length.

The 1.5mm diameter plunger is compatible with the Low Passive and Swiss Loc NG attachment systems.



Plunger Loc: titanium locking plunger for Overdenture Bars. May be rigid or resilient. 7.5mm long, may be cut to 6, 5, or 4mm. 1.5mm Ø plunger Plunger Diameter compatible with Swiss-Loc & Low Passive



Universal Plunger Loc Attachment	
Item #	Product
1776000	Universal Plunger Loc Attachment
1776005	Universal Plunger Loc Processing Jig
1776010	Universal Plunger Loc Male Only
1776015	Universal Plunger Loc Twist Drill
1776030	Universal Plunger Loc Delrin Inserts

Sphero Reconstruction Kits: designed to salvage and restore worn mini-implants, direct placement or post coping balls, Sagital Balls (like Bredent VKS and Sagix) and o-ring attachments



Reconstruction
Kit



Reconstruction
Kit




0800266

Concave Sphere Retrofit Kits

Item #	Product
0800142	Sphero Reconstruction Kit 1.8mm Ø 2.7mm* vertical, 4.3mm Ø, 2 hollow TiN spheres, 1 strip holder, 2 pink females, 1 insertion tool
0800280	1.8mm Stainless Cap Housing (2-Pack)
0800144	Sphero Reconstruction Kit 2.2mm 2.7mm* vertical, 4.3mm Ø, 2 hollow TiN spheres, 1 strip holder, 2 pink females, 1 insertion tool
0800146	Sphero Reconstruction Kit 2.5mm Ø 2.7mm* vertical, 4.8mm Ø, 2 hollow TiN spheres, 1 strip holder, 2 pink females, 1 insertion tool
0800180	2.5mm Stainless Cap Housing (2-Pack)
0800266	1.8mm Retention Caps, White (6-Pack)

13

Concave Sphere Retrofit Kits



Introducing PreatLoc™

The Next Generation in Overdenture Attachments

Backed by 45 years of Preat overdenture expertise and exceptional support, **PreatLoc™** helps you deliver reliable performance and long-lasting smiles for their patients.

**Compatible with the leading
overdenture attachment system!**



The PreatLoc™ Advantage



Compatible With All Major Implant Systems

Compatible with 18 implant systems and 35 implant connections.



Tissue-Friendly ZrN Coating for Healthier Smiles

Less biofilm and less plaque adhesion reduce inflammation.



Trusted, Long-Lasting Retention

In 5,000-cycle internal testing, PreatLoc™ inserts retain reliably over time.




Experience PreatLoc™
Priced up to 33% less than
leading attachment systems!

SYSTEM OVERVIEW

PreatLoc™ is a complete overdenture attachment system designed for reliable performance and seamless integration.


Tissue-Friendly ZrN Coating

- Biocompatible surface with reduced plaque adhesion^{1,2}
- Has shown greater bacterial resistance than TiN²



Reliable Retention

- Standard pivot inserts accommodate up to 20° divergence
- Extended pivot inserts for up to 40° divergence
- Demonstrated reliable and long-lasting performance³




Anodized Aesthetic Denture Cap

- Enhances aesthetics in thin acrylic for an improved appearance without compromising prosthesis contours



Highly Compatible. Zero Learning Curve

- Compatible with the leading overdenture attachment system
- Compatible with all major implant systems, including NobelActive™, Straumann® Bone Level, ZimVie Tapered Screw-Vent® and more.



1. Chih-Ho Lai a, Yin-Yu Chang b,*, Heng-Li Huang c, Ho-Yi Kao Characterization and antibacterial performance of ZrCN/amorphous carbon coatings deposited on titanium implants Available online 25 June 2011

2. Brunello G, Brun P, Gardin C, Ferroni L, Bressan E, Meneghello R, et al. (2018) Biocompatibility and antibacterial properties of zirconium nitride coating on titanium abutments: An in vitro study. PLoS ONE 13(6): e0199591. <https://doi.org/10.1371/journal.pone.0199591>



3. Internal data on file

Specifications are subject to change without notice. All product names, logos, and brands are property of their respective owners. All company, product and service names are for identification purposes only. Use of these third-party names, logos, and brands does not imply, nor should it be understood as, any indication that the respective third-party brand owners have endorsed or sponsored, or are otherwise affiliated with or connected to, Preat or Preat's goods or services.

SPECIFICATIONS

PreatLoc™ Retention Inserts

16


Standard Pivot with Dual Retention For use with two or more abutments. Accommodates up to 20° divergence. 			Try-In For denture try-ins and protection. Not for long-term use.			Extended Range For use when divergence is greater than 20°. Accommodates up to 40° divergence. 		
								
clear strong retention 5 lbs*	pink medium retention 3 lbs*	blue light retention 1.5 lbs*	gray zero retention 0 lbs	red light retention 1 lb*	orange medium retention 2 lbs*	green strong retention 4 lbs*		

The selection of the retention insert depends on the individually desired strength of the anchorage or retention.

Always start with the inserts with the least retention!

Dimensional Specifications: Abutment and Denture Cap


PREATLOC™ ABUTMENTS




PREATLOC™ BAR ABUTMENT




ABUTMENT FOR MULTI-UNIT



DENTURE CAP






PREATLOC™ ACCESSORIES

Item #	Product
7100001-4	PreatLoc™ Gray Replacement Insert (4 Pack)
7100001-20	PreatLoc™ Gray Replacement Insert (20 Pack)
7100002-4	PreatLoc™ Blue Replacement Insert (4 Pack)
7100002-20	PreatLoc™ Blue Replacement Insert 20 Pack
7100003-4	PreatLoc™ Pink Replacement Insert 4 Pack
7100003-20	PreatLoc™ Pink Replacement Insert 20 Pack
7100004-4	PreatLoc™ Clear Replacement Insert 4 Pack
7100004-20	PreatLoc™ Clear Replacement Insert 20 Pack
7100005-4	PreatLoc™ Red Replacement Insert 4 Pack
7100005-20	PreatLoc™ Red Replacement Insert 20 Pack
7100006-4	PreatLoc™ Orange Replacement Insert 4 Pack
7100006-20	PreatLoc™ Orange Replacement Insert 20 Pack
7100007-4	PreatLoc™ Green Replacement Insert 4 Pack
7100007-20	PreatLoc™ Green Replacement Insert 20 Pack
7100008-4	PreatLoc™ Black Processing Insert 4 Pack
7100008-20	PreatLoc™ Black Processing Insert 20 Pack
7100009-4	PreatLoc™ Cap Assembly - Anodized Cap 4 pack
7100009-10	PreatLoc™ Cap Assembly - Anodized Cap 10 Pack
7100010-4	PreatLoc™ 4mm Lab Analog 4 Pack
7100010-20	PreatLoc™ 4mm Lab Analog 20 Pack

ATTACHMENTS

OVERDENTURE STUD

18



PREATLOC™ ACCESSORIES

Item #	Product
7100011-4	PreatLoc™ Impression Coping 4 Pack
7100011-20	PreatLoc™ Impression Coping 20 Pack
7100012-4	PreatLoc™ Yellow Processing Insert for Bar, 4 Pack
7100012-20	PreatLoc™ Yellow Processing Insert for Bar, 20 Pack
7100013	PreatLoc™ Universal Instrument 4-piece with PEEK Abutment Holder Sleeve
7100014	PreatLoc™ Abutment Driver for Latch Handpiece with PEEK Abutment Holder Sleeve
7100015-2	PreatLoc™ Standard Pivot Processing Pack Assembly - Anodized Cap 2 Pack
7100015-10	PreatLoc™ Standard Pivot Processing Pack Assembly - Anodized Cap 10 Pack
7100016-2	PreatLoc™ Extended Pivot Processing Pack Assembly - Anodized Cap 2 Pack
7100016-10	PreatLoc™ Extended Pivot Processing Pack Assembly - Anodized Cap 10 Pack
7100017-2	PreatLoc™ Processing Pack for Bar - Anodized Cap, 2 Pack
7100017-10	PreatLoc™ Processing Pack for Bar - Anodized Cap, 10 Pack
7100018-2	PreatLoc™ Abutment for Bar, 2.0mm Thread, 2 Pack
7100018-10	PreatLoc™ Abutment for Bar, 2.0mm Thread, 10 Pack
7100019-20	PreatLoc™ Block-Out Spacer 20 Pack
7100020-4	PreatLoc™ Processing Spacer 4 Pack
7100020-20	PreatLoc™ Processing Spacer 20 Pack
7100021-4	PreatLoc™ Parallel Post 4 Pack

PREATLOC™ ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.0MM	0.75MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
BIOHORIZONS®								
Tapered Internal 3.0mm			710521	710522	710523	710524	710525	710526
Tapered Internal 3.5mm	710540 (0.5mm)		710541	710542	710543	710544	710545	
Tapered Internal 4.5mm	710550 (0.5mm)		710551	710552	710553	710554	710555	
Tapered Internal 5.7mm			710561	710562	710563	710564	710565	
BIOMET 3I®								
Certain 3.4mm			710161	710162	710163	710164	710165	
Certain 4.1mm			710171	710172	710173	710174	710175	
DENTSPLY SIRONA®								
Ankylos™ C/X (Conical Connection)				710482	710483	710484		
Astra™ Aqua			710211	710212	710213	710214		
Astra™ Lilac			710221	710222	710223	710224	710225	
Astra™ EV 3.6mm			710651	710652	710653	710654		
Astra™ EV 4.2mm			710661	710662	710663	710664	710665	
Astra™ EV 4.8mm			710671	710672	710673	710674	710675	
MIS® C1 NP			710741	710742	710743	710744		
MIS® C1 SP			710751	710752	710753	710754		
HIOSEN™								
Mini			710501	710502	710503	710504	710505	710506
Regular (Standard)			710511	710512	710513	710514	710515	710516
KEYSTONE®								
TiLobe™ 3.5mm				710572		710574		
TiLobe™ 4.1mm				710582		710584		
TiLobe™ 5.0mm				710592		710594		
MEGAGEN®								
AnyRidge™			710441	710442	710443	710444	710445	710446

ATTACHMENTS



OVERDENTURE STUD

20

PREATLOC™ ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.0MM	0.75MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
NEODENT®								
Grand Morse®		710711 (0.8mm)	710712 (1.5mm)	710713 (2.5mm)	710714 (3.5mm)	710715 (4.5mm)	710716 (5.5mm)	
NOBELBIOCARE™/NOBEL BRÅNEMARK								
Brånemark RP		710101		710102	710103	710104	710105	710106
NobelBiocare™ Tri-Lobe NP			710041	710042	710043	710044	710045	
NobelBiocare™ Tri-Lobe RP			710051	710052	710053	710054	710055	710056
NobelActive™/Conical 3.0mm			710471	710472	710473	710474	710475	710476
NobelActive™/Conical NP			710081	710082	710083	710084	710085	710086
NobelActive™/Conical RP			710091	710092	710093	710094	710095	710096
STRAUMANN®								
Bone Level NC			710151	710152	710153	710154	710155	710156
Bone Level RC			710141	710142	710143	710144	710145	710146
BLX RB/WB			710701	710702	710703	710704	710705	710706
Tissue Level RN			710121	710122	710123			
ZIMVIE®								
Tapered Screw-Vent® 3.5mm	710010 (0.5mm)		710011	710012 (2.5mm)	710013 (3.5mm)	710014 (4.5mm)	710015 (5.5mm)	710016 (6.5mm)
Tapered Screw-Vent® 4.5mm	710020 (0.5mm)		710021	710022 (2.5mm)	710023 (3.5mm)	710024 (4.5mm)	710025 (5.5mm)	
Tapered Screw-Vent® 5.7mm			710031	710032 (2.5mm)	710033 (3.5mm)	710034 (4.5mm)		

To determine the correct size O-ring needed, measure the retention zone (red line) – or the area on the male post where the O-ring seats. Select the O-ring from the size chart with the closest matching internal diameter to assure proper fit. The O-ring Combo pack is a great option as well.



		Size	Internal Diameter (mm)	Outside Diameter (mm)	Corresponding Metal Retainer (aka ring)
	Black #2 (12-Pack 1500603)	Micro	1.00mm	3.6mm	Metal Retainer #2 (1500602), Closed Metal Housing #2 MH-2 (1500601)
	Black #2M (12-Pack 1500197)	1.8 Mini Implant	1.20mm	3.60mm	
	White #3 (12-Pack 1500504)	Standard	1.40mm	4.50mm	Metal Retainer #3 (1500502). For Mini Implants Only: Closed Metal Housing MH-1 (1500501)
	Black #3 (12-Pack 1500503)	Standard	1.40mm	4.50mm	
	Black #4 (6-Pack 1500194)	Large	2.00mm	5.20mm	Metal Retainer #4 (1500204)
	Black #5 (6-Pack 1500195)	Extra Large	2.40mm	6.25mm	Metal Retainer #5 (1500205)
	Combo Pack (10-Pack 1500200)	Micro (2) 1.8 Mini Implant (2) Standard (2) Large (2) Extra Large (2)	2 of each size (Black)		See above
	J-Ring Servicing Kit 1500202)	Micro (2) 1.8 Mini Implant (2) Standard (2) Large (2) Extra Large (2) Insertion Tool (1)	2 of each size (Black) Insertion Tool		See above

ATTACHMENTS

OVERDENTURE STUD

O-Ring Attachment: features simple techniques, low cost, versatility, and is the most resilient option to protect weaker abutments 3.20mm tall 3.00mm Ø tooth surface 4.50mm Ø housing



O-Ring



1500200



4008049



1500602



1500603



1500501

O-Ring Attachment

Item #	Product
1500200	O-Ring Black O-Ring Combo Pack (10-Pack) <i>2 of each size of O-Ring (10 total): #2, #2M, #3, #4, and #5</i>
4008049	O-Ring Insertion Tool
1500202	O-Ring Servicing Kit <i>O-Ring Combo pack plus O-Ring Insertion Tool</i>
1500632	O-Ring Overdenture Complete #2 w/ housing <i>1 Castable Male, 3 black final rings, 1 yellow processing ring, 1 metal retainer ring, 1 processing analog</i>
1500630	O-Ring Castable Male #2 Micro
1500602	O-Ring Metal Retainer #2 & 2M
1500501	O-Ring Closed Metal Housing MH-1
1500601	O-Ring Closed Metal Housing #2 MH-2
1500605	O-Ring Yellow Processing #2 Micro (12-Pack)
1500603	O-Ring Black Rings #2 Micro (12-Pack)
1500606	O-Ring Overdenture Analog #2 Micro



1500502



1500505



1500503




1500504



1500506



1600025



O-Ring

O-Ring Attachment

Item #	Product
1500197	O-Ring Black Rings #2M (12-Pack)
1500532	O-Ring Overdenture Complete #3 w/ housing 1 Castable Male, 3 white final rings, 1 red processing ring, 1 metal retainer ring, 1 processing analog
1500530	O-Ring Castable Male Standard #3
1500502	O-Ring Metal Retainer Standard #3
1500505	O-Ring Red Processing Rings #3 (12-Pack)
1500503	O-Ring Black Rings Standard #3 (12-Pack)
1500504	O-Ring White Rings Standard #3 (12-Pack)
1500506	O-Ring Overdenture Analog Standard #3
1600025	O-Ring Implant Analog
1500204	O-Ring Metal Retainer #4
1500194	O-Ring Black Rings #4 (6-Pack)
1500214	O-Ring Overdenture Analog #4
1500205	O-Ring Metal Retainer #5
1500195	O-Ring Black Rings #5 (6-Pack)

23

O-Ring Attachment

O-Ring Mini Implant Attachment: compatible with mini implant systems featuring a square base and 1.8mm ball top.



1500610



1500212

O-Ring Mini Implant Attachment

Item #	Product
1500610	O-Ring Mini Implant Servicing Kit 1 impression coping, 1 model analog, 1 MH-2 closed housing, 1 yellow processing O-ring, 2 black O-rings
1500212	O-Ring Mini Implant Impression Copings
1500210	O-Ring Mini Implant Analog #2M

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

O-Ring Attachment: the extremely resilient titanium O-Ring attachment requires 4.0mm vertical, 5.0mm Ø All abutments include metal retainer ring, 1 red processing ring, and 3 black rubber rings



24

O-Ring Abutment



O-Ring Abutment

Item #	Product
1600801C	O-Ring Abutment Zimmer TSV 3.5 X 1mm
1600802C	O-Ring Abutment Zimmer TSV 3.5 X 2mm
1600803C	O-Ring Abutment Zimmer TSV 3.5 X 3mm
1600804C	O-Ring Abutment Zimmer TSV 3.5 X 4mm
1600805C	O-Ring Abutment Zimmer TSV 3.5 X 5mm
1600806C	O-Ring Abutment Zimmer TSV 3.5 X 6mm
1600811C	O-Ring Abutment Zimmer TSV 4.5 X .5mm
1600812C	O-Ring Abutment Zimmer TSV 4.5 X 2mm
1600813C	O-Ring Abutment Zimmer TSV 4.5 X 3mm
1600814C	O-Ring Abutment Zimmer TSV 4.5 X 4mm
1600815C	O-Ring Abutment Zimmer TSV 4.5 X 5mm
1600816C	O-Ring Abutment Zimmer TSV 4.5 X 6mm
1600821C	O-Ring Abutment Branemark RP x 1mm
1600822C	O-Ring Abutment Branemark RP x 2mm
1600823C	O-Ring Abutment Branemark RP x 3mm



O-Ring Abutment	
Item #	Product
1600824C	O-Ring Abutment Branemark RP x 4mm
1600825C	O-Ring Abutment Branemark RP x 5mm
1600826C	O-Ring Abutment Branemark RP x 6mm
1500502	O-Ring Metal Retainer Standard #3 (each)
1500505	O-Ring Red Processing Rings #3 (12-Pack)
1500503	O-Ring Black Rings Standard #3 (12-Pack)
1500504	O-Ring White Rings Standard #3 (12-Pack)
1600025	O-Ring Implant Analog
1600030	O-Ring/Clx Latch Driver
1600035	Latch to Square Driver Insert
4008049	O-Ring Insertion Tool

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

Clix® Attachment: 2.25mm spherical attachment with slotted plastic female to engage a larger surface area of the sphere and reduce wear. This unique female makes an audible click to ensure prosthesis security. Housing may be paralleled to correct up to 30° divergence requires: 4.0mm vertical, 4.0mm Ø

26

Clix®



Clix®

Item #

Product

1608001C Clix Complete Ball Abutment for Nobel Active NP x 1mm

1608002C Clix Complete Ball Abutment for Nobel Active NP x 2mm

1608003C Clix Complete Ball Abutment for Nobel Active NP x 3mm

1608004C Clix Complete Ball Abutment for Nobel Active NP x 4mm

1608005C Clix Complete Ball Abutment for Nobel Active NP x 5mm

1608006C Clix Complete Ball Abutment for Nobel Active NP x 6mm

1600901C Clix Complete Ball Abutment for Nobel Active RP x 1mm

1600902C Clix Complete Ball Abutment for Nobel Active RP x 2mm

1600903C Clix Complete Ball Abutment for Nobel Active RP x 3mm


1600904C Clix Complete Ball Abutment for Nobel Active RP x 4mm

1600905C Clix Complete Ball Abutment for Nobel Active RP x 5mm

1600906C Clix Complete Ball Abutment for Nobel Active RP x 6mm

1601500C Clix Complete Ball Abutment for Straumann Bone Level NC x 1mm

1601501C Clix Complete Ball Abutment for Straumann Bone Level NC x 2mm



Clix®

Item #	Product
1601502C	Clix Complete Ball Abutment for Straumann Bone Level NC x 3mm
1601503C	Clix Complete Ball Abutment for Straumann Bone Level NC x 4mm
1601504C	Clix Complete Ball Abutment for Straumann Bone Level NC x 5mm
1601505C	Clix Complete Ball Abutment for Straumann Bone Level NC x 6mm
1601400C	Clix Complete Ball Abutment for Straumann Bone Level RC x 1mm
1601401C	Clix Complete Ball Abutment for Straumann Bone Level RC x 2mm
1601402C	Clix Complete Ball Abutment for Straumann Bone Level RC x 3mm
1601403C	Clix Complete Ball Abutment for Straumann Bone Level RC x 4mm
1601404C	Clix Complete Ball Abutment for Straumann Bone Level RC x 5mm
1601405C	Clix Complete Ball Abutment for Straumann Bone Level RC x 6mm
1600100C	Clix Complete Ball Abutment for Zimmer TSV NP x 1mm
1600102C	Clix Complete Ball Abutment for Zimmer TSV NP x 2mm
1600103C	Clix Complete Ball Abutment for Zimmer TSV NP x 3mm
1600104C	Clix Complete Ball Abutment for Zimmer TSV NP x 4mm
1600105C	Clix Complete Ball Abutment for Zimmer TSV NP x 5mm
1600106C	Clix Complete Ball Abutment for Zimmer TSV NP x 6mm

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

28

Clix®



Clix®

Item #	Product
1600200C	Clix Complete Ball Abutment for Zimmer TSV RP x 1mm
1600202C	Clix Complete Ball Abutment for Zimmer TSV RP x 2mm
1600203C	Clix Complete Ball Abutment for Zimmer TSV RP x 3mm
1600204C	Clix Complete Ball Abutment for Zimmer TSV RP x 4mm
1600205C	Clix Complete Ball Abutment for Zimmer TSV RP x 5mm
1600206C	Clix Complete Ball Abutment for Zimmer TSV RP x 6mm
1600402C	Clix Complete Ball Abutment for Tri-Lobe NP x 2mm
1600403C	Clix Complete Ball Abutment for Tri-Lobe NP x 3mm
1600404C	Clix Complete Ball Abutment for Tri-Lobe NP x 4mm
1600405C	Clix Complete Ball Abutment for Tri-Lobe NP x 5mm
1600406C	Clix Complete Ball Abutment for Tri-Lobe NP x 6mm
1600500C	Clix Complete Ball Abutment for Tri-Lobe RP x 0mm



Clix®




Item #	Product
1600502C	Clix Complete Ball Abutment for Tri-Lobe RP x 2mm
1600503C	Clix Complete Ball Abutment for Tri-Lobe RP x 3mm
1600504C	Clix Complete Ball Abutment for Tri-Lobe RP x 4mm
1600505C	Clix Complete Ball Abutment for Tri-Lobe RP x 5mm
1600506C	Clix Complete Ball Abutment for Tri-Lobe RP x 6mm
1602102C	Clix Complete Ball Abutment for Astra Aqua 3.5/4.0mm x 2mm
1602104C	Clix Complete Ball Abutment for Astra Aqua 3.5/4.0mm x 4mm
1602202C	Clix Complete Ball Abutment for Astra Lilac 4.5/5.0mm x 2mm
1602204C	Clix Complete Ball Abutment for Astra Lilac 4.5/5.0mm x 4mm
1601602C	Clix Complete Ball Abutment for Biomet 3i Certain 3.4mm x 2mm
1601604C	Clix Complete Ball Abutment for Biomet 3i Certain 3.4mm x 4mm
1601702C	Clix Complete Ball Abutment for Biomet 3i Certain 4.1mm x 2mm
1601704C	Clix Complete Ball Abutment for Biomet 3i Certain 4.1mm x 4mm

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

30


Clix®



Clix®

Item #	Product
1600303C	Clix Complete Ball Abutment for Straumann RN x 2mm
1600304C	Clix Complete Ball Abutment for Straumann RN x 3mm
1600305C	Clix Complete Ball Abutment for Straumann RN x 4mm
1601002C	Clix Complete Ball Abutment for Neoss 3.5-6.0mm x 2mm
1601004C	Clix Complete Ball Abutment for Neoss 3.5-6.0mm x 4mm
1600020	Clix Implant Ball Analog
1600030	Clix Latch Driver
1600035	Latch to Square Driver Insert
011222	Clix Female Insertion Tool
0101240	Clix Complete Female TI housing and yellow female
0101232W-6	Clix White Females- Reduced (6-Pack) white = 1.7lbs retention
0101232W-50	Clix White Females- Reduced (50-Pack)
0101231-6	Clix Yellow Females- Standard (6-Pack) yellow = 2.5 lbs retention
0101231-50	Clix Yellow Females- Standard (50-Pack)
0101233-6	Clix Red Females- Increased (6-Pack) red = 3.3 lbs retention
0101233-50	Clix Red Females- Increased (50-Pack)
0112340	Clix Combo Females 2 white, 2 yellow, 2 red Clix Females

LOCATOR®: the most popular stud attachment in the world due to its angle correction, ease of use, lowprofile 9 different retention levels to meet patient's needs, and patented PIVOT between the housing and nylon male that reduces potential wear on the abutment. *Beware imitations that do not pivot* 2.73mm vertical, 5.45mm Ø



LOCATOR® Self Aligning

Item #	Product
4008260	LOCATOR® 35n Square Driver (15mm) <i>For use with Straumann torque wrenches</i>
4008280	LOCATOR® 35n Square Driver (21mm) <i>For use with Straumann torque wrenches</i>
4008913	LOCATOR® Insert Driver (Latch 23mm)
4008914	LOCATOR® Insert Driver (Latch 29mm)
4008926	LOCATOR® Square Insert Driver (15mm)
4008927	LOCATOR® Square Insert Driver (21mm)
4008393	LOCATOR® Core Tool
4008390	LOCATOR® Gold Driver
4008397	LOCATOR® Male Removal Tool
4008394	LOCATOR® Abutment Carrier (4-Pack)
4009006	LOCATOR® Demo Model
4009014	LOCATOR® Root Model
4009015	LOCATOR® Implant Model w/External Hexed Implants
4008505	LOCATOR® Impression Coping (4-Pack)
4008505-20	LOCATOR® Impression Coping (20-Pack)
4008530	LOCATOR® Female Analog (4.0mm)(4-Pack) <i>The 4mm ȳ analog is for Narrow and Regular Platform Implants</i>
4008530-20	LOCATOR® Female Analog (4.0mm)(20-Pack)
4008514	LOCATOR® Block Out Spacer (20-Pack)
4008515	LOCATOR® Processing Replacement Male (Black)(4-Pack)

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS


32

LOCATOR® Self Aligning



LOCATOR® Self Aligning

Item #	Product
4008515-20	LOCATOR® Processing Replacement Male (Black) (20-Pack)
4008516	LOCATOR® Female Analog (5.0mm)(4-Pack) <i>The 5mm Ø analog is for Wide Platform Implants & Natural Teeth</i>
4008516-20	LOCATOR® Female Analog (5.0mm)(20-Pack)
4008517	LOCATOR® Parallel Post (4-Pack)
4008524	LOCATOR® Replacement Male (Clear)(4-Pack)
4008524-20	LOCATOR® Replacement Male (Clear)(20-Pack)
4008527	LOCATOR® Light Retention Replacement Male (Pink) (4-Pack)
4008527-20	LOCATOR® Light Retention Replacement Male (Pink) (20-Pack)
4008529	LOCATOR® Extra Light Retention Replacement Male (Blue)(4-Pack)
4008529-20	LOCATOR® Extra Light Retention Replacement Male (Blue)(20-Pack)
4008547	LOCATOR® Extended Range Replacement Male (Green)(4-Pack)
4008547-20	LOCATOR® Extended Range Replacement Male (Green)(20-Pack)
4008548	LOCATOR® Extra Light Retention Extended Range Replacement Male (Red)(4-Pack)
4008548-20	LOCATOR® Extra Light Retention Extended Range Replacement Male (Red)(20-Pack)
4008915	LOCATOR® Extended Range Replacement Male (Orange)(4-Pack)
4008915-20	LOCATOR® Extended Range Replacement Male (Orange)(20-Pack)
4008558	LOCATOR® Zero Retention Extended Range Replacement Male (Gray)(4-Pack)
4008558-20	LOCATOR® Zero Retention Extended Range Replacement Male (Gray)(20-Pack)
4008559	LOCATOR® High Retention Replacement Male (Purple) (4-Pack)
4008026	LOCATOR® Bar Processing Replacement Male (Yellow) (4-Pack)
4008026-20	LOCATOR® Bar Processing Replacement Male (Yellow) (20-Pack)
4008569	LOCATOR® Processing Spacer (4-Pack)



LOCATOR® Self Aligning

Item #	Product
4008569-20	LOCATOR® Processing Spacer (20-Pack)
4008510-4	LOCATOR® Replacement Denture Cap Male Assembly (4-Pack)
4008510-10	LOCATOR® Replacement Denture Cap Male Assembly (10-Pack)
4008519-2	LOCATOR® Male Processing Package (2-Pack)
4008519-10	LOCATOR® Male Processing Package (10-Pack)
4008540-2	LOCATOR® Extended Range Male Processing Package (2-Pack)
4008540-10	LOCATOR® Extended Range Male Processing Package (10-Pack)
4008550-2	LOCATOR® Male Processing Package SS (2-Pack) <i>Stainless Steel housing ideal for welding to CoCr framework, SS housing, blackout spacer, blue / pink / clear males</i>
4008550-10	LOCATOR® Male Processing Package SS (10-Pack)
4008028-2	LOCATOR® Bar Male Processing (2-Pack) <i>Ti housing. w/ yellow, 2 block. spacers, blue / pink / clear</i>
4008028-10	LOCATOR® Bar Male Processing (10-Pack)
4008010	LOCATOR® Refill Kit <i>Four 0° females, 2 parallel posts, 4 complete males</i>
4008011	LOCATOR® Starter Kit <i>Spotface Drill, Pilot Drill, Four 0° females, 4 complete males, 2 parallel posts</i>
4008050	LOCATOR® Refill Coping Kit <i>4 Cast-to females w/parallel posts, 4 complete males</i>
4008051	LOCATOR® Coping Kit <i>Spotface Drill, Pilot Drill, 4 Cast-to females w/parallel posts, 4 complete males</i>
4008070	LOCATOR® Refill Mixed Kit <i>Two 0°, One 10°, One 20° females, 2 parallel posts, 4 complete males</i>

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

34

LOCATOR® Self Aligning



LOCATOR® Self Aligning

Item #	Product
4008071	LOCATOR® Mixed Kit <i>Spotface Drill, Pilot Drill, Two 0°, One 10°, One 20° females, 2 parallel posts, 4 complete males</i>
4008922	LOCATOR® Spotface Diamond Bur
4008924	LOCATOR® Pilot Drill
4008520-2	LOCATOR® 0 Degree Female (2-Pack)
4008520-10	LOCATOR® 0 Degree Female (10-Pack)
4008521-2	LOCATOR® 10 Degree Female (2-Pack)
4001747-2	LOCATOR® Retrofit Kit for Worn ERA Root Attachments (2-Pack) <i>Includes LOCATOR® Retrofit, Parallel Post, Weep Hole Alignment Pin</i>
4008521-10	LOCATOR® 10 Degree Female (10-Pack)
4008522-2	LOCATOR® 20 Degree Female (2-Pack)
4008522-10	LOCATOR® 20 Degree Female (10-Pack)
4008528-2	LOCATOR® 0 Degree Cast-to Female (2-Pack)
4008528-10	LOCATOR® 0 Degree Cast-to Female (10-Pack)
4008586-2	LOCATOR® Cast-To Bar Female (2-Pack)
4008586-10	LOCATOR® Cast-To Bar Female (10-Pack)
4008587-2	LOCATOR® Bar Female 2-56 Thread (2-Pack)
4008587-10	LOCATOR® Bar Female 2-56 Thread (10-Pack)
4008588-2	LOCATOR Laser Bar Female, Titanium (2-Pack)
4008588-10	LOCATOR® Laser Bar Female, Titanium (10-Pack)
4008589-2	LOCATOR® Bar Female, 2.0mm Thread (2-Pack)



4008014



4008016



4009102



4009104



4009105



4009107



4009260-2



4008014

LOCATOR® Self Aligning

Item #	Product
4008589-10	LOCATOR® Bar Female, 2.0mm Thread (10-Pack)
4008590-2	LOCATOR® Laser Bar Female, Stainless Steel (2-Pack)
4008590-10	LOCATOR® Laser Bar Female, Stainless Steel (10-Pack)
4008013	LOCATOR® 2-56 Castable sleeve (10-Pack)
4008014	LOCATOR® 2mm Castable Threaded Inserts (10-Pack)
4008016	LOCATOR® Tap Holder
4009020	30N-cm Torque and Thumb Knob Wrench Kit
4009102	1.7mm Bar Drill (2.0mm Thread)
4009103	1.8mm Bar Drill (2-56 Thread)
4009104	2.0mm Bar Tap (2.0mm Thread)
4009105	LOCATOR® 2-56 Bar Tap (2-56 Thread)
4009107	LOCATOR® Paralleling Mandrel for Bar Female
4008909-2	LOCATOR® for Multi-Unit Abutment TI (2-Pack)
4009260-2	LOCATOR Scan Body (2-Pack)
4009260-4	LOCATOR® Scan-body (4-Pack)
4009260-10	LOCATOR Scan Body (10-Pack)

ATTACHMENTS


LOCATOR® RETENTION INSERTS

Attachments




Dimensional Specifications


LOCATOR® IMPLANT ATTACHMENTS




LOCATOR® DENTURE CAP & BLACK PROCESSING MALE



LOCATOR® BAR ATTACHMENTS




LOCATOR® ROOT ATTACHMENTS



LOCATOR® PIVOTING FUNCTION



LOCATOR® Multi-Unit



IMPLANT ABUTMENTS/ATTACHMENTS

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.0MM	0.73MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
AB DENTAL®								
3.75mmD (NP)	8661		8662	8663 2.5mm	8664 3.5mm	8665 4.5mm	8625 5.5mm	8626 6.5mm
4.5mmD (WP)	8671		8672	8673	8674	8675	8627	8628
ALMITECH®								
3.3mm/4.1mm/4.8mm	1438		1439	1440	1441	1442	1443	
ALPHA BIO IMPLANT®								
3.5mmD	8661		8662	8663 2.5mm	8664 3.5mm	8665 4.5mm	8625 5.5mm	8626 6.5mm
3.75mmD		1192	1193	1194	1195	1196	1197	1198
ANTHOGYR®								
3.4mm - Short			1697 1.55mm					
3.4mm - Medium				1698 2.5mm				
3.4mm - Tall					1699 3.5mm			
ASTRA®								
TX 3.5mm/4.0mmD			8115	8116	8117	8118	8119	
TX 4.5mm/5.0mmD		8686 0.5mm		8687	8688	8689	8690	
M2 Thread (old style) 3.5mmD		8636 0.81mm		8637	8638	8639	8640	
M2 Thread (old style) 4.0mmD ST		8656 0.5mm		8657	8658	8659	8660	
EV 3.6mmD			2790	2791	2792	2793	2794	
EV 4.2mmD			2795	2796	2797	2798	2799	
EV 4.8mmD			2800	2801	2802	2803	2804	
BICON®								
2.0 Post			8746	8747	8748	8749	8750	
2.5mm			2278	2279	2280	2281	2282	
3.0 Post			8930	8931	8932	8933	8934	
BIOHORIZONS®								
External 3.5mmD (Ex Hex)			8531 1.25mm	8532	8533	8534	8535	
External 4.0mmD (Ex Hex)			8541	8542	8543	8544	8545	
External 5.0mmD (Ex Hex)			8551 1.5mm	8552 2.5mm	8553 3.5mm	8554 4.5mm	8555 5.5mm	
Tapered Internal 3.0mmD			1927	1928	1929	1930		
Tapered Internal 3.5mm/3.8mmD	8661		8662	8663 2.5mm	8664 3.5mm	8665 4.5mm	8625 5.5mm	8626 6.5mm
Tapered Internal 4.0mm/4.5mm/4.6mmD	8671		8672	8673	8674	8675	8627	8628
Tapered Internal 5.7mm/5.8mmD			8190	8191	8192	8193	8194	8195

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.0MM	0.73MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
BIOMET 3I®								
Certain® 3.4mmD			8431	8432	8433	8434	8435	8436
Certain® 4.1mmD (Compatible with 5.0mm and 6.0mm)		8630		8631	8632	8633	8634	8635
MicroMiniPlant 3.4mmD (Ex Hex)			8441	8442	8443	8444	8445	8446
MiniPlant 3.75mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
Cylindrical 4.0mm Diameter (Ex Hex)		8681		8682	8683	8684	8685	8696
External 5.0mm Diameter (Ex Hex)		8820 0.6mm		8821	8822	8823	8824	
TG (Internal Connection)			8621 1.25mm	8622 2.25mm	8623 3.25mm	8624 4.25mm		
BIOTECH DENTAL KONTACT®								
			1792		1794		1796	
BLUE SKY BIO®								
Internal Hex 3.7mmD	8661		8662	8663 2.5mm	8664 3.5mm	8665 4.5mm	8625 5.5mm	8626 6.5mm
Internal Hex 4.7mmD	8671		8672	8673	8674	8675	8627	8628
Trilobe 3.5mmD (NP)	8756		8757	8758	8759	8760	8734	8735
Trilobe 4.3mmD (RP)	8761		8762	8763	8764	8765	8753	8754
Trilobe 5.0mmD (WP)			8766	8767	8768	8769	8770	
One Stage Regular		8950		8951	8952	8953	8954	8955
BREDDENT®								
4.1mm			1827	1828	1829	1830	1831	1832
BTI - BIOTECHNOLOGY INSTITUTE SL®								
Internal Hex - 4.1mm			8291	8292	8293	8294	8295	8296
External Hex 4.1mmD		8681		8682	8683	8684	8685	8696
External Hex 5.5mmD				8302	8303	8304	8305	8306
DENTEGRIZ DENTAL TECH®								
3.75mm	8052		8053	8054	8055	8056	8057	
4.5mm		8058 0.63mm		8059	8060	8061	8062	8063
4.8mm Green/Gray/Blue		8950		8951	8952	8953	8954	8955
DENTIUM USA®								
Implantium / Superline 3.4mm/3.8mm, 4.3mm/4.8mm, 4.8mmD		8656 0.5mm		8657	8658	8659	8660	
DENTSPLY / FRIADENT®								
Ankylos® C/X (Conical Connection)				1729	1730	1731	1732	1733
Frialit-2 / Xive 3.4mmD			8805	8806	8807	8808	8809	
Frialit-2 / Xive 3.8mmD			8810	8811	8812	8813	8814	

IMPLANT ABUTMENTS/ATTACHMENTS

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.0MM	0.73MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
DENTSPLY® / FRIADENT® cont'd								
Frialit-2 / Xive 4.5mmD			8815	8816	8817	8818	8819	
Frialit-2 / Xive 5.5mmD			8139 1.3mm	8140	8141	8142	8143	
DIO®								
Dio SM 3.1mm SM (Narrow)			1351	1352	1353	1354	1355	1356
Dio SM 3.5mm (Regular & Wide)			1361	1362	1363	1364	1365	1366
IMTEC®								
Endure® 3.5mm/4.3mm (Internal Hex)		8994 0.63mm		8995	8996	8997	8998	
Self Tapping 3.75mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
Universal 3.4mm/4.0mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
Mini 3.3mmD (Ex Hex)		8856 0.91mm		8857	8858	8859	8860	
IMPLANT DIRECT®								
InterActive / SwishActive™ 3.0mmD (NP)			2069	2070	2071	2072	2073	2074
InterActive / SwishActive™ 3.4mmD (RP)			2075	2076	2077	2078	2079	2080
Legacy™ 3.0mmD (Internal Hex)		2308 0.75mm		2309	2310	2311	2312	2313
Legacy™ 3.5mmD (Internal Hex)	8661		8662	8663 2.5mm	8664 3.5mm	8665 4.5mm	8625 5.5mm	8626 6.5mm
Legacy™ 4.5mmD (Internal Hex)	8671		8672	8673	8674	8675	8627	8628
Legacy™ 5.7mmD (Internal Hex)			1471	1472	1473	1474	1475	1476
RePlus / RePlant / ReActive™ 3.5mmD (Trilobe NP)	8756		8757	8758	8759	8760	8734	8735
RePlus / RePlant / ReActive™ 4.3mmD (Trilobe RP)	8761		8762	8763	8764	8765	8753	8754
RePlus / RePlant / ReActive™ 5.0mmD (Trilobe WP)			8766	8767	8768	8769	8770	
ScrewPlant™ 3.7mmD			1371	1372	1373	1374	1375	1376
ScrewPlant™ 4.7mmD			1381 1.5mm	1382	1383	1384	1385	1386
ScrewPlant™ 5.7mmD				1391	1392	1393	1394	1395
ScrewInDirect™ 5.0mmD				1250	1251	1252	1253	1254
SwishTapered / SwishPlus™ 3.7mmD		1821 0.4mm	1822	1823	1824	1825	1826	
SwishTapered / SwishPlus™ 4.8mmD		8950		8951	8952	8953	8954	8955
IMZ®								
3.3mm/4.0mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
3.3mmD (Non-Hex)	8020		8021	8022	8023	8024	8025	
4.0mmD (Non-Hex)	8030 0.35mm		8031	8032	8033	8034	8035	
INTRA-LOCK SYSTEM INTERNATIONAL®								
Threaded 3.75mm/4.0mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
Threaded Conic 4.3mmD (Ex Hex)			8781	8782	8783	8784	8785	

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.0MM	0.73MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
INTRA-LOCK SYSTEM INTERNATIONAL® cont'd								
Threaded 4.75mmD (Ex Hex)			8727 1.42mm	8728	8729	8730	8731	
Threaded Conic 5.0mmD (Ex Hex)			8727 1.42mm	8728	8729	8730	8731	
KEYSTONE® / LIFECORE BIOMEDICAL®								
3.3mm/4.0mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
Stage-1 Regular Diameter		8950		8951	8952	8953	8954	8955
Renova 3.75mmD (Internal Connection)	8090		8091	8092	8093	8094	8095	8096
Revona 4.5mm/4.75mmD (Internal Connection)		8125 0.63mm		8126	8127	8128	8129	8130
Restore SD 3.3mm/3.4mmD (Ex Hex)			8196	8197	8198	8199	8200	
Restore RD 3.75mm/4.0mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
Restore WD 5.0mm/6.0mmD (Ex Hex)			8727 1.42mm	8728	8729	8730	8731	
Sustain 4.2mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
Sustain 4.0mm/4.7mmD	8651		8652	8653	8654	8655		
KYOCERA POI®								
4.2mm (Non-Hexed)		8716 0.35mm	8717	8718-1	8719	8720 4.5mm	8721	
LASAK®								
3.7mm	1207		1208	1209	1210	1211	1212	
MEGAGEN® / INTEGRATED DENTAL SYSTEMS®								
AnyRidge® (Internal Connection)			2100	2101	2102	2103	2104	2105
MIS - MEDICAL IMPLANT SYSTEM®								
Internal Hex 3.75mm/4.2mmD (SP)	8661		8662	8663 2.5mm	8664 3.5mm	8665 4.5mm	8625 5.5mm	8626 6.5mm
Internal Hex 5.0mmD (WP)	8671		8672	8673	8674	8675	8627	8628
External SP (Ex Hex)		8681		8682	8683	8684	8685	8696
External WP (Ex Hex)		8820 0.6mm		8821	8822	8823	8824	
Mistral Internal Octagon (Stage One)		8950		8951	8952	8953	8954	8955
SEVEN 3.3 (Internal Hex 3.3mm NP)	2062 0.1mm		2063	2064	2065	2066	2067	2068
NEOSS®								
4.0mmD (Non-Hex)			8451	8452	8453	8454	8455	8456
NEODENT®								
Cone Morse				2168	2169	2170	2171	2172
Grand Morse			2112	2113	2114	2115	2116	2117
NOBEL BIO CARE® / BRANEMARK® / STERIOSS®								
NobelActive™ / Conical Connection 3.0mmD			2212	2213	2214	2215	2216	2217
NobelActive™ / Conical Connection NP (3.5mmD)			2069	2070	2071	2072	2073	2074
NobelActive™ / Conical Connection RP (4.3mm/5.0mmD)			2075	2076	2077	2078	2079	2080

IMPLANT ABUTMENTS/ATTACHMENTS

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.0MM	0.73MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
NOBEL BIO CARE® / BRANEMARK® / STERIOSS® cont'd								
Steri-Oss Flat-Top 3.8mmD	8651		8652	8653	8654	8655		
Replace Select NP (3.5mmD, Trilobe)	8756		8757	8758	8759	8760	8734	8735
Replace Select RP (4.3mmD, Trilobe)	8761		8762	8763	8764	8765	8753	8754
Replace Select WP (5.0mmD, Trilobe)			8766	8767	8768	8769	8770	
Replace Select 6.0mmD (Trilobe)			8766	8767	8768	8769	8770	
Branemark NP (3.3mmD, Ex Hex)		8105		8106	8107	8108	8109	
Branemark RP (3.75mm/4.0mm/5.0mmD, Ex Hex)		8681		8682	8683	8684	8685	8696
Branemark WP (5.0mmD, Ex Hex)			8825 1.1mm	8826	8827	8828	8829	
Nobel Speedy NP (3.5mmD, Ex Hex)		8105		8106	8107	8108	8109	
Nobel Speedy RP (4.3mmD, Ex Hex)		8681		8682	8683	8684	8685	8696
Replace Ex Hex NP (3.5mmD)			8676	8677	8678	8679	8680	
Replace Ex Hex RP (4.3mmD)			8781	8782	8783	8784	8785	
Replace Ex Hex WP (5.0mmD)			8835 1.37mm	8836	8837	8838	8839	
Hex-Lock (HL) 3.25mmD (Ex Hex)			8646 1.5mm	8647 2.5mm	8648 3.5mm	8649 4.5mm	8650 5.5mm	
Hex-Lock (HL) 3.8mm/4.5mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
Hex-Lock (HL) 5.0mmD (Ex Hex)			8835 1.37mm	8836	8837	8838	8839	
OSSTEM® / HIOSSEN USA® / DENTIS®								
SS Internal / HS Regular 4.8mmD		8950		8951	8952	8953	8954	8955
ET/ HG/ GS Mini 3.5mmD (Internal Connection)			1561	1562	1563	1564	1565	1566
ET/ HG/ GS Regular 4.0mmD (Internal Connection)			1551	1552	1553	1554	1555	1556
US Mini 3.3mm / 3.5mmD (Ex Hex)			8105	8106	8107	8108	8109	
US Regular/ HU Regular 3.75mm / 4.0mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
OSTEO-IMPLANT CORPORATION®								
3.25mm SD, 3.75mm ST, 4.0mm DT		8681		8682	8683	8684	8685	8696
PHIBO® / IMPLADENT®								
TSA S3 – 3.7mm			1755	1756	1757	1758	1759	1760
TSA S4 – 4.7mm			1761 1.4mm	1762	1763	1764	1765	1766
TSA Advanced S5 – 6.0mm				1767	1768	1769	1770	1771
TSH S2 – 3.3mm			1916	1917	1918	1919	1920	1921
TSH S3 & S4 – 4.0mm			1942	1943	1944	1945	1946	1947
TSH S5 – 5.0mm				1922	1923	1924	1925	1926
RIEMSER®								
Arzneimittel Revois – 3.8/4.3/5.0mm			1213	1214	1215	1216	1217	1218
Arzneimittel Revois II – 3.8/4.3/5.0mm	1303 0.2mm		1304	1305	1306	1307	1308	

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.0MM	0.73MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
RITTER IMPLANTS®								
Spiral/QSI NP (Internal Hex)		2308 0.75mm		2309	2310	2311	2312	2313
Spiral/QSI RP (Internal Hex)			2284	2285	2286	2287	2288	2289
SOUTHERN IMPLANT USA®								
Tri-Nex 3.5mmD	8756		8757	8758	8759	8760	8734	8735
Tri-Nex 4.3mmD	8761		8762	8763	8764	8765	8753	8754
Tri-Nex 5.0mmD			8766	8767	8768	8769	8770	
Tri-Nex 6.0mmD			8766	8767	8768	8769	8770	
IT 4.8mmD (Internal Octagon)		8950		8951	8952	8953	8954	8955
STERNGOLD IMPLAMED®								
Hex Cylinder 3.3mm/4.0mmD		8681		8682	8683	8684	8685	8696
Hex Screw 3.75mm/4.0mmD		8681		8682	8683	8684	8685	8696
STRAUMANN®								
Tissue Level 3.5mmD (NNC)			2056	2057	2058	2059	2060	2061
Tissue Level 4.8mmD (RN)			8742	8743	8744	8745	8751	8752
Tissue Level 6.5mmD (WN)			8870	8871	8872	8873	8874	
Bone Level 3.3mmD (NC)				8422	8423	8424	8425	8426
Bone Level 4.1mmD (RC)			8170	8171	8172	8173	8174	8175
BLX Regular Base (RB) & Wide Base (WB)			2326	2327	2328	2329	2330	2331
SYBRON® / INNOVA®								
Endopore 3.5mmD (0.9mm Ex Hex)			8851	8852	8853	8854	8855	
Endopore 4.1mmD (0.7mm Ex Hex)		8681		8682	8683	8684	8685	8696
Endopore 5.0mmD (0.9mm Ex Hex)			8691 1.35mm	8692	8693	8694	8695	
Endopore 4.1mm/5.0mmD (One Stage)		8950		8951	8952	8953	8954	8955
Endopore 4.1mmD (Internal Hex)	8145		8146	8147	8148	8149	8150	
Endopore 5.0mmD (Internal Hex)			8164	8165	8166	8167	8168	8169
Entegra 3.25mm/4.0mmD (0.7mm Ex Hex)		8681		8682	8683	8684	8685	8696
Entegra 3.25mmD (0.9mm Ex Hex)			8851	8852	8853	8854	8855	
Entegra 4.7mmD (0.9mm Ex Hex)			8691 1.35mm	8692	8693	8694	8695	
Pitt Easy 3.25mm/3.3mmD			8311	8312	8313	8314	8315	8316
Pitt Easy 3.75mm/3.8mmD			8321	8322	8323	8324	8325	8326
Pitt Easy 4.0mm/4.1mmD			8331	8332	8333	8334	8335	8336
Pitt Easy 4.9mmD			8341	8342	8343	8344	8345	8346
Pro Series XRT Bone Level 3.5mmD			1631	1632	1633	1634	1635	1636
Pro Series XRT Bone Level 4.1mmD			1637	1638	1639	1640	1641	1642
Pro Series XRT Bone Level 4.8mmD		8950		8951	8952	8953	8954	8955

IMPLANT ABUTMENTS/ATTACHMENTS

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.0MM	0.73MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
SYBRON® / INNOVA® CONT'D								
Pro Series Tissue Level (TL) 4.1mm/4.8mmD		8950		8951	8952	8953	8954	8955
Hex Plus 3.5mmD		8105		8106	8107	8108	8109	
Hex Plus 4.1mmD		8681		8682	8683	8684	8685	8696
Hex Plus 5.1mmD			8825 1.1mm	8826	8827	8828	8829	
SARGON®								
External Hex 3.8mmD			8722 1.8mm		8723	8724	8725	8726
TEKKA®								
I-Kone		1937 .70mm	1938 1.5mm		1939		1940	1941 7.0mm
Krestal – 3.2mm			1691	1692	1693	1694	1695	1696
Krestal – 3.4mm	1220		1221	1222	1223	1224	1225	
Krestal – 4.0mm			1226	1227	1228	1229	1230	1231
Krestal – 4.5mm			1232	1233	1234	1235	1236	1237
Krestal – 5.0mm			1238	1239	1240	1241	1242	1243
Krestal – 6.0mm			1244	1245	1246	1247	1248	1249
ZIMMER® / CENTERPULSE® / CALCITEK® / PARAGON®								
Advent 3.7mm/4.7mmD (4.5mmD Platform)		8841 0.71mm		8842	8843	8844	8845	8846
Screw-Vent/ Tapered Screw Vent 3.3mm/3.7mm/4.1mmD (3.5mmD)	8661		8662	8663 2.5mm	8664 3.5mm	8665 4.5mm	8625 5.5mm	8626 6.5mm
Screw-Vent / Tapered Screw Vent 4.7mmD (4.5mmD)	8671		8672	8673	8674	8675	8627	8628
Screw-Vent/ Tapered Screw Vent 6.0mm (5.7mmD)			8180	8181	8182	8183	8184	8185
Calcitek: Integral SD 3.25mmD	8611		8612	8613	8614	8615		
Calcitek: Integral 4.0mmD	8651		8652	8653	8654	8655		
Calcitek: Omniloc 3.25mmD	8611		8612	8613	8614	8615		
Calcitek: Omniloc 4.0mmD	8666		8667	8668	8669	8670		
Calcitek: Spline 3.25mmD			8701	8702	8703	8704	8705	
Calcitek: Spline 3.75mm/4.0mmD			8711	8712	8713	8714	8715	
Calcitek: Spline 5.0mmD		8736 0.86mm		8737	8738	8739	8740	
Calcitek: Swede-Vent 4.0mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
Calcitek: Swede-Vent TL (External Hex)		8681		8682	8683	8684	8685	8696
Calcitek: Taperlock (Ex Hex)		8681		8682	8683	8684	8685	8696
Calcitek: Threadlock 4.0mmD (Ex Hex)		8681		8682	8683	8684	8685	8696
Eztetic 3.1mmD			2218	2219	2220	2221	2222	2223
Paragon: Bio-Vent 3.5mmD (3.5mmD Platform)	8661		8662	8663 2.5mm	8664 3.5mm	8665 4.5mm	8625 5.5mm	8626 6.5mm
Paragon: Bio-Vent 4.5mmD (4.5mmD Platform)	8671		8672	8673	8674	8675	8627	8628

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.0MM	0.73MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
ZIMMER® / CENTERPULSE® / CALCITEK® / PARAGON® CONT'D								
Paragon: Core-Vent 4.3mmD (3.5mmD Platform)	8661		8662	8663 2.5mm	8664 3.5mm	8665 4.5mm	8625 5.5mm	8626 6.5mm
Paragon: Core-Vent 5.3mmD (4.5mmD Platform)	8671		8672	8673	8674	8675	8627	8628
Paragon: Micro-Vent 3.25mmD (3.5mmD Platform)	8661		8662	8663 2.5mm	8664 3.5mm	8665 4.5mm	8625 5.5mm	8626 6.5mm
Paragon: Micro-Vent 4.25mmD (4.5mmD Platform)	8671		8672	8673	8674	8675	8627	8628
Paragon: Swiss Plus 3.8mmD			8935	8936	8937	8938		
Paragon: Swiss Plus 4.8mmD		8950		8951	8952	8953	8954	8955

44


LOCATOR® Abutments

Avoid Challenges When Processing Overdentures and Overdenture Attachments




For predictable processing into full and partial overdentures

- » Simple
- » Versatile
- » Patient-Friendly




Ideal for blocking our undercuts

- » Easy Handling
- » Multiple Uses



LOCATOR® R-Tx: Corrects up to 60° of implant divergence, provides an easy self-aligning path of insertion, esthetic pink housings, and new Duratec abutment coating for longer life and greater wear resistance. 2.85mm vertical, 5.55mm Ø



LOCATOR® R-Tx	
Item #	Product
30001-01	R-Tx Zero Retention Insert, Grey (4-Pack)
30002-01	R-Tx Low Retention Insert, Blue (4-Pack)
30003-01	R-Tx Medium Retention Insert, Pink (4-Pack)
30004-01	R-Tx High Retention Insert, White (4-Pack)
30027-01	R-Tx Low Retention, Limited Range, Blue (4-Pack)
30028-01	R-Tx Medium Retention, Limited Range, Pink (4-Pack)
30029-01	R-Tx High Retention, Limited Range, Yellow (4-Pack)
30030-01	R-Tx Attachment System, Mixed Retention Ltd Range Insert (2-Pack)
30012-01	R-Tx Processing Insert (4-Pack)
30013-01	R-Tx Denture Attachment Processing Assembly (4-Pack)
30014-01	R-Tx 3.35mm Abutment Analog (4-Pack) <i>R-Tx analogs for use with narrow diameter implants</i>
30015-01	R-Tx 4mm Abutment Analog (4-Pack) <i>R-Tx analogs for use with regular platform implants</i>
30016-01	R-Tx 5mm Abutment Analog (4-Pack) <i>R-Tx analogs for use with wide diameter implants</i>
30017-01	R-Tx Impression Coping (4-Pack)
30018-01	R-Tx Processing Spacer (4-Pack)
30021-01	R-Tx Retention Insert Tool
4008514	R-Tx Block Out Spacer (20-Pack)
4009665	R-Tx Patient Education Model

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.5MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
AB DENTAL®							
3.75mmD (NP)	30200-00	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mmD (WP)	30201-00	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
ACE®							
Bone Level NP		30900-01	30900-02	30900-03	30900-04	30900-05	30900-06
Bone Level RP		30901-01	30901-02	30901-03	30901-04	30901-05	30901-06
Tissue Level RP		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
3.5mm NP		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
3.5mm NP Tri-Channel	30500-00 0.25mm	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
3.5mm NP Internal Connection	30200-00	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.1mm RP		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
4.3mm RP Tri-Channel	30501-00	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
4.5mm RP Internal Connection	30201-00	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
5.0mm WP Tri-Channel		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06
BIOHORIZONS®							
Internal 3.0mmD		30600-01	30600-02	30600-03	30600-04	30600-05	30600-06
Internal 3.5mm/3.8mmD	30200-00	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
Internal 4.0mm/4.5mm/4.6mmD	30201-00	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
Internal 5.7mm/5.8mmD		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06
BIO-LOK®							
4.0mm (Ex Hex & Internal Connection)		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
BIOMET 3I®							
Certain 3.4mmD		30300-01	30300-02	30300-03	30300-04	30300-05	30300-06
Certain 4.1mmD		30301-01	30301-02	30301-03	30301-04	30301-05	30301-06
MicroMiniPlant 3.4mmD (Ex Hex)		30305-01	30305-02	30305-03	30305-04	30305-05	30305-06
MiniPlant 3.75mmD (Ex Hex)		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
Cylindrical 4.0mmD (Ex Hex)		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
BLUE SKY BIO®							
MicroMiniPlant 3.4mmD (Ex Hex)		30305-01	30305-02	30305-03	30305-04	30305-05	30305-06
BTI – Biotechnology Institute SL®		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
Cylindrical 4.0mmD (Ex Hex)		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
One Stage Regular		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
One Stage WP		30907-01	30907-02	30907-03	30907-04	30907-05	30907-06
BTI – BIOTECHNOLOGY INSTITUTE SL®							
External Hex 4.1mmD		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.5MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
DENTSPLY® / ASTRA®							
Ankylos C/X (Conical Connection)		30715-01	30715-02	30715-03	30715-04	30715-05	30715-06
EV 3.6mmD		30700-01	30700-02	30700-03	30700-04	30700-05	30700-06
EV 4.2mmD		30701-01	30701-02	30701-03	30701-04	30701-05	30701-06
EV 4.8mmD		30702-01	30702-02	30702-03	30702-04	30702-05	30702-06
TX 3.5mm/4.0mmD		30706-01	30706-02	30706-03	30706-04	30706-05	30706-06
4.5mm/5.0mmD		30707-01	30707-02	30707-03	30707-04	30707-05	30707-06
Frialit-2 / Xive 3.4mmD		30711-01	30711-02	30711-03	30711-04	30711-05	30711-06
Frialit-2 / Xive 3.8mmD		30712-01	30712-02	30712-03	30712-04	30712-05	30712-06
Frialit-2 / Xive 4.5mmD		30713-01	30713-02	30713-03	30713-04	30713-05	30713-06
Frialit-2 / Xive 5.5mmD		30714-01	30714-02	30714-03	30714-04	30714-05	30714-06
IMTEC®							
Self Tapping 3.75mmD (Ex Hex)		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
Universal 3.4mm/4.0mmD (Ex Hex)		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
IMPLANT DIRECT®							
Legacy 3.5mmD (Internal Hex)	30200-00	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
Legacy 4.5mmD (Internal Hex)	30201-00	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
Legacy 5.7mmD (Internal Hex)		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06
RePlus / RePlant / ReActive 3.5mmD (Trilobe NP)	30500-00	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
RePlus / RePlant / ReActive 4.3mmD (Trilobe RP)	30501-00	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
RePlus / RePlant / ReActive 5.0mmD (Trilobe WP)		30502-01	30502-02	30502-03	30502-04	30502-05	
SwishTapered / SwishPlus 4.8mmD		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
IMZ®							
3.3mmD (Ex Hex)		30306-01 1.2mm	30306-02	30306-03	30306-04	30306-05	30306-06
4.0mmD (Ex Hex)		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
INTRA-LOCK SYSTEM INTERNATIONAL®							
Threaded 3.75mm/4.0mmD (Ex Hex)		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
KEYSTONE / LIFECORE BIOMEDICAL®							
3.3mm/4.0mmD (Ex Hex)		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
5.0/5.5/6.0mmD (Ex Hex)		31202-01	31202-02	31202-03	31202-04	31202-05	31202-06
Stage-1 Regular		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
Restore RD 3.75mm/4.0mmD (Ex Hex)		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
MEGAGEN® / INTEGRATED DENTAL SYSTEMS®							
AnyRidge 3.5mmD (Internal Connection)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.5MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
MIS - MEDICAL IMPLANT SYSTEM®							
Internal Hex 3.75mm/4.2mmD (SP)	30200-00	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
Internal Hex 5.0mmD (WP)	30201-00	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
External SP (Ex Hex)		30306-01	30306-02	30306-03	30306-04	30306-05	30306-06
Mistral Internal Octagon (Stage One)		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
NEODENT®							
Cone Morse		32301-01	32301-02	32301-03	32301-04	32301-05	32301-06
Grand Morse		32700-01	32700-02	32700-03	32700-04	32700-05	32700-06
NOBEL BIO CARE® / BRANEMARK® / STERIOSS®							
Active / Conical Connection NP		30506-01	30506-02	30506-03	30506-04	30506-05	30506-06
Active / Conical Connection RP		30507-01	30507-02	30507-03	30507-04	30507-05	30507-06
Active / Conical Connection WP		30508-01	30508-02	30508-03	30508-04	30508-05	30508-06
Replace Select NP (3.5mmD, Trilobe)	30500-00	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
Replace Select RP (4.3mmD, Trilobe)	30501-00	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
Replace Select WP (5.0mmD, Trilobe)		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06
Branemark NP (3.3mmD, Ex Hex)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
Branemark RP (3.75mm/4.0mm/5.0mmD, Ex Hex)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Branemark WP (5.0mmD, Ex Hex)		30512-01	30512-02	30512-03	30512-04	30512-05	30512-06
Nobel Speedy NP (3.5mmD, Ex Hex)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
Nobel Speedy RP (4.3mmD, Ex Hex)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Replace Ex Hex NP (3.5mmD)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
Replace Ex Hex RP (4.3mmD)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Replace Ex Hex WP (5.0mmD)		30512-01	30512-02	30512-03	30512-04	30512-05	30512-06
Hex-Lock (HL) 3.8mm/4.5mmD (Ex Hex)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
OCO BIOMEDICAL®							
3.0mmD		31500-01 0.75mm	31500-02	31500-03	31500-04	31500-05	31500-06
3.5mmD	30200-00	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mmD	30201-00	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
OSSTEM® / HIOSSEN USA®							
SS Internal / HS Regular 4.8mmD		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
SS Internal Octagon RP		31606-01 1.5mm	31606-02	31606-03	31606-04	31606-05	31606-06
US Mini 3.3mm / 3.5mmD (Ex Hex)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
US Regular/ HU Regular 3.75mm / 4.0mmD (Ex Hex)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
ET Mini		31600-01	31600-02	31600-03	31600-04	31600-05	31600-06
ET Regular		31601-01	31601-02	31601-03	31601-04	31601-05	31601-06
SSII/III 4.5/4.8		31606-01	31606-02	31606-03	31606-04	31606-05	31606-06
OSTEO-IMPLANT CORPORATION®							
3.25mm SD, 3.75mm ST, 4.0mm DT		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06

IMPLANT ABUTMENTS/ATTACHMENTS

ABUTMENT TISSUE CUFF HEIGHT

IMPLANT ABUTMENT DESCRIPTION	0.5MM	1.0MM	2.0MM	3.0MM	4.0MM	5.0MM	6.0MM
SOUTHERN IMPLANT USA®							
Tri-Nex 3.5mmD	30500-00	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
Tri-Nex 4.3mmD	30501-00	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
Tri-Nex 5.0mmD		30502-01	30502-02	30502-03	30502-04	30502-05	
IT 4.8mmD (Internal Octagon)		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
STERNGOLD IMPLAMED®							
Hex Cylinder 3.3mm/4.0mmD		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Hex Screw 3.75mm/4.0mmD		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
STRAUMANN®							
Tissue Level 3.5mmD (NNC)		30405-01	30405-02	30405-03	30405-04	30405-05	30405-06
Tissue Level 4.8mmD (RN)		30406-01	30406-02	30406-03	30406-04	30406-05	30406-06
Tissue Level 6.5mmD (WN)		30407-01 1.5mm	30407-02	30407-03	30407-04	30407-05	30407-06
Bone Level 3.3mmD (NC)		30401-01	30401-02	30401-03	30401-04	30401-05	30401-06
Bone Level 4.1mmD (RC)		30402-01	30402-02	30402-03	30402-04	30402-05	30402-06
BLX Regular Base (RB) & Wide Base (WB)		30410-01	30410-02	30410-03	30410-04	30410-05	30410-06
SYBRON® / INNOVA®							
Endopore 4.1mmD (0.7mm Ex Hex)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Endopore 4.1mm/5.0mmD (One Stage)		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
Entegra 3.25mm/4.0mmD (0.7mm Ex Hex)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Pro Series XRT Bone Level 4.8mmD		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
Pro Series Tissue Level (TL) 4.1mm/4.8mmD		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
Hex Plus 4.1mmD		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
ZIMMER® / CENTERPULSE® / CALCITEK® / PARAGON®							
Screw-Vent/ Tapered Screw Vent 3.3mm/3.7mm/4.1mm Diameter (3.5mm Dia.)	30200-00	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
Screw Vent / Tapered Screw Vent 4.7mm Diameter (4.5mm Platform Dia.)	30201-00	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
Screw-Vent/ Tapered Screw Vent 6.0mm Diameter (5.7mm Platform Dia.)		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06
Calcitek: Spline 4.1mm (Ex Hex)		30207-01	30207-02	30207-03	30207-04	30207-05	30207-06
Calcitek: Swede-Vent 4.0mm Diameter (Ex Hex)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Calcitek: Swede-Vent TL (External Hex)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Calcitek: Taperlock (Ex Hex)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Calcitek: Threadlock 4.0mm (Ex Hex)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Paragon: Bio-Vent 3.5mm Diameter (3.5mm Platform Dia.)	30200-00	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
Paragon: Bio-Vent 4.5mm Diameter (4.5mm Platform Dia.)	30201-00	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
Paragon: Core-Vent 4.3mm Diameter (3.5mm Platform Dia.)	30200-00	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
Paragon: Core-Vent 5.3mm Diameter (4.5mm Platform Dia.)	30201-00	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
Paragon: Micro-Vent 3.25mm Diameter (3.5mm Platform Dia.)	30200-00	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
Paragon: Micro-Vent 4.25mm Diameter (4.5mm Platform Dia.)	30201-00	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
Paragon: Swiss Plus 4.8mm Diameter		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
Internal 3.5mm/3.8mmD	30200-00	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
Internal 4.0mm/4.5mm/4.6mmD	30201-00	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
Internal 5.7mm/5.8mmD		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

LOCATOR® Overdenture Implant (LODI): designed to provide a specialty overdenture implant /attachment option where narrow ridges, minimized surgical invasiveness and financial considerations play a role in the implant choice and where a superior overdenture attachment can be utilized with no compromise. The best of both worlds—a narrow diameter implant utilizing the proven LOCATOR® Attachment!

50

LOCATOR Overdenture Implant System



4007351



4007352



4007564



4007561



4007450



4007451





4007452



4007455

LODI NDI Implant System

Item #	Product
4007351	LODI Attachment, 2.5mm
4007352	LODI Attachment, 4.0mm
4007318	LODI Attachment, 3.0mm
4007319	LODI Attachment, 5.0mm
4007320	LODI Attachment, 6.0mm
4007365	LODI Direction Indicator
4007538	LODI Drill 1.6mm
4007539	LODI Drill 2.1mm
4007541	LODI Drill 2.6mm
4007374	LODI Drill Extender
4007540	LODI Drill 2.4mm
4007560	LODI Surgical Tray
4007561	LODI Standard Surgical Kit
4007564	LODI Premium Surgical Kit
4007450	LODI Implant 2.4mm x 10mm, 2.5mm Cuff
4007451	LODI Implant 2.4mm x 12mm, 2.5mm Cuff
4007452	LODI Implant 2.4mm x 14mm, 2.5mm Cuff



LODI NDI Implant System

Item #	Product
4007455	LODI Implant 2.4mm x 10mm, 4.0mm Cuff
4007456	LODI Implant 2.4mm x 12mm, 4.0mm Cuff
4007457	LODI Implant 2.4mm x 14mm, 4.0mm Cuff
4007460	LODI Implant 2.9mm x 10mm, 2.5mm Cuff
4007461	LODI Implant 2.9mm x 12mm, 2.5mm Cuff
4007462	LODI Implant 2.9mm x 14mm, 2.5mm Cuff
4007465	LODI Implant 2.9mm x 10mm, 4.0mm Cuff
4007466	LODI Implant 2.9mm x 12mm, 4.0mm Cuff
4007467	LODI Implant 2.9mm x 14mm, 4.0mm Cuff
4007470	LODI Implant 3.5mm x 10mm, 2.5mm Cuff
4007477	LODI Implant 3.5mm x 14mm, 4.0mm Cuff
4007476	LODI Implant 3.5mm X 12mm, 4.0mm Cuff
4007475	LODI Implant 3.5mm x 10mm, 4.0mm Cuff
4007471	LODI Implant 3.5mm x 12mm, 2.5mm Cuff
4007472	LODI Implant 3.5mm x 14mm, 2.5mm Cuff
4007362	LODI Torque Indicating Ratchet Wrench
4007357	LODI Implant Latch Driver

Call for additional implant diameter options!

ATTACHMENTS

IMPLANT ABUTMENTS/ATTACHMENTS

OT Equator is the smallest dimensional overdenture attachment, in terms of size, available on the market today. The maximum vertical dimension (attachment, retentive cap and housing) is only 2.1mm with a diameter (attachment, retentive cap and housing) of 4.5mm.

The attachment system includes retentive caps in various colors corresponding to different levels of retention as well as the Smartbox titanium housing for retentive caps. The innovative Smartbox design allows a passive insertion of the cap even under extreme divergence conditions up to 50° thanks to a tilting mechanism with a pivoting fulcrum. The Smartbox is compatible with the OT Equator direct attachment.

Call for information on OT Equator attachments which are available with heights ranging from 0.5mm to 7.0mm for a variety of implant systems and may even be custom-made if needed.

52






OT Equator® & OT Equator® Smartbox



OT Equator® & OT Equator® Smartbox

Item #	Product
2239GSF	Equator 1.6mm Threaded Sleeve
2192ECE	Equator Cap Assortment Kit
244CAIN	Equator Impression Copings (2-Pack)
20144AE	Equator Laboratory Analogs (2-Pack)
20760CE	Equator Latch-Type Driver
2140CEN	Equator Processing Caps - Black (4-Pack)
2140CET	Equator Retentive Caps - Clear, Standard (4-Pack)
2140CER	Equator Retentive Caps - Pink, Soft (4-Pack)
2140CEV	Equator Retentive Caps - Violet, Strong (4-Pack)
2140CEG	Equator Retentive Caps - Yellow, Extra Soft (4-Pack)
2335SBC	Equator SmartBox
2330SBE	Equator SmartBox Housing
2335CSB	Equator SmartBox Processing Caps - Black (4-Pack)
2141CAE	Equator SS Housings (2-Pack)
2141CTE	Equator Ti Housings (2-Pack)

Sphero Flex is a self-parallel ball attachment correcting up to 15° of divergence for use with 2.5mm ball abutments. Sphero Block is a fixed sphere available in either 1.8 or 2.5mm ball diameters.




Sphero Flex and Block Ball Attachments

Item #	Product
0800161	2.5mm Retention Caps—Light Green (6-Pack) <i>For worn spheres, Rigid Retention (2.2mm)</i>
0800162	2.5mm Retention Caps—Pink (6-Pack) <i>Standard- 1.65lbs</i>
0800163	2.5mm Retention Caps—Green (6-Pack) <i>Soft, gummy for patient comfort</i>
0800164	2.5mm Retention Caps—Yellow (6-Pack) <i>Soft- 1.21lbs</i>
0800166	2.5mm Retention Caps—White (6-Pack) <i>Strong- 2.87lbs</i>
0800167	2.5mm Retention Caps—Orange (6-Pack) <i>For worn spheres, Elastic Retention (2.2mm)</i>
0800194	Cap Insertion Tool 2.5mm/1.8mm
0800197	Sphero 2.5mm Analog (2-Pack)
0800180	2.5mm Stainless Cap Housing (2-Pack)
0800280	1.8mm Stainless Cap Housing (2-Pack)
0800720	Sphero 1.8mm Analog (2-Pack)
0800264	1.8mm Retention Caps—Yellow (6-Pack) - Soft - 1.00lbs
0800262	1.8mm Retention Caps—Pink (6-Pack) - Standard - 1.10lbs
0800266	1.8mm Retention Caps—White (6-Pack) - Strong - 2.20 lbs
0800263	1.8mm Retention Caps—Green (6-Pack) <i>Soft, gummy for patient comfort</i>
0800169	2.5mm Retention Caps—Black (6-Pack) Processing Cap
0800269	1.8mm Retention Caps—Black (6-Pack) Processing Cap
0800267	1.6mm Retention Caps—Orange (6-Pack) <i>For worn spheres, Elastic Retention (1.6mm)</i>
0800261	1.6mm Retention Caps—Light Green (6-Pack) <i>For worn spheres, Rigid Retention (1.6mm)</i>
0800250	Sphero Directional Rings

ATTACHMENTS

REMOVABLE PARTIAL DENTURE ATTACHMENTS


The demand for partial dentures is expected to grow at 6.8% year to year over the next decade. Partner with the company that can help you provide the best patient specific solutions.



AMERICAN MINI

*Kennedy Class I
semi-precision attachment
Mini-Dalbo/Swiss Mini*

- » 3mm tall, male 2.5mm Ø,
female 3.5mm Ø



IC ATTACHMENT

*Kennedy Class III
Hannes Anchor*

- » male: 5.2mm length,
2mm Ø, 1mm nose
- » female: 2mm Ø,
0.8mm deep



CEKA REVAX

*Kennedy Class I Precision
Resilient attachment*


- » M2 (small): 3.80mm tall,
3.40mm Ø
- » M3 (standard): 4.35mm
tall, 4.40mm Ø



VERTIX

Kennedy Class IV

- » 6mm (may be reduced
to 3mm) tall, 3.3mm Ø



SAGIX

*Ideal for Kennedy Class I
VKS-SG, OT Strategy, Allegro*


- » 1.7mm ball: 3.10mm tall,
3.10mm Ø
- » 2.2mm ball: 4.10mm tall,
4.20mm Ø



O-RING DISTAL EXTENSION

Kennedy Class I

- » 5.2mm (may be reduced
to 2.5mm) tall, 2.00mm Ø



MAYS UNILATERAL

Ideal for Kennedy Class II


- » 3mm tall, 2mm Ø



SPHERO RETROFIT

- » replacing worn extracoronal
attachments, such as the
ERA or Ceka
- » 1.5mm tall ball, 1.8mm Ø

Sphero Retrofit Kits: designed to retrofit worn ERA, Ceka, and other worn extracoronaral and overdenture attachments
Use Ceka Site or metal-to-metal bonding composite. retrofit: 2.7mm* vertical, 4.3mm Ø 2 TiN coated spheres, 1 mandrel/ key, 2 pink females, 2 protector disks. **metal housings 0800280 available separately**






Sphero Retrofit Kit

Item #	Product
0800140	Sphero Retrofit Kit

Sphero Retrofit Kit

55

Sagix Attachment: Sagital 1.7 and 2.2mm Ball attachment with female that engages the full sphere for maximum retention. Unique female engages entire sphere for more retention and less attachment wear Compatible with Bredent VKS-SG 1.7 and 2.2mm balls 1.7mm: 3.1mm vertical, 3.1mm Ø 2.2mm: 4.1mm vertical, 4.2mm Ø



Sagix

Item #	Product
0101305	Sagix 1.7 Complete Plastic w/ housing 1 plastic male, 1 each yellow, white and red females, 1 processing pink female, 1 metal housing
0101356	Sagix 1.7 Complete Metal w/ housing 1 metal male, 1 each yellow, white and red females, 1 processing pink female, 1 metal housing
0101357	Sagix 1.7 Complete Metal 1 metal male, 1 each yellow, white and red females, 1 processing pink female
0101325-2	Sagix 1.7 Metal Housing (2-Pack) CoCr
0101331-2	Sagix 1.7 Model Analog (2-Pack)
0101316-6	Sagix 1.7 Female Dupl Dummy (6-Pack)
0101316-50	Sagix 1.7 Female Dupl Dummy (50-Pack)
0101312-6	Sagix 1.7 White Female- Reduced (6-Pack) white = 1.5lbs retention
0101312-50	Sagix 1.7 White Female- Reduced (50-Pack)

Sagix

ATTACHMENTS

EXTRACORONAL ATTACHMENTS

Sagix

56



0101311-6



0101313-6



0101314



0101300P




0101323



0101358



0101307



0101375-2



0101362-6




0101274

Sagix (cont'd)

Item # Product


0101311-6	Sagix 1.7 Yellow Female- Standard (6-Pack) <i>yellow = 2.3lbs retention</i>
0101311-50	Sagix 1.7 Yellow Female- Standard (50-Pack)
0101313-6	Sagix 1.7 Red Female- Increased (6-Pack) <i>red = 3.1lbs retention</i>
0101313-50	Sagix 1.7 Red Female- Increased (50-Pack)
0101314	Sagix 1.7 Female Combo pack <i>2 white, 2 yellow, 2 red 1.7mm females</i>
0101300P	Sagix Universal Parallel Mandrel
0101323	Sagix Universal Insertion Tool
0101307	Sagix 2.2 Complete Plastic w/ housing <i>1 plastic male, 1 each yellow, white and red females, 1 processing pink female, 1 metal housing.</i>
0101358	Sagix 2.2 Complete Metal w/ housing <i>1 metal male, 1 each yellow, white and red females, 1 processing pink female, 1 metal housing</i>
0101274	Universal Threaded Ball 2.2mm (M2 x .040)
0101375-2	Sagix 2.2 Metal Housing (2-Pack) CoCr
0101381-2	Sagix 2.2 Model Analog (2-Pack)
0101366-6	Sagix 2.2 Female Dupl Dummy (6-Pack)
0101366-50	Sagix 2.2 Female Dupl Dummy (50-Pack)
0101362-6	Sagix 2.2 White Female- Reduced (6-Pack) <i>white = 1.7 lbs retention</i>
0101362-50	Sagix 2.2 White Female- Reduced (50-Pack)



Sagix (cont'd)

Item #	Product
0101361-6	Sagix 2.2 Yellow Female- Standard (6-Pack) <i>yellow = 2.5 lbs retention</i>
0101361-50	Sagix 2.2 Yellow Female- Standard (50-Pack)
0101363-6	Sagix 2.2 Red Female- Increased (6-Pack) <i>red = 3.3 lbs retention</i>
0101363-50	Sagix 2.2 Red Female- Increased (50-Pack)
0101364	Sagix 2.2 Female Combo pack


Vertix® Attachment: extracoronral slide attachment providing limited resilience and excellent lateral stability and esthetics with castable housing: 3.0mm vertical, 3.3mm Ø with metal housing: 3.0mm vertical, 4.1mm Ø



Vertix

Item #	Product
0101799	Vertix P Combo Kit (2-Pack) <i>2 castable P males, 2 castable housings, 1 insertion tool, 1 paralleling mandrel, 6 female riders (2 white, 2 yellow, 2 red)</i>

American Mini®: (Mini Dalbo): extracoronral attachment allowing hinge movements with excellent stability mini: 3.0mm vertical, male 2.5mm Ø, female 3.5mm Ø



American Mini (Mini Dalbo)

Item #	Product
1800000	American Mini Complete <i>Stainless Female and Delrin Male</i>

ATTACHMENTS

EXTRACORONAL ATTACHMENTS

O-Ring Distal Extension: provides free movement in all planes for maximum abutment protection. A proven retentive system utilizing an easily replaceable rubber O ring. Resilient attachment with vertical and hinge stress-breaking action for free-end and bounded partials.

OSO Distal Extension



O-Ring



1500150



1500153


OSO Distal Extension

Item #	Product
1500150	O-Ring Distal Extension Attachment #3 <i>Burn-out plastic pattern, 3 white o-rings, 1 red o-ring</i>
1500153	O-Ring Distal Ext. Paralleling Mandrel #3
1500502	O-Ring Metal Retainer Standard #3

58

Mays® Unilateral: designed specifically for the unilateral distal extension. Does not require a lingual or palatal arm. Metal framework is not required; parallelism is not necessary, even on bilateral cases. Requires: 3mm tall, 2mm Ø

Mays® Unilateral



Mays


Mays® Unilateral

Item #	Product
0200060	Mays Unilateral Attach UR/LL

IC Attachment: a popular spring loaded plunger attachment that provides resilient movement for abutment protection without requiring an abutment crown. The IC attachment requires a 180 degree reciprocal lingual arm if free end.
Universal Male: 5.2mm length, 2mm Ø, 1mm nose Universal Female: 2mm Ø, 0.8mm deep



IC



IC

Universal IC Attachment

Item #	Product
1800000	Universal IC Attachment

Universal IC Attachment

Ceka® and Ceka Revax®: adjustable precision attachment allowing vertical, hinge, & rotational movements M2 Revax: 3.8mm vertical, 3.4mmØ M3 Revax: 4.35mm vertical, 4.4mmØ M3 Classic: 3.55mm vertical, 4.3mmØ M2 female attachments only




Ceka




2RP0075



2RP0061



2RQ0032



2RQ0031

Ceka® and Ceka Revax


Item #	Product
M2 Retention Caps	
2RP0061	M2 Pallax Solder Cap (2-Pack)
2RP0075	M2 Acrylic Retention Caps (2-Pack)
2RP0095	M2 Axial Acrylic Retention Cap (2-Pack)
M2 Male Spring Pins	
2RQ0031	M2 Spring Pins Standard (2-Pack)
2RQ0032	M2 Spring Pins 195 (2-Pack)
2RQ0033	M2 Spring Pins 202 (2-Pack)
2RP0096	M2 Space Maintainer

Ceka® and Ceka Revax


EXTRACORONAL ATTACHMENTS



2RT00A1



3X694C2



3694C19

Ceka® and Ceka Revax

Item #	Product
M2 Attachment Tools	
2RT00A1	A1 Activating Tool & Driver
M3 Male Spring Pins	
3X694C2	M3 694C Standard Spring Pin (2-Pack)
3694C19	M3 694C Spring Pin 195 (2-Pack)
3694C20	M3 694C Spring Pin 202 (2-Pack)
3694C29	M3 694C Spring Pin 209 (2-Pack)
3694CL8	Ceka Sizer
3TL03H1	M3 Working Dummy Pin (2-Pack)
3X694B0	M3 Space Maintainer 694B
M3 Retention Caps	
3694KS2	M3 Acrylic Retention Cap (2-Pack) 694 AKS
3X694A2	M3 Solder Retention Cap (2-Pack) 694 AS
3X6AKS2	M3 694 AKS2 Acrylic Retention (2-Pack)

Splintastik®: a very simple plastic male and female rod and tube for use in anterior segmented bridgework. The female features a smooth ceramic insert to assure accurate castings. Small 0.8mm diameter Large 1.2mm diameter.



0200001



0200003


Splintastik

Item #	Product
0200001	Splintastik Starter Kit 2 each: Large Male, Female, and ceramic insert 2 each: Small Male, Female, and ceramic insert 1 each: Small and Large Mandrels, and Precision Reamers
0200003	Splintastik Refill LG (6-Pack) Large Male, Female, and Ceramic insert, size 1.2mm Ø


Splintastik

61


Better Plastic Dovetail (BPD): a precision stress breaker ideal for posterior bridges as well as Intracoronal semi-precision slide attachment. Rigid attachment (always use a retentive clasp arm) with friction (not adjustable) retention. Ideal for nonparallel abutments and segmented units. Inside of female rounded to facilitate devesting and avoid attachment wear. Built in paralleling mandrels for both male and female.



BPD



BPD Female



BPD Male

Double Mandrel BPD

Item #	Product
0200040-6	BPD (Better Plastic Dovetail)(6-Pack)
0200040-12	BPD (Better Plastic Dovetail)(12-Pack)
0200040-50	BPD (Better Plastic Dovetail)(50-Pack)
0200040-100	BPD (Better Plastic Dovetail)(100-Pack)

Double Mandrel BPD

TRILOR®

HIGH PERFORMANCE POLYMER

TriLor is a millable fiber-reinforced resin material designed for suprastructures for full arch restorations and bridge frameworks.

This super-strong resin flexes and bends under stress, transferring less energy to the implants, making it an ideal support material for implant restorations.

- Lightweight and Esthetic
- Full Covalent Bond
- Time and Cost Savings



**Fabricated at the bench in under
30 minutes with a handpiece!**

THE CHOICE IS YOURS

DIGITAL

Trilor 98 and 95
CAD/CAM Pucks



ANALOG

3 Different
Thicknesses



See What our Customers Say


"Predictable outcome is predictable income. I use TriLor as a tool to make my implant retained cases predictable."

– Jeremy Wohlers, CDT
Owner of Esthetics in Function Dental Lab

"TriLor has improved workflow in my laboratory and chairside while still achieving a very esthetic outcome"

– Eugene Royzengurt
Owner of Apple Dental Lab

TriLor® is a millable fiber-reinforced resin material designed for suprastructures for full arch restorations and bridge frameworks. This super-strong resin flexes and bends under stress, transferring less energy to the implants, making it an ideal support material for implant restorations. The TriLor® Arch uniquely allows the Dental Professional to design the framework using a lathe or handpiece. CAD / CAM TriLor® Pucks also available.




Fiberglass Technology

Item #	Product
50000-35	TriLor Arch, 3.5mm (3-Pack)
50000-55	TriLor Arch, 5.5mm (3-Pack)
50000-75	TriLor Arch, 7.5mm (3-Pack)
50000-16	TriLor 98 CAD Puck, 16mm
50000-16P	TriLor 98 CAD Puck Pink, 16mm
50000-20	TriLor 98 CAD / CAM Puck, 20mm
50000-20P	TriLor 98 CAD / CAM Puck Pink, 20mm
50000-25	TriLor 98 CAD / CAM Puck, 25mm
50005-16	TriLor 95 CAD / CAM Puck, 16mm
50005-20	TriLor 95 CAD / CAM Puck, 20mm
50005-25	TriLor 95 CAD / CAM Puck, 25mm

TRILOR® INDICAITONS FOR SUPRASTRUCTURE DESIGN

- A. Minimum must be 7mm² (3.5mm x 2mm)
- B. Minimum thickness around cylinders is 0.8mm
- C. Maximum cantilever is 10mm



TRILOR® MATERIAL COMPARISON

Properties	TriLor®	Trinia™	PEKKTON®	JUVORA® PEEK	Zirconia + Ittrio	CrCo
Tensile Strength	380 MPa	169 MPa	115 MPa	100 MPa	348 MPa	695 MPa
Flexural Strength	540 MPa	393 MPa	200 MPa	170 MPa	1200 MPa	600 MPa
Tensile Elongation	2%	2.7%	4.4%	20%	-	8%
Flexural Modulus	26 GPa	18.8 GPa	5.0 GPa	4.0 GPa	210 GPa	175 GPa
Compressive Strength	530 MPa	347 MPa	246 MPa	118 MPa	2000 MPa	448 MPa

ELEVATE YOUR RESTORATIONS

With G-CAM

AVAILABLE IN 9 SHADES & MULTIPLE THICKNESSES

Compatible with universal &
Zirkonzahn® CAD/CAM systems



Experience the next generation of durability and aesthetics with our graphene-reinforced biopolymer disc, designed for superior resilience, patient comfort and long-term patient satisfaction.

Strength & Durability

Engineered for resilience with high fracture resistance, flexibility and masticatory load absorption.

Multichromatic Aesthetics



Multichromatic structure mimics natural translucency and color gradients for stunning restorations.

Patient Comfort First

Lightweight, low water absorption and gentle on opposing dentition for long-term wearability.

Versatile & Compatible

Ideal for temporary and permanent restorations – crowns to full arches.




G-CAM Multichroma 95mm

Item #	Product
GC00089	G-CAM Multichroma 95mm, 14mm x A1
GC00090	G-CAM Multichroma 95mm, 14mm x A2
GC00157	G-CAM Multichroma 95mm, 14mm x A3
GC00158	G-CAM Multichroma 95mm, 14mm x A3'5
GC00091	G-CAM Multichroma 95mm, 14mm x B1
GC00092	G-CAM Multichroma 95mm, 14mm x B2
GC00094	G-CAM Multichroma 95mm, 14mm x BL2
GC00095	G-CAM Multichroma 95mm, 14mm x C2
GC00096	G-CAM Multichroma 95mm, 16mm x A1
GC00097	G-CAM Multichroma 95mm, 16mm x A2
GC00098	G-CAM Multichroma 95mm, 16mm x A3
GC00099	G-CAM Multichroma 95mm, 16mm x A3'5
GC00100	G-CAM Multichroma 95mm, 16mm x B1
GC00101	G-CAM Multichroma 95mm, 16mm x B2
GC00103	G-CAM Multichroma 95mm, 16mm x BL2
GC00104	G-CAM Multichroma 95mm, 16mm x C2
GC00112	G-CAM Multichroma 95mm, 20mm x A1
GC00113	G-CAM Multichroma 95mm, 20mm x A2
GC00114	G-CAM Multichroma 95mm, 20mm x A3
GC00115	G-CAM Multichroma 95mm, 20mm x A3'5
GC00116	G-CAM Multichroma 95mm, 20mm x B1
GC00117	G-CAM Multichroma 95mm, 20mm x B2
GC00119	G-CAM Multichroma 95mm, 20mm x BL2
GC00120	G-CAM Multichroma 95mm, 20mm x C2
GC00130	G-CAM Multichroma 95mm, 24mm x A1
GC00131	G-CAM Multichroma 95mm, 24mm x A2
GC00132	G-CAM Multichroma 95mm, 24mm x A3
GC00133	G-CAM Multichroma 95mm, 24mm x A3'5
GC00134	G-CAM Multichroma 95mm, 24mm x B1
GC00135	G-CAM Multichroma 95mm, 24mm x B2
GC00137	G-CAM Multichroma 95mm, 24mm x BL2
GC00138	G-CAM Multichroma 95mm, 24mm x C2

ATTACHMENTS

GRAPHENE TECHNOLOGY

G-CAM Multichroma 98.5mm




G-CAM Multichroma 98.5mm

Item #	Product
GC00250	G-CAM Multichroma 98.5mm, 14mm x A1
GC00251	G-CAM Multichroma 98.5mm, 14mm x A2
GC00317	G-CAM Multichroma 98.5mm, 14mm x A3
GC00252	G-CAM Multichroma 98.5mm, 14mm x A3'5
GC00253	G-CAM Multichroma 98.5mm, 14mm x B1
GC00254	G-CAM Multichroma 98.5mm, 14mm x B2
GC00255	G-CAM Multichroma 98.5mm, 14mm x BL1
GC00256	G-CAM Multichroma 98.5mm, 14mm x BL2
GC00257	G-CAM Multichroma 98.5mm, 14mm x C2
GC00318	G-CAM Multichroma 98.5mm, 16mm x A1
GC00258	G-CAM Multichroma 98.5mm, 16mm x A2
GC00259	G-CAM Multichroma 98.5mm, 16mm x A3
GC00260	G-CAM Multichroma 98.5mm, 16mm x A3'5
GC00261	G-CAM Multichroma 98.5mm, 16mm x B1
GC00262	G-CAM Multichroma 98.5mm, 16mm x B2
GC00263	G-CAM Multichroma 98.5mm, 16mm x BL1
GC00264	G-CAM Multichroma 98.5mm, 16mm x BL2
GC00265	G-CAM Multichroma 98.5mm, 16mm x C2
GC00266	G-CAM Multichroma 98.5mm, 18mm x A1
GC00319	G-CAM Multichroma 98.5mm, 18mm x A2
GC00320	G-CAM Multichroma 98.5mm, 18mm x A3
GC00267	G-CAM Multichroma 98.5mm, 18mm x A3'5
GC00268	G-CAM Multichroma 98.5mm, 18mm x B1
GC00269	G-CAM Multichroma 98.5mm, 18mm x B2
GC00270	G-CAM Multichroma 98.5mm, 18mm x BL1
GC00271	G-CAM Multichroma 98.5mm, 18mm x BL2
GC00272	G-CAM Multichroma 98.5mm, 18mm x C2
GC00273	G-CAM Multichroma 98.5mm, 20mm x A1
GC00274	G-CAM Multichroma 98.5mm, 20mm x A2
GC00275	G-CAM Multichroma 98.5mm, 20mm x A3
GC00276	G-CAM Multichroma 98.5mm, 20mm x A3'5
GC00277	G-CAM Multichroma 98.5mm, 20mm x B1
GC00278	G-CAM Multichroma 98.5mm, 20mm x B2

G-CAM Multichroma 98.5mm	
Item #	Product
GC00279	G-CAM Multichroma 98.5mm, 20mm x BL1
GC00280	G-CAM Multichroma 98.5mm, 20mm x BL2
GC00281	G-CAM Multichroma 98.5mm, 20mm x C2
GC00282	G-CAM Multichroma 98.5mm, 22mm x A1
GC00283	G-CAM Multichroma 98.5mm, 22mm x A2
GC00284	G-CAM Multichroma 98.5mm, 22mm x A3
GC00285	G-CAM Multichroma 98.5mm, 22mm x A3'5
GC00286	G-CAM Multichroma 98.5mm, 22mm x B1
GC00287	G-CAM Multichroma 98.5mm, 22mm x B2
GC00288	G-CAM Multichroma 98.5mm, 22mm x BL1
GC00289	G-CAM Multichroma 98.5mm, 22mm x BL2
GC00290	G-CAM Multichroma 98.5mm, 22mm x C2
GC00291	G-CAM Multichroma 98.5mm, 24mm x A1
GC00292	G-CAM Multichroma 98.5mm, 24mm x A2
GC00293	G-CAM Multichroma 98.5mm, 24mm x A3
GC00294	G-CAM Multichroma 98.5mm, 24mm x A3'5
GC00295	G-CAM Multichroma 98.5mm, 24mm x B1
GC00296	G-CAM Multichroma 98.5mm, 24mm x B2
GC00297	G-CAM Multichroma 98.5mm, 24mm x BL1
GC00298	G-CAM Multichroma 98.5mm, 24mm x BL2

G-CAM Multichroma 98.5mm	
Item #	Product
GC00299	G-CAM Multichroma 98.5mm, 24mm x C2
GC00300	G-CAM Multichroma 98.5mm, 26mm x A1
GC00301	G-CAM Multichroma 98.5mm, 26mm x A2
GC00302	G-CAM Multichroma 98.5mm, 26mm x A3
GC00303	G-CAM Multichroma 98.5mm, 26mm x A3'5
GC00304	G-CAM Multichroma 98.5mm, 26mm x B1
GC00305	G-CAM Multichroma 98.5mm, 26mm x B2
GC00306	G-CAM Multichroma 98.5mm, 26mm x BL1
GC00307	G-CAM Multichroma 98.5mm, 26mm x BL2
GC00308	G-CAM Multichroma 98.5mm, 26mm x C2
GC00309	G-CAM Multichroma 98.5mm, 30mm x A1
GC00310	G-CAM Multichroma 98.5mm, 30mm x A2
GC00311	G-CAM Multichroma 98.5mm, 30mm x A3
GC00312	G-CAM Multichroma 98.5mm, 30mm x A3'5
GC00313	G-CAM Multichroma 98.5mm, 30mm x B1
GC00314	G-CAM Multichroma 98.5mm, 30mm x B2
GC00161	G-CAM Multichroma 98.5mm, 30mm x BL1
GC00315	G-CAM Multichroma 98.5mm, 30mm x BL2
GC00316	G-CAM Multichroma 98.5mm, 30mm x C2
1960	G-TOOLS Polishing Kit



CLINICAL SPECIALTIES

Chairside®: pink dual-cure Bis-Acryl chairside attachment pickup material, odorless and tasteless. Working time 2 minutes and 15 seconds, curing time 7 minutes (self cure). 4mL syringe enables the pickup of about 30 attachments




4009566



4009566



4009566



4009565



4009574

Clinical Specialties

Item

Product

4009566

CHAIRSIDE Pickup Material 4mL
Syringe and Plunger, 15 Mixing Tips, 15 Angled Tips

4009565

CHAIRSIDE Pickup Material 18mL
Syringe and Plunger, 10 Mixing Tips, 10 Angled Tips
Low viscosity provides excellent control during chairside pickup. Ideal for picking up Temporary Cylinders in hybrid dentures.

Optimal handling characteristics:

- Tissue-colored material bonds to itself.
- Dual cure – self or light cure.
- Will set in approximately 5 minutes. The working time is approximately 2 minutes and 15 seconds.
- The material will self cure in approximately 7 minutes. Curing time can be shortened by light curing.

4009125-15

CHAIRSIDE Short Mixing Tips (15-Pack)

4009549


Denture Removal Tool (10-Pack)

4009574

CHAIRSIDE Block Out Material-4 Syringes
(1.2mL), 20 Mixing Tips

A soft, light-cured resin material designed to protect a patient's tissue from undercuts during the chairside overdenture processing. The optimal handling and low viscosity ensure precise application and easy removal, while the contrasting blue pigment provides better visibility during clinical procedures

Perma Ret: a simple solution for Permanent Retention of acrylic and composite denture teeth in both full and partial dentures. Simply drill a small hole in the tooth with the provided twist drill. Cut a length of wire and insert in the wire holder. Thread the wire into the tooth. Bend the wire into desired position. Then process the denture base material. Wire is now anodized for improved esthetics.




Perma Ret



0109000




0109002



Perma Ret



Perma Ret



6000005

Perma Ret

Item #	Product
0109000	Perma Ret Kit 3 twist drills, 1 wire holder, and 13 pieces of 80mm long threaded silver wire
0109002	Perma Ret Wire Refill 13 pieces of 80mm long threaded silver wire
0109002-6	Perma Ret Wire Refill (6-Pack)
0109004	Perma Ret Twist Drills (3)
0109006	Perma Ret Tool
6000005	Denture Stippling Mandrel The 1/4 inch shaft will fit most standard chucks on a dental lathe. The user should place 4 of either a #2, #4, or #6 round straight handpiece or FG steel or carbide burs into the stippling mandrel and tighten securely with the Allen wrench that is included with the stippling mandrel. Best results are achieved using dull or used burs. *Burs are not included with the mandrel*

Perma Ret

69

Pro Blue: The finest single step polishing compound available. Water soluble for easy cleaning and a high shine on acrylic, composites, thermoplastics, and all alloys.



0601251

Pro Blue NTR Polish

Item #	Product
0601251	Pro Blue NTR Polishing Compound 150g

Pro Blue NTR Polish

ATTACHMENTS

CAD-CAM BAR SOLUTIONS




4008589



4008028



0101702




1776000



0101274



0101380



0101280

CAD-CAM Bar Solutions

Item #	Product
4008589-2	LOCATOR Bar Female threaded 2mm (2-Pack)
4008589-10	LOCATOR Bar Female threaded 2mm (10-Pack)
4008028-2	LOCATOR Bar Processing Package (2-Pack)
4008028-10	LOCATOR Bar Processing Package (10-Pack)
0101702	Milled Hader Bar Attachment
1776000	Universal Plunger Loc Attachment
0101274	Universal Threaded Ball 2.2mm <i>Based on the universally accepted thread: M2 x .040 Our existing drill (4009102) and tap (4009104) may be used with the Universal Threaded Ball.</i>
0101380	Milled Bar Sagix Complete Attachment
0101280	Milled Bar Clix Complete Attachment
2339SFE	OT Equator Threaded 1.6mm Abutment
239SFE2	OT Equator Threaded 2mm Abutment
2239GSF	OT Equator 1.6mm Threaded Sleeve

ALLOY CHART

PALLAX	STAINLESS STEEL
White - Au 2 - Ag 37 - Pt 9.5 - Pd 37 - Cu 12.5 - Co 2	White - Fe 66 - Cr 17 - Ni 11 -
Melting range: 1055-1130 °C (1931-2066 °F)	Mo 2 - other 4 - Do not heat!
IRAX - Precious	POLYPROPYLENE
White - Au 60 - Pt 24 - Pd 15 - Ir 1	Hader Clips
Melting range: 1400-1460 °C (2552-2660 °F)	
NOPRAX - Non-precious	POLYPROPYLENE
White - Au 2 - Ag 37 - Pt 9.5 - Pd 37 - Cu 12.5 - Co 2	Clix, Sagix Females
Melting range: 1055-1130 °C (1931-2066 °F)	
ORAX	THERMOPLASTIC ELASTOMER
Yellow - Au 67 - Pt 8.6 - Ag 13.2 - Cu 11 - Zn 0.Y 2	# 1201 B
Melting range: 910-995 °C (1670-1823°F)	
Heat treatment: 60 min at 400 °C (752 °F)	
PLATIRAX	TITANAX
White - Pt 85 - Ir 15	White - Ti 90 - Al 6 - V 4 - Do not heat!
Melting range: 1790-1820 °C (3254-3308°F)	

Restoring Confidence.

PREAT ABUTMENTS

will meet and/or exceed fit and function of the OEM abutments*

- » OEM Implant Warranty Matched
- » Manufactured in the US
- » Lifecycle Product and Technical Support



SAMPLE TESTING RESULT

MATING FEATURE 1

Screw head clearance to the abutment

MATING FEATURE 2

Screw shaft clearance to the abutment

MATING FEATURE 3

Hexagonal connection clearance

MATING FEATURE 4


Abutment tapered seat

MATING FEATURE 5

Abutment hexagonal clearance to the implant


1

Big Implant Company #2 Mating Feature 1




2

Big Implant Company #3 Mating Feature 2




3

Big Implant Company #4 Mating Feature 3




4

Big Implant Company #5 Mating Feature 4



5

Big Implant Company #3 Mating Feature 5



*Test results by independant third party, 2018



IMPLANT PROSTHETICS

Preat **CHOICE LIBRARY**

We've expanded the Preat Choice Library to give you even more implant system choices. With these additions, the Preat Choice Library is compatible with all major implant systems – **40 platforms in total**.

For each system platform, the scan body gives you access to 18 different restorative options, including engaging and non-engaging, for fixture-level restorations.

» **ANGLED SCREW CHANNEL (ASC) TITANIUM BASES**

designed for optimal esthetics and function with 0° to 30° channel angulation, cement groove to prevent overflow, plus retention grooves and blasted surface for enhanced cement adhesion

» **STRAIGHT TITANIUM BASES** for digital or traditional workflow

» **VERIFICATION CYLINDERS** for provisional restorations with variable height options

» **PLUS 9 DIFFERENT RESTORATIVE SOLUTIONS** for full arch cases

» **PATIENT-SPECIFIC ABUTMENTS** tailored to your patients' needs, offering up to 30° of angulation




Quick Reference Chart

Manufacturer	Compatible Implant/ Interface	Healing Abutment	Scan Body	Transfer	Digital Analog ²	Analog ³	Verification Cylinder		Patient- Specific Abutment	ASC Ti Base	
							Engaging	Non- Engaging		Engaging	Non-Engaging
BioHorizons®	Tapered Internal 3.0mm		✓		✓	✓			✓		
	Tapered Internal 3.5mm	*	✓	*	✓	✓	*	*	✓	*See Zimmer Tapered Screw-Vent® 3.5mm	
	Tapered Internal 4.5mm	**	✓	**	✓	✓	**	**	✓	**See Zimmer Tapered Screw-Vent® 4.5mm	
Branemark®	RP 3.75/4.1mm		✓			✓	✓	✓			
CAMLOG®	Screw-Line 3.3mm		✓			✓	✓	✓			
	Screw-Line 3.8mm		✓			✓	✓	✓			
	Screw-Line 4.3mm		✓			✓	✓	✓			
	Screw-Line 5.0mm		✓			✓	✓	✓			
	Screw-Line 6.0mm		✓			✓	✓	✓			
Dentsply Sirona	Ankylos®		✓			✓	✓	✓			
	Astra TX 3.0mm		✓		✓	✓	✓		✓		
	Astra TX Aqua 3.5/4.0mm	✓	✓	✓	✓	✓	✓	✓	✓		
	Astra TX Lilac 4.5/5.0mm	✓	✓	✓	✓	✓	✓	✓	✓		
	Astra Tech EV 3.0mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Astra Tech EV 3.6mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Astra Tech EV 4.2mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Astra Tech EV 4.8mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Astra Tech EV 5.4mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Hiossen®	Mini	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Regular	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Implant Direct	Legacy™ 3.0mm		✓		✓	✓			✓		
	Legacy™ 3.5mm	*	✓	*	✓	✓	*	*	✓	*See Zimmer Tapered Screw-Vent® 3.5mm	
	Legacy™ 4.5mm	**	✓	**	✓	✓	**	**	✓	**See Zimmer Tapered Screw-Vent® 4.5mm	
	Legacy™ 5.7mm	***	✓	***	✓	✓	***	***	✓	***See Zimmer Tapered Screw-Vent® 5.7mm	
Keystone Dental	PrimaConnex™ TiLobe® 3.5mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	PrimaConnex™ TiLobe® 4.1mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	PrimaConnex™ TiLobe® 5.0mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MegaGen	AnyRidge®	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Unit		✓	✓	✓	✓	✓	N/A	✓	N/A	N/A	✓
Neodent®	Helix GM™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Neoss®	Regular Platform 3.5/4.0/4.5/5.0/5.5mm		✓		✓	✓			✓		
NobelBiocare™	NobelActive™ 3.0mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	NobelActive™ NP	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	NobelActive™ RP	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	NobelActive™ WP		✓	✓	✓	✓	✓	✓	✓		
	Tri-Lobe NP		✓	✓	✓	✓	✓	✓	✓	✓	✓
	Tri-Lobe RP		✓	✓	✓	✓	✓	✓	✓	✓	✓
	Tri-Lobe WP		✓	✓	✓	✓	✓	✓	✓		
	Tri-Lobe 6.0mm		✓		✓	✓	✓	✓	✓		
Straumann®	BLX RB/WB	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	BLX WB	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Bone Level NC	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Bone Level RC	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Tissue Level NN 3.5mm					✓					
	Tissue Level RN 4.8mm		✓		✓	✓	✓	✓	✓	✓	✓
	Tissue Level WN 6.5mm		✓		✓	✓	✓	✓	✓		
ZimVie (formerly Zimmer Dental)	Certain® 3.4mm		✓		✓	✓	✓	✓	✓	✓	✓
	Certain® 4.1mm		✓		✓	✓	✓	✓	✓	✓	✓
	Certain® 5.0mm		✓	Call	✓	✓	✓	✓	✓		
	Certain® 6.0mm		✓	Call	✓	✓	✓	✓	✓		
	Tapered Screw-Vent® 3.5mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Tapered Screw-Vent® 4.5mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Tapered Screw-Vent® 5.7mm		✓		✓	✓	✓	✓	✓	✓	✓

Quick Reference Chart (cont'd)

Manufacturer	Compatible Implant/ Interface	Ti Base			Multi- Unit Abutment			Esthetic Abutment		Gold/Plastic UCLA		Plastic UCLA	
		Engaging	6mm Engaging	Non- Engaging	Straight	17°	30°	Straight	15°	Engaging	Non- Engaging	Engaging	Non- Engaging
BioHorizons®	Tapered Internal 3.0mm												
	Tapered Internal 3.5mm	*See Zimmer Tapered Screw-Vent® 3.5mm											
	Tapered Internal 4.5mm	**See Zimmer Tapered Screw-Vent® 4.5mm											
Branemark®	RP 3.75/4.1mm	✓	✓	✓	✓	✓	✓					✓	✓
CAMLOG®	Screw-Line 3.3mm	✓	✓	✓				✓	✓				✓
	Screw-Line 3.8mm	✓	✓	✓	✓	✓	✓	✓	✓				✓
	Screw-Line 4.3mm	✓	✓	✓	✓	✓	✓	✓	✓				✓
	Screw-Line 5.0mm	✓	✓	✓	✓			✓	✓				✓
	Screw-Line 6.0mm	✓	✓	✓									✓
Dentsply Sirona	Ankylos®	✓	✓	✓	✓			✓	✓				✓
	Astra TX 3.0mm	✓	✓					✓	✓			✓	✓
	Astra TX Aqua 3.5/4.0mm	✓	✓	✓	✓	✓	✓	✓	✓			✓	✓
	Astra TX Lilac 4.5/5.0mm	✓	✓	✓	✓	✓	✓	✓	✓			✓	✓
	Astra Tech EV 3.0mm	✓	✓	✓									
	Astra Tech EV 3.6mm	✓	✓	✓									
	Astra Tech EV 4.2mm	✓	✓	✓									
	Astra Tech EV 4.8mm	✓	✓	✓									
	Astra Tech EV 5.4mm	✓	✓	✓									
Hiossen®	Mini	✓	✓	✓	✓	✓	✓	✓	✓			✓	✓
	Regular	✓	✓	✓	✓	✓	✓	✓	✓			✓	✓
Implant Direct	Legacy™ 3.0mm												
	Legacy™ 3.5mm	*See Zimmer Tapered Screw-Vent® 3.5mm											
	Legacy™ 4.5mm	**See Zimmer Tapered Screw-Vent® 4.5mm											
	Legacy™ 5.7mm	***See Zimmer Tapered Screw-Vent® 5.7mm											
Keystone Dental	PrimaConnex™ TiLobe® 3.5mm	✓	✓	✓									
	PrimaConnex™ TiLobe® 4.1mm	✓	✓	✓									
	PrimaConnex™ TiLobe® 5.0mm	✓	✓	✓									
MegaGen	AnyRidge®	✓	✓	✓	✓	✓	✓	✓	✓			✓	✓
Multi-Unit		N/A	N/A	✓	N/A	N/A	N/A	N/A	N/A	N/A	✓	N/A	N/A
Neodent®	Helix GM™	✓	✓	✓									
Neoss®	Regular Platform 3.5/4.0/4.5/5.0/5.5mm												
NobelBiocare™	NobelActive™ 3.0mm	✓	✓	✓	✓			✓	✓	✓	✓	✓	✓
	NobelActive™ NP	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	NobelActive™ RP	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	NobelActive™ WP	✓	✓	✓				✓		✓	✓	✓	✓
	Tri-Lobe NP	✓	✓	✓	✓	✓		✓	✓	✓		✓	✓
	Tri-Lobe RP	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓
	Tri-Lobe WP	✓	✓	✓	✓			✓	✓	✓		✓	✓
	Tri-Lobe 6.0mm	✓	✓	✓				✓		✓		✓	✓
Straumann®	BLX RB/WB	✓	✓	✓									
	BLX WB	✓	✓	✓									
	Bone Level NC	✓	✓	✓	✓	✓		✓	✓			✓	✓
	Bone Level RC	✓	✓	✓	✓	✓	✓	✓	✓			✓	✓
	Tissue Level NN 3.5mm												
	Tissue Level RN 4.8mm	✓	✓	✓						✓	✓		
	Tissue Level WN 6.5mm	✓	✓	✓						✓	✓		
ZimVie (formerly Zimmer Dental)	Certain® 3.4mm	✓	✓	✓	✓	✓		✓	✓			✓	✓
	Certain® 4.1mm	✓	✓	✓	✓	✓	✓	✓	✓			✓	✓
	Certain® 5.0mm	✓	✓	✓	✓	✓		✓	✓			✓	✓
	Certain® 6.0mm	✓	✓	✓				✓				✓	✓
	Tapered Screw-Vent® 3.5mm	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓
	Tapered Screw-Vent® 4.5mm	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓
	Tapered Screw-Vent® 5.7mm	✓	✓	✓								✓	✓

IMPLANT PROSTHETICS



Compatible with
3i Certain®



9001645



9001600



9001603



9001630



9001635



9001690




9001795

3i Certain® 3.4mm Compatible

Item #	Product
9001645	3.4mm Preat Choice Titanium Scan Body
9001600	3.4mm Analog
9001600-10	3.4mm Analog (10-Pack)
9001600-25	3.4mm Analog (25-Pack)
9001603	3.4mm Digital Analog
9001603-10	3.4mm Digital Analog (10-Pack)
9001630	3.4mm Engaging Verification Cylinder
9001630-10	3.4mm Engaging Verification Cylinder (10-Pack)
9001635	3.4mm Non-Engaging Verification Cylinder
9001635-10	3.4mm Non-Engaging Verification Cylinder (10-Pack)
9001690	ASC 3.4mm x 9mm Engaging Ti Base
9001795	ASC 3.4mm x 9mm Non-Engaging Ti Base
9001653	ASC 3.4mm/4.1mm/5.0mm/6.00mm Titanium Screw

3i® Certain 3.4mm Compatible

76



3i Certain® 3.4mm Compatible (cont'd)

Item #	Product
9001659	3.4mm Non-Engaging Ti Base x 4.5mm
9001661	3.4mm Engaging Ti Base x 4.5mm
9001662	3.4mm Engaging Ti Base x 6mm
9001650	<u>3.4mm / 4.1mm / 5.0mm / 6.0mm</u> Titanium Screw
9001650-10	<u>3.4mm / 4.1mm / 5.0mm / 6.0mm</u> Titanium Screw (10-Pack)
9001651	<u>3.4mm / 4.1mm / 5.0mm / 6.0mm</u> Guide Pin 20mm
8001789	<u>3.4mm / 4.1mm / 5.0mm / 6.0mm</u> Multi-Unit Titanium Screw
8001680	3.4mm 17° Multi-Unit Abutment x 2mm
8001681	3.4mm 17° Multi-Unit Abutment x 3mm
8001675	3.4mm Straight Multi-Unit Abutment x 1mm
8001676	3.4mm Straight Multi-Unit Abutment x 2mm
8001677	3.4mm Straight Multi-Unit Abutment x 3mm
8001670	3.4mm Esthetic Abutment 15° Angle, Anterior
8001671	3.4mm Esthetic Abutment Straight, Anterior
8001610	3.4mm Engaging Plastic UCLA
8001615	3.4mm Non-Engaging Plastic UCLA
9001667	3.4mm 3i Certain®-compatible Titanium Blank
9001665	3.4mm ASC 3i Certain®-compatible Titanium Blank

IMPLANT PROSTHETICS



9001745



9001703



9001730



9001735



9001790



9001795

3i Certain® 4.1mm

Item # Product

9001745	4.1mm Preat Choice Titanium Scan Body
9001700	4.1mm Analog
9001700-10	4.1mm Analog (10-Pack)
9001700-25	4.1mm Analog (25-Pack)
9001703	4.1mm Digital Analog
9001703-10	4.1mm Digital Analog (10-Pack)
9001730	4.1mm Engaging Verification Cylinder
9001730-10	4.1mm Engaging Verification Cylinder (10-Pack)
9001735	4.1mm Non-Engaging Verification Cylinder
9001735-10	4.1mm Non-Engaging Verification Cylinder (10-Pack)
9001790	ASC 4.1mm x 9mm Engaging Ti Base
9001795	ASC 4.1mm x 9mm Non-Engaging Ti Base
9001759	4.1 mm Non-Engaging Ti Base x 4.5mm
9001761	4.1 mm Engaging Ti Base x 4.5mm
9001762	4.1mm Engaging Ti Base x 6.0mm

3i Certain® 4.1mm Compatible

78




8001775




9001759




9001762



8001780



8001787



8001777



8001772




9001765

3i Certain® 4.1mm (cont'd)

Item #	Product
	3i® Screws and Guide Pins - See page 73
8001780	4.1mm 17° Multi-Unit Abutment x 2mm
8001781	4.1mm 17° Multi-Unit Abutment x 3mm
8001786	4.1mm 30° Multi-Unit Abutment x 4mm
8001787	4.1mm 30° Multi-Unit Abutment x 5mm
8001775	4.1mm Straight Multi-Unit Abutment x 1mm
8001776	4.1mm Straight Multi-Unit Abutment x 2mm
8001777	4.1mm Straight Multi-Unit Abutment x 3mm
8001779	4.1mm Straight Multi-Unit Abutment x 5mm
8001770	4.1mm Esthetic Abutment 15° Angle, Anterior
8001772	4.1mm Esthetic Abutment 15° Angle, Posterior
8001771	4.1mm Esthetic Abutment Straight, Anterior
8001773	4.1mm Esthetic Abutment Straight, Posterior
8001710	4.1mm Engaging Plastic UCLA
8001715	4.1mm Non-Engaging Plastic UCLA
9001765	4.1mm ASC 3i Certain®-compatible Titanium Blank
9001767	4.1mm 3i Certain®-compatible Titanium Blank

IMPLANT PROSTHETICS




3i® Certain 5.0mm Compatible


80

3i® Certain 5.0mm


Item #	Product
8001841	5.0mm Titanium Scan Body (for use with Ti Bases)
9001845	5.0mm Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9001800	5.0mm Analog
9001800-10	5.0mm Analog (10-Pack)
9001800-25	5.0mm Analog (25-Pack)
9001803	5.0mm Digital Analog
9001803-10	5.0mm Digital Analog (10-Pack)
9001835	5.0mm Non-Engaging Verification Cylinder
9001835-10	5.0mm Non-Engaging Verification Cylinder (10-Pack)
9001830	5.0mm Engaging Verification Cylinder
9001830-10	5.0mm Engaging Verification Cylinder (10-Pack)
8001859	5.0mm Non-Engaging Ti Base x 4.5mm
8001861	5.0mm Engaging Ti Base x 4.5mm
8001862	5.0mm Engaging Ti Base x 6mm
	Screws and Guide Pins - See page 69
8001880	5.0mm 17° Multi-Unit Abutment x 2mm
8001881	5.0mm 17° Multi-Unit Abutment x 3mm
8001875	5.0mm Straight Multi-Unit Abutment x 1mm
8001876	5.0mm Straight Multi-Unit Abutment x 2mm



8001859




8001872




8001881



8001875



8001877



8001872



8001871



9001867

3i® Certain 5.0mm (cont'd)

Item #	Product
8001877	5.0mm Straight Multi-Unit Abutment x 3mm
8001870	5.0mm Esthetic Abutment 15° Angle, Anterior
8001872	5.0mm Esthetic Abutment 15° Angle, Posterior
8001871	5.0mm Esthetic Abutment Straight, Anterior
8001873	5.0mm Esthetic Abutment Straight, Posterior
8001810	5.0mm Engaging Plastic UCLA
8001815	5.0mm Non-Engaging Plastic UCLA
9001865	5.0mm ASC 3i Certain®-compatible Titanium Blank
9001867	5.0mm 3i Certain®-compatible Titanium Blank

IMPLANT PROSTHETICS

3i® Certain 6.0mm Compatible

82



9001945




9001900



8001930



8001959



8001961



9001903



8001910



8001915

3i® Certain 6.0mm

Item #	Product
8001941	6.0mm Scan Body (for use with Ti Bases)
9001945	6.0mm Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9001900	6.0mm Analog
9001900-10	6.0mm Analog (10-Pack)
9001900-25	6.0mm Analog (25-Pack)
9001903	6.0mm Digital Analog
9001903-10	6.0mm Digital Analog (10-Pack)
8001930	6.0mm Engaging Verification Cylinder
8001935	6.0mm Non-Engaging Verification Cylinder
8001959	6.0 mm Non-Engaging Ti Base x 4.5mm
8001961	6.0 mm Engaging Ti Base x 4.5mm
8001962	6.0mm Engaging Ti Base x 6mm
	Screws and Guide Pins - See page 69
8001973	6.0mm Esthetic Abutment Straight, Posterior
8001910	6.0mm Engaging Plastic UCLA
8001915	6.0mm Non-Engaging Plastic UCLA
9001965	6.0mm ASC 3i Certain®-compatible Titanium Blank
9001967	6.0mm 3i Certain®-compatible Titanium Blank

3i® External Hex


3i External Hex 3.75 / 4.0mm	<i>Please see Nobel Branemark RP</i>
------------------------------	--------------------------------------

Ace Surgical Implant Systems

ACE Infinity Internal 3.7 / 4.1 (3.5mmD)	<i>Please see Zimmer TSV 3.5mm</i>
ACE Infinity Internal 4.7 / 5.1 (4.5mmD)	<i>Please see Zimmer TSV 4.5mm</i>
ACE Infinity TRI-CAM 3.5mm	<i>Please see Nobel Tri-Lobe NP</i>
ACE Infinity TRI-CAM 4.3mm	<i>Please see Nobel Tri-Lobe RP</i>
ACE Infinity TRI-CAM 5.0mm	<i>Please see Nobel Tri-Lobe WP</i>
ACE Infinity X-hex 3.7 / 4.0	<i>Please see Nobel Branemark RP</i>
ACE Infinity X-hex 4.75	<i>Please see Nobel Branemark WP</i>
ACE Octagon Bone Level 3.3mm	<i>Please see Nobel Active / Conical NP</i>
ACE Octagon Bone Level 4.1mm	<i>Please see Nobel Active / Conical RP</i>

IMPLANT PROSTHETICS

Ankylos® C/X (indexed) & C (non-indexed) Implant 3.5 / 4.5 / 5.5mm Compatible



Compatible with
Dentsply Ankylos®



Ankylos® C/X (indexed) & C (non-indexed) Implant 3.5 / 4.5 / 5.5mm	
Item #	Product
8004841	Titanium Scan Body <i>(for use with Ti Bases)</i>
8004800	3.5mm Analog
8004800-10	3.5mm Analog (10-Pack)
8004805	4.5mm Analog
8004805-10	4.5mm Analog (10-Pack)
8004806	5.5mm Analog
8004806-10	5.5mm Analog (10-Pack)
8004806-25	5.5mm Analog (25-Pack)
8004861	3.5mm Engaging Ti Base x 4.5mm
8004857	3.5mm Non-Engaging Ti Base x 4.5mm
8004863	3.5mm Engaging Ti Base x 6mm
8004858	4.5mm Non-Engaging Ti Base x 4.5mm
8004859	5.5mm Non-Engaging Ti Base x 4.5mm
8004851	Guide Pin



8004859



8004851



8004850



8004876



8004870



8004871



8004830




8004815

Ankylos® C/X (indexed) & C (non-indexed) Implant 3.5 / 4.5 / 5.5mm

Item #	Product
8004850	Titanium Screw
8004875	Straight Multi-Unit Abutment x 1.5mm
8004876	Straight Multi-Unit Abutment x 2.5mm
8004877	Straight Multi-Unit Abutment x 3.5mm
8004878	Straight Multi-Unit Abutment x 4.5mm
8004870	Esthetic Abutment 15° Angle, Anterior
8004872	Esthetic Abutment 15° Angle, Posterior
8004871	Esthetic Abutment Straight, Anterior
8004873	Esthetic Abutment Straight, Posterior
8004830	Engaging Verification Cylinder
8004835	Non- Engaging Verification Cylinder
8004815	Non-Engaging Plastic UCLA

IMPLANT PROSTHETICS



Compatible with
Astra® TX



8002030



9002000



8002061



8002062



8002050



8002051



8002070



8002071





8002015



9002067

Astra® 3.0mm

Item #	Product
8002041	3.0mm Titanium Scan Body (for use with Ti Bases)
9002045	3.0mm Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9002000	3.0mm Analog
9002000-10	3.0mm Analog (10-Pack)
9002000-25	3.0mm Analog (25-Pack)
9002003	3.0mm Digital Analog
9002003-10	3.0mm Digital Analog (10-Pack)
8002030	3.0mm Engaging Verification Cylinder
8002030-10	3.0mm Engaging Verification Cylinder (10-Pack)
8002061	3.0mm Engaging Ti Base x 4.5mm
8002062	3.0mm Engaging Ti Base x 6mm
8002050	3.0mm Titanium Screw
9002053	ASC 3.0mm Titanium Screw
8002050-10	3.0mm Titanium Screw (10-Pack)
8002051	3.0mm Guide Pin 20mm
8002070	3.0mm Esthetic Abutment 15° Angle, Anterior
8002071	3.0mm Esthetic Abutment Straight, Anterior
8002010	3.0mm Engaging Plastic UCLA
8002015	3.0mm Non-Engaging Plastic UCLA
9002065	3.0mm ASC Astra®-compatible Titanium Blank
9002067	3.0mm Astra®-compatible Titanium Blank



Astra® Aqua 3.5 / 4.0mm

Item #	Product
8002123	3.5 / 4.0mm Healing Abutment 4 x 5mm
8002104	3.5 / 4.0mm Closed Tray Impression Tray Coping 3.8 x 9mm
8002141	3.5 / 4.0mm Titanium Scan Body (for use with Ti Bases)
9002145	3.5 / 4.0mm Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9002100	3.5 / 4.0mm Analog
9002100-10	3.5 / 4.0mm Analog (10-Pack)
9002100-25	3.5 / 4.0mm Analog (25-Pack)
9002103	3.5/4.0mm Digital Analog
9002103-10	3.5/4.0mm Digital Analog (10-Pack)
8002130	3.5 / 4.0mm Engaging Verification Cylinder
8002130-10	3.5 / 4.0mm Engaging Verification Cylinder (10-Pack)
8002135	3.5 / 4.0mm Non-Engaging Verification Cylinder
8002135-10	3.4 / 4.0mm Non-Engaging Verification Cylinder (10-Pack)
8002159	3.5 / 4.0mm Non-Engaging Ti Base x 4.5mm
8002161	3.5 / 4.0mm Engaging Ti Base x 4.5mm
8002162	3.5 / 4.0mm Engaging Ti Base x 6mm
8002150	3.5 / 4.0mm Titanium Screw
9002153	ASC Aqua 3.5mm/5.0mm Titanium Screw
8002151	3.5 / 4.0mm Guide Pin 20mm
8002189	3.5 / 4.0mm Multi-Unit Titanium Screw
8002175	3.5 / 4.0mm Straight Multi-Unit Abutment x 1.5mm

IMPLANT PROSTHETICS

Astra® Aqua 3.5 / 4.0mm Compatible


88



8002150




8002151



8002177



8002186



8002172




8002110



9002165

Astra® Aqua 3.5 / 4.0mm (cont'd)

Item #	Product
8002176	3.5 / 4.0mm Straight Multi-Unit Abutment x 2.5mm
8002177	3.5 / 4.0mm Straight Multi-Unit Abutment x 3.5mm
8002180	3.5 / 4.0mm 17° Multi-Unit Abutment x 2.5mm
8002181	3.5 / 4.0mm 17° Multi-Unit Abutment x 3.5mm
8002185	3.5 / 4.0mm 30° Multi-Unit Abutment x 3.5mm
8002186	3.5 / 4.0mm 30° Multi-Unit Abutment x 4.5mm
8002170	3.5 / 4.0mm Esthetic Abutment 15° Angle, Anterior
8002172	3.5 / 4.0mm Esthetic Abutment 15° Angle, Posterior
8002171	3.5 / 4.0mm Esthetic Abutment Straight, Anterior
8002173	3.5 / 4.0mm Esthetic Abutment Straight, Posterior
8002110	3.5 / 4.0mm Engaging Plastic UCLA
8002115	3.5 / 4.0mm Non-Engaging Plastic UCLA
9002165	3.5/4.0mm ASC Astra®-compatible Aqua Titanium Blank
9002167	3.5/4.0mm Astra®-compatible Aqua Titanium Blank




Astra® Lilac 4.5 / 5.0mm

Item #	Product
8002223	4.5 / 5.0mm Healing Abutment 5 x 5mm
8002204	4.5 / 5.0mm Closed Tray Impression Tray Copping 4.5 x 9mm
8002241	4.5 / 5.0mm Titanium Scan Body (for use with Ti Bases)
9002245	4.5 / 5.0mm Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9002200	4.5 / 5.0mm Analog
9002200-10	4.5 / 5.0mm Analog (10-Pack)
9002200-25	4.5 / 5.0mm Analog (25-Pack)
9002203	4.5/5.0mm Digital Analog
9002203-10	4.5/5.0mm Digital Analog (10-Pack)
8002230	4.5 / 5.0mm Engaging Verification Cylinder
8002230-10	4.5 / 5.0mm Engaging Verification Cylinder (10-Pack)
8002235	4.5 / 5.0mm Non-Engaging Verification Cylinder
8002235-10	4.5 / 5.0mm Non-Engaging Verification Cylinder (10-Pack)
8002259	4.5 / 5.0mm Non-Engaging Ti Base x 4.5mm
8002261	4.5 / 5.0mm Engaging Ti Base x 4.5mm
8002262	4.5 / 5.0mm Engaging Ti Base x 6mm
8002250	4.5 / 5.0mm Titanium Screw
8002250-10	4.5 / 5.0mm Titanium Screw (10-Pack)
9002253	ASC Lilac 4.5mm/6.0mm Titanium Screw
8002251	4.5 / 5.0mm Guide Pin 20mm
8002289	4.5 / 5.0mm Multi-Unit Titanium Screw
8002280	4.5 / 5.0mm 17° Multi-Unit Abutment x 2.5mm

IMPLANT PROSTHETICS

Astra® Lilac 4.5 / 5.0mm Compatible

90



8002285



8002276



8002270




8002210



8002215




9002265




8002273

Astra® Lilac 4.5 / 5.0mm (cont'd)



Item #	Product
8002281	4.5 / 5.0mm 17° Multi-Unit Abutment x 3.5mm
8002285	4.5 / 5.0mm 30° Multi-Unit Abutment x 3.5mm
8002286	4.5 / 5.0mm 30° Multi-Unit Abutment x 4.5mm
8002275	4.5 / 5.0mm Straight Multi-Unit Abutment x 1.5mm
8002276	4.5 / 5.0mm Straight Multi-Unit Abutment x 2.5mm
8002277	4.5 / 5.0mm Straight Multi-Unit Abutment x 3.5mm
8002278	4.5 / 5.0mm Straight Multi-Unit Abutment x 4.5mm
8002270	4.5 / 5.0mm Esthetic Abutment 15° Angle, Anterior
8002272	4.5 / 5.0mm Esthetic Abutment 15° Angle, Posterior
98002271	4.5 / 5.0mm Esthetic Abutment Straight, Anterior
8002273	4.5 / 5.0mm Esthetic Abutment Straight, Posterior
8002210	4.5 / 5.0mm Engaging Plastic UCLA
8002215	4.5 / 5.0mm Non-Engaging Plastic UCLA
9002265	4.5/5.0mm ASC Astra®-compatible Lilac Titanium Blank
9002267	4.5/5.0mm Astra®-compatible Lilac Titanium Blank



Compatible with
Astra® Tech EV



Astra® Tech EV 3.0mm	
Item #	Product
9006445	3.0mm Preat Choice Titanium Scan Body
9006400	3.0mm Analog
9006400-10	3.0mm Analog (10-Pack)
9006400-25	3.0mm Analog (25-Pack)
9006401	3.0mm Open-Tray Transfer
9006404	3.0mm Closed-Tray Transfer
9006452	3.0mm Closed-Tray Pin
9006403	3.0mm Digital Analog
9006403-10	3.0mm Digital Analog (10-Pack)
9006421	3.0mm x 2mm Healing Abutment (Wide)
9006422	3.0mm x 2mm Healing Abutment (Wide)
9006423	3.0mm x 4mm Healing Abutment
9006424	3.0mm x 4mm Healing Abutment (Wide)
9006430	3.0mm Engaging Verification Cylinder
9006430-10	3.0mm Engaging Verification Cylinder (10-Pack)
9006435	3.0mm Non-Engaging Verification Cylinder
9006435-10	3.0mm Non-Engaging Verification Cylinder (10-Pack)
9006490	ASC 3.0mm x 9mm Engaging Ti Base
9006495	ASC 3.0mm x 9mm Non-Engaging Ti Base



Astra® Tech EV 3.0mm (cont'd)	
Item #	Product
9006453	ASC 3.0mm Titanium Screw
9006459	3.0mm Non-Engaging Ti Base x 4.5mm
9006461	3.0mm Engaging Ti Base x 4.5mm
9006462	3.0mm Engaging Ti Base x 6mm
9006450	3.0mm Titanium Screw
9006450-10	3.0mm Titanium Screw (10-Pack)
9006451	3.0mm Guide Pin
9006465	3.0mm ASC Astra® EV-compatible Titanium Blank
9006467	3.0mm Astra® EV-compatible Titanium Blank



Astra® Tech EV 3.6mm

Item #	Product
9006545	3.6mm Preat Choice Titanium Scan Body
9006500	3.6mm Analog
9006500-10	3.6mm Analog (10-Pack)
9006500-25	3.6mm Analog (25-Pack)
9006501	3.6mm Open-Tray Transfer
9006504	3.6mm Closed-Tray Transfer
9006552	3.6mm Closed-Tray Pin
9006503	3.6mm Digital Analog
9006503-10	3.6mm Digital Analog (10-Pack)
9006521	3.6mm X 2mm Healing Abutment
9006522	3.6mm X 2mm Healing Abutment (Wide)
9006523	3.6mm X 4mm Healing Abutment
9006524	3.6mm X 4mm Healing Abutment (Wide)
9006530	3.6mm Engaging Verification Cylinder
9006530-10	3.6mm Engaging Verification Cylinder (10-Pack)
9006535	3.6mm Non-Engaging Verification Cylinder
9006535-10	3.6mm Non-Engaging Verification Cylinder (10-Pack)
9006590	ASC 3.6mm x 9mm Engaging Ti Base
9006595	ASC 3.6mm x 9mm Non-Engaging Ti Base



Astra® Tech EV 3.6mm (cont'd)	
Item #	Product
9006553	ASC 3.6mm Titanium Screw
9006559	3.6 X 4.5mm Non-Engaging Ti Base
9006561	3.6mm x 4.5mm Engaging Ti Base
9006562	3.6mm x 6mm Engaging Ti Base
9006550	3.6mm Titanium Screw
9006550-10	3.6mm Titanium Screw (10-Pack)
9006551	3.6mm Guide Pin
9006565	3.6mm ASC Astra® EV-compatible Titanium Blank
9006567	3.6mm Astra® EV-compatible Titanium Blank



9006645



9006600



9006630



9006635



9006659



9006603

Astra® Tech EV 4.2mm

Item #	Product
9006645	4.2mm Preat Choice Titanium Scan Body
9006600	4.2mm Analog
9006600-10	4.2mm Analog (10-Pack)
9006600-25	4.2mm Analog (25-Pack)
9006601	4.2mm Open-Tray Transfer
9006604	4.2mm Closed-Tray Transfer
9006652	4.2mm Closed-Tray Pin
9006603	4.2mm Digital Analog
9006603-10	4.2mm Digital Analog (10-Pack)
9006621	4.2mm x 2mm Healing Abutment
9006622	4.2mm x 2mm Healing Abutment (Wide)
9006623	4.2mm x 4mm Healing Abutment
9006624	4.2mm x 4mm Healing Abutment (Wide)
9006630	4.2mm Engaging Verification Cylinder
9006630-10	4.2mm Engaging Verification Cylinder (10-Pack)
9006635	4.2mm Non-Engaging Verification Cylinder
9006635-10	4.2mm Non-Engaging Verification Cylinder (10 Pack)
9006690	4.2mm x 9mm ASC Engaging Ti Base
9006695	4.2mm x 9mm ASC Non-Engaging Ti Base



Astra® Tech EV 4.2mm (cont'd)	
Item #	Product
9006653	ASC 4.2mm Titanium Screw
9006659	4.2mm x 4.5mm Non-Engaging Ti Base
9006661	4.2mm x 4.5mm Engaging Ti Base
9006662	4.2mm x 6mm Engaging Ti Base
9006650	4.2mm Titanium Screw
9006650-10	4.2mm Titanium Screw (10-Pack)
9006651	4.2mm Guide Pin
9006665	4.2mm ASC Astra® EV-compatible Titanium Blank
9006667	4.2mm Astra® EV-compatible Titanium Blank



9006745



9006700



9006701



9006703



9006735




9006730

Astra® Tech EV 4.8mm

Item #	Product
9006745	4.8mm Preat Choice Titanium Scan Body
9006700	4.8mm Analog
9006700-10	4.8mm Analog (10-Pack)
9006700-25	4.8mm Analog (25-Pack)
9006701	4.8mm Open-Tray Transfer
9006704	4.8mm Closed-Tray Transfer
9006752	4.8mm Closed-Tray Pin
9006703	4.8mm Digital Analog
9006703-10	4.8mm Digital Analog (10-Pack)
9006721	4.8mm X 2mm Healing Abutment
9006722	4.8mm X 2mm Healing Abutment (Wide)
9006723	4.8mm X 4mm Healing Abutment
9006724	4.8mm X 4mm Healing Abutment (Wide)
9006730	4.8mm Engaging Verification Cylinder
9006730-10	4.8mm Engaging Verification Cylinder (10-Pack)
9006735	4.8mm Non-Engaging Verification Cylinder
9006735-10	4.8mm Non-Engaging Verification Cylinder (10-Pack)

IMPLANT PROSTHETICS

Astra® Tech EV 4.8mm Compatible



Astra® Tech EV 4.8mm (cont'd)	
Item #	Product
9006790	ASC 4.8mm x 9mm Engaging Ti Base
9006795	ASC 4.8mm x 9mm Non-Engaging Ti Base
9006753	ASC 4.8mm/5.4mm Titanium Screw
9006759	4.8mm x 4.5mm Non-Engaging Ti Base
9006761	4.8mm x 4.5mm Engaging Ti Base
9006762	4.8mm x 6mm Engaging Ti Base
9006750	4.8mm / 5.4mm Titanium Screw
9006750-10	4.8mm / 5.4mm Titanium Screw (10-Pack)
9006751	4.8mm / 5.4mm Guide Pin
9006765	4.8mm ASC Astra® EV-compatible Titanium Blank
9006767	4.8mm Astra® EV-compatible Titanium Blank



9006845



9006800



9006830



9006835



9006859



9006803

Astra® Tech EV 5.4mm

Item #	Product
9006845	5.4mm Preat Choice Titanium Scan Body
9006800	5.4mm Analog
9006800-10	5.4mm Analog (10-Pack)
9006800-25	5.4mm Analog (25-Pack)
9006801	5.4mm Open-Tray Transfer
9006804	5.4mm Closed-Tray Transfer
9006852	5.4mm Closed-Tray Pin
9006803	5.4mm Digital Analog
9006803-10	5.4mm Digital Analog (10-Pack)
9006821	5.4mm X 2mm Healing Abutment
9006822	5.4mm X 2mm Healing Abutment (Wide)
9006823	5.4mm X 4mm Healing Abutment
9006824	5.4mm X 4mm Healing Abutment (Wide)
9006830	5.4mm Engaging Verification Cylinder
9006830-10	5.4mm Engaging Verification Cylinder (10-Pack)
8006850-10	5.4 Titanium Screw (10-Pack)
9006835	5.4mm Non-Engaging Verification Cylinder
9006835-10	5.4mm Non-Engaging Verification Cylinder (10-Pack)

IMPLANT PROSTHETICS



9006890



9006867

Astra® Tech EV 5.4mm (cont'd)

Item #	Product
9006890	ASC 5.4mm x 9mm Engaging Ti Base
9006895	ASC 5.4mm x 9mm Non-Engaging Ti Base
9006753	ASC 4.8mm/5.4mm Titanium Screw
9006859	5.4mm x 4.5mm Non-Engaging Ti Base
9006861	5.4mm x 4.5mm Engaging Ti Base
9006862	5.4mm x 6mm Engaging Ti Base
9006865	5.4mm ASC Astra® EV-compatible Titanium Blank
9006867	5.4mm Astra® EV-compatible Titanium Blank

Screws and Guide Pins - [See page 90](#)

BIO-LOK Implant Systems


BIO-LOK External Hex	<i>Please see Nobel Branemark RP</i>
----------------------	--------------------------------------

Blue Sky Bio Implant Systems

Blue Sky Bio Internal Hex 3.5mm Green	<i>Please see Zimmer TSV 3.5mm</i>
Blue Sky Bio Internal Hex 4.5mm Purple	<i>Please see Zimmer TSV 4.5mm</i>
Blue Sky Bio Tri-Lobe NP	<i>Please see Nobel Tri-Lobe NP</i>
Blue Sky Bio Tri-Lobe RP	<i>Please see Nobel Tri-Lobe RP</i>
Blue Sky Bio Tri-Lobe WP	<i>Please see Nobel Tri-Lobe WP</i>
Blue Sky Bio Tri-Lobe 6.0	<i>Please see Nobel Tri-Lobe 6.0</i>
Blue Sky Bio MAX	<i>Please see NobelActive™ / Conical NP</i>
Blue Sky Bio Three	<i>Please see Astra 3.0mm</i>
Blue Sky Bio CONUS 12 3.5 / 4.0	<i>Please see Astra Aqua 3.5 / 4.0mm</i>
Blue Sky Bio CONUS 12 4.5 / 5.0	<i>Please see Astra Lilac 4.5 / 5.0mm</i>
Blue Sky Bio One Stage RN (4.1mm)	<i>Please see Straumann Tissue Level RN</i>

Astra® Tech EV 5.4mm Compatible

100



Compatible with
BioHorizons®



BioHorizons® 3.0mm

Item #	Product
9005245	3.0mm Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9005200	3.0mm Analog
9005200-10	3.0mm Analog (10-Pack)
9005200-25	3.0mm Analog (25-Pack)
9005203	3.0mm Digital Analog
9005253	ASC 3.0mm Titanium Screw
9005265	3.0mm ASC BioHorizons®-compatible Titanium Blank
9005267	3.0mm BioHorizons®-compatible Titanium Blank

BioHorizons® 3.5mm

Item #	Product
9005445	3.5mm Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9005400	3.5mm Analog
9005400-10	3.5mm Analog (10-Pack)
9005400-25	3.5mm Analog (25-Pack)
9005403	3.5mm Digital Analog
9005453	ASC 3.5mm/4.5mm Titanium Screw
9005465	3.5mm ASC BioHorizons®-compatible Titanium Blank
9005467	3.5mm BioHorizons®-compatible Titanium Blank

See Zimmer TSV for additional options

BioHorizons® 4.5mm

Item #	Product
9005545	4.5mm Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9005500	4.5mm Analog
9005500-10	4.5mm Analog (10-Pack)
9005500-25	4.5mm Analog (25-Pack)
9005503	4.5mm Digital Analog
9005565	4.5mm ASC BioHorizons®-compatible Titanium Blank
9005567	4.5mm BioHorizons®-compatible Titanium Blank

See Zimmer TSV for additional options

BioHorizons® 3.0 Implant Compatible

101

Compatible with Brånemark®



9001000



8001035



8001030




8001062



8001051



8001089




8001087



8001091

Brånemark® RP 3.75 / 4.1mm

Item #	Product
8001041	Titanium Scan Body <i>(for use with Ti Bases)</i>
9001000	Analog
9001000-10	Analog (10-Pack)
9001000-25	Analog (25-Pack)
8001035	Non-Engaging Verification Cylinder
8001035-10	Non-Engaging Verification Cylinder (10-Pack)
8001030	Engaging Verification Cylinder
8001030-10	Engaging Verification Cylinder (10-Pack)
8001059	Non-Engaging Ti Base x 4.5mm
8001061	Engaging Ti Base x 4.5mm
8001062	Engaging Ti Base x 6mm
9001050	Titanium Screw
9001050-10	Titanium Screw (10-Pack)
8001051	Guide Pin 20mm
8001089	Multi-Unit Titanium Screw
8001085	17° Multi-Unit Abutment x 2mm
8001086	17° Multi-Unit Abutment x 3mm
8001087	17° Multi-Unit Abutment x 4mm
8001091	30° Multi-Unit Abutment x 4mm
8001092	30° Multi-Unit Abutment x 5mm



Brånemark® RP 3.75 / 4.1mm (cont'd)	
Item #	Product
8001081	Straight Multi-Unit Abutment x 1mm
8001082	Straight Multi-Unit Abutment x 2mm
8001083	Straight Multi-Unit Abutment x 3mm
8001010	Engaging Plastic UCLA
8001015	Non-Engaging Plastic UCLA

IMPLANT PROSTHETICS

Compatible with
CAMLOG® Screw-Line



8003241



8003200



8003230



8003235



8003259



8003261



8003262



8003270




8003271



8003215

CAMLOG Screw-Line 3.3mm

Item #	Product
8003241	3.3mm Titanium Scan Body (for use with Ti Bases)
8003200	3.3mm Analog
8003200-10	3.3mm Analog (10-Pack)
8003230	3.3mm Engaging Verification Cylinder
8003230-10	3.3mm Engaging Verification Cylinder (10-Pack)
8003235	3.3mm Non-Engaging Verification Cylinder
8003235-10	3.3mm Non-Engaging Verification Cylinder (10-Pack)
8003259	3.3mm x 4.5mm Non-Engaging Ti Base
8003261	3.3mm Engaging Ti Base x 4.5mm
8003262	3.3mm x 6mm Engaging Ti Base
8003250	<u>3.3mm / 3.8mm / 4.3mm</u> Titanium Screw
8003250-10	<u>3.3mm / 3.8mm / 4.3mm</u> Titanium Screw (10-Pack)
8003251	<u>3.3mm / 3.8mm / 4.3mm</u> Guide Pin 20mm
8003270	3.3mm Esthetic Abutment 15° Angle, Anterior
8003271	3.3mm Esthetic Abutment Straight, Anterior
8003215	3.3mm Non-Engaging Plastic UCLA



CAMLOG Screw-Line 3.8mm

Item #	Product
8003341	3.8mm Titanium Scan Body (for use with Ti Bases)
8003300	3.8mm Analog
8003300-10	3.8mm Analog (10-Pack)
8003335	3.8mm Non-Engaging Verification Cylinder
8003335-10	3.8mm Non-Engaging Verification Cylinder (10-Pack)
8003330	3.8mm Engaging Verification Cylinder
8003330-10	3.8mm Engaging Verification Cylinder (10-Pack)
8003359	3.8mm x 4.5mm Non-Engaging Ti Base
8003361	3.8mm x 4.5mm Engaging Ti Base
8003362	3.8mm x 6mm Engaging Ti Base
	Screws and Guide Pins - See page 96
8003489	3.8mm / 4.3mm Multi-Unit Titanium Screw
8003375	3.8mm x 1mm Straight Multi-Unit Abutment
8003376	3.8mm x 2mm Straight Multi-Unit Abutment
8003377	3.8mm x 3mm Straight Multi-Unit Abutment
8003380	3.8mm x 2mm 17° Multi-Unit Abutment



CAMLOG Screw-Line 3.8mm (cont'd)	
Item #	Product
8003381	3.8mm x 3mm 17° Multi-Unit Abutment
8003385	3.8mm x 3mm 30° Multi-Unit Abutment
8003387	3.8mm x 5mm 30° Multi-Unit Abutment
8003370	3.8mm Esthetic Abutment 15° Angle, Anterior
8003371	3.8mm Esthetic Abutment Straight, Anterior
8003315	3.8mm Non-Engaging Plastic UCLA



8003441



8003400




8003430



8003435



8003459



8003461



8003462



8003480



8003487

CAMLOG Screw-Line 4.3mm

Item #	Product
8003441	4.3mm Titanium Scan Body (for use with Ti Bases)
8003400	4.3mm Analog
8003400-10	4.3mm Analog (10-Pack)
8003430	4.3mm Engaging Verification Cylinder
8003430-10	4.3mm Engaging Verification Cylinder (10-Pack)
8003435	4.3mm Non-Engaging Verification Cylinder
8003435-10	4.3mm Non-Engaging Verification Cylinder (10-Pack)
8003459	4.3mm x 4.5mm Non-Engaging Ti Base
8003461	4.3mm x 4.5mm Engaging Ti Base
8003462	4.3mm x 6mm Engaging Ti Base
	Screws and Guide Pins - See page 96
8003480	4.3mm 17° Multi-Unit Abutment x 2mm
8003481	4.3mm 17° Multi-Unit Abutment x 3mm
8003483	4.3mm 17° Multi-Unit Abutment x 5mm
8003485	4.3mm 30° Multi-Unit Abutment x 3mm
8003487	4.3mm 30° Multi-Unit Abutment x 5mm




8003476



8003470



8003471



8003473



8003415

CAMLOG Screw-Line 4.3mm (cont'd)

Item #	Product
8003475	4.3mm Straight Multi-Unit Abutment x 1mm
8003476	4.3mm Straight Multi-Unit Abutment x 2mm
8003477	4.3mm Straight Multi-Unit Abutment x 3mm
8003479	4.3mm Straight Multi-Unit Abutment x 5mm
8003470	4.3mm Esthetic Abutment 15° Angle, Anterior
8003472	4.3mm Esthetic Abutment 15° Angle, Posterior
8003471	4.3mm Esthetic Abutment Straight, Anterior
8003473	4.3mm Esthetic Abutment Straight, Posterior
8003415	4.3mm Non-Engaging Plastic UCLA




8003541



8003500



8003530



8003535



8003562



8003555




8003551



8003575



8003576



8003571



8003515

CAMLOG Screw-Line 5.0mm

Item #	Product
8003541	5.0mm Titanium Scan Body (for use with Ti Bases)
8003500	5.0mm Analog
8003500-10	5.0mm Analog (10-Pack)
8003530	5.0mm Engaging Verification Cylinder
8003530-10	5.0mm Engaging Verification Cylinder (10-Pack)
8003535	5.0mm Non-Engaging Verification Cylinder
8003535-10	5.0mm Non-Engaging Verification Cylinder (10-Pack)
8003559	5.0mm Non-Engaging Ti Base x 4.5mm
8003561	5.0mm Engaging Ti Base x 4.5mm
8003562	5.0mm Engaging Ti Base x 6mm
8003555	<u>5.0mm / 6.0mm</u> Titanium Screw
8003555-10	<u>5.0mm / 6.0mm</u> Titanium Screw (10-Pack)
8003551	<u>5.0mm / 6.0mm</u> Guide Pin 20mm
8003575	5.0mm Straight Multi-Unit Abutment x 1mm
8003576	5.0mm Straight Multi-Unit Abutment x 2mm
8003577	5.0mm Straight Multi-Unit Abutment x 3mm
8003571	5.0mm Esthetic Abutment 15° Angle, Posterior
8003573	5.0mm Esthetic Abutment Straight, Posterior
8003515	5.0mm Non-Engaging Plastic UCLA

IMPLANT PROSTHETICS

CAMLOG Screw-Line 6.0mm Compatible



8003641



8003600



8003630



8003661




8003615

CAMLOG Screw-Line 6.0mm

Item #	Product
8003641	6.0mm Titanium Scan Body (for use with Ti Bases)
8003600	6.0mm Analog
8003600-10	6.0mm Analog (10-Pack)
8003635	6.0mm Non-Engaging Verification Cylinder
8003630	6.0mm Engaging Verification Cylinder
8003659	6.0mm Non-Engaging Ti Base x 4.5mm
8003661	6.0mm Engaging Ti Base x 4.5mm
8003662	6.0mm Engaging Ti Base x 6mm
	Screws and Guide Pins - See page 101
8003615	6.0mm Non-Engaging Plastic UCLA

Glidewell Laboratories Implant Systems

Glidewell Inclusive 3.5mm	Please see Zimmer TSV 3.5mm
Glidewell Inclusive 3.5mm	Please see Zimmer TSV 4.5mm
Glidewell Hahn Tapered 3.0	Please see Nobel Active / Conical 3.0
Glidewell Hahn Tapered 3.5 / 4.3	Please see Nobel Active / Conical NP
Glidewell Hahn Tapered 5.0	Please see Nobel Active / Conical RP
Glidewell Hahn Tapered 7.0	Please see Nobel Active / Conical WP



**Compatible with
Hiossen® ETIII**



Hiossen Mini 3.5mm

Item #	Product
8005021	HG Mini Healing Abutment 4 x 3mm
8005023	HG Mini Healing Abutment 4 x 5mm
8005004	HG Mini Closed Tray Impression Coping 3.8 x 9mm
9005045	HG Mini Preat Choice Titanium Scan Body
9005100	HG Mini Analog
9005100-10	HG Mini Analog (10-Pack)
9005100-25	HG Mini Analog (25-Pack)
9005103	HG Mini Digital Analog
9005103-10	HG Mini Digital Analog (10-Pack)
9005035	HG Mini Non-Engaging Verification Cylinder
9005035-10	HG Mini Non-Engaging Verification Cylinder (10-Pack)
9005030	HG Mini Engaging Verification Cylinder
9005030-10	HG Mini Engaging Verification Cylinder (10-Pack)
9005090	ASC HG Mini x 9mm Engaging Ti Base
9005095	ASC HG Mini x 9mm Non-Engaging Ti Base




9005061




9005059



9005051



8005080




8005086




8005076



8005070



8005071



8005073




8005010



9005065

Hiossen Mini 3.5mm (cont'd)

Item #	Product
9005053	ASC HG Mini Titanium Screw
9005061	HG Mini Engaging Ti Base x 4.5mm
9005062	HG Mini Engaging Ti Base x 6mm
9005050	HG Mini Titanium Screw
9005050-10	HG Mini Titanium Screw (10-Pack)
9005051	HG Mini Guide Pin
9005089	HG Mini Multi Unit Titanium Screw
8005080	HG Mini 17° Multi-Unit Abutment x 2.5mm
8005081	HG Mini 17° Multi-Unit Abutment x 3.5mm
8005085	HG Mini 30° Multi-Unit Abutment x 3.5mm
8005086	HG Mini 30° Multi-Unit Abutment x 4.5mm
8005075	HG Mini Straight Multi-Unit Abutment x 1.5mm
8005076	HG Mini Straight Multi-Unit Abutment x 2.5mm
8005077	HG Mini Straight Multi-Unit Abutment x 3.5mm
8005070	HG Mini Esthetic Abutment 15° Angle, Anterior
8005072	HG Mini Esthetic Abutment 15° Angle, Posterior
8005071	HG Mini Esthetic Abutment Straight, Anterior
8005073	HG Mini Esthetic Abutment Straight, Posterior
8005010	HG Mini Engaging Plastic UCLA
8005015	HG Mini Non-Engaging Plastic UCLA
9005065	ASC Hiossen® HG-compatible Mini Titanium Blank
9005067	Hiossen® HG-compatible Mini Titanium Blank



Hiossen Regular 4.0 / 4.5 / 5.0mm

Item #	Product
8005118	HG Regular Healing Abutment 5 x 3mm
8005119	HG Regular Healing Abutment 5 x 5mm
8005121	HG Regular Healing Abutment 6 x 3mm
8005123	HG Regular Healing Abutment 6 x 5mm
8005104	HG Regular Closed Tray Impression Coping 4.5 x 9mm
9005145	HG Regular Preat Choice Titanium Scan Body
9005000	HG Regular Analog
9005000-10	HG Regular Analog (10-Pack)
9005003	HG Regular Digital Analog
9005003-10	HG Regular Digital Analog (10-Pack)
9005130	HG Regular Engaging Verification Cylinder
9005130-10	HG Regular Engaging Verification Cylinder (10-Pack)
9005135	HG Regular Non-Engaging Verification Cylinder
9005135-10	HG Regular Non-Engaging Verification Cylinder (10-Pack)
9005190	ASC HG Regular x 9mm Engaging Ti Base
9005195	ASC HG Regular x 9mm Non-Engaging Ti Base
9005153	ASC Hiossen® HG-compatible Regular Titanium Screw
9005159	HG Regular Non-Engaging Ti Base x 4.5mm
9005161	HG Regular Engaging Ti Base x 4.5mm
9005162	HG Regular Engaging Ti Base x 6mm

Preat Angled Screw Channel Titanium Bases

EXPANDED CUSTOMIZATION. SIMPLIFIED INVENTORY.

Designed with a 9 mm height,
customizable to 7 mm or 5 mm
for optimal fit


Blasted surface and
retention grooves improve
cement adhesion

Gold hue anodization for
improved aesthetics

With 29 interfaces, Preat's ASC Ti Bases provide broad compatibility
with the leading implants on the market. Available in engaging and non-
engaging implant-level and multi-unit options.

Variable screw channel
provides **0-30°** angulation for
easy adjustment to optimize
esthetics and accommodate
challenging cases










Cement groove helps
prevent cement overflow
when adhering crown



SCAN OR [CLICK HERE](#)
to Explore the Collection

Hiossen Regular 4.0 / 4.5 / 5.0mm (cont'd)


Item #	Product
9005150	HG Regular Titanium Screw
9005150-10	HG Regular Titanium Screw (10-Pack)
9005151	HG Regular Guide Pin
8005189	HG Regular Multi Unit Titanium Screw
8005180	HG Regular 17° Multi-Unit Abutment x 2.5mm
8005181	HG Regular 17° Multi-Unit Abutment x 3.5mm
8005185	HG Regular 30° Multi-Unit Abutment x 3.5mm
8005186	HG Regular 30° Multi-Unit Abutment x 4.5mm
8005175	HG Regular Straight Multi-Unit Abutment x 1.5mm
8005176	HG Regular Straight Multi-Unit Abutment x 2.5mm
8005177	HG Regular Straight Multi-Unit Abutment x 3.5mm
8005178	HG Regular Straight Multi-Unit Abutment x 4.5mm
8005170	HG Regular Esthetic Abutment 15° Angle, Anterior
8005172	HG Regular Esthetic Abutment 15° Angle, Posterior
8005171	HG Regular Esthetic Abutment Straight, Anterior
8005173	HG Regular Esthetic Abutment Straight, Posterior
8005115	HG Regular Non-Engaging Plastic UCLA
8005110	HG Regular Engaging Plastic UCLA
9005165	ASC Hiossen® HG-compatible Regular Titanium Blank
9005167	Hiossen® HG-compatible Regular Titanium Blank



Hiossen Implant Systems

Hiossen Plus Wide 5.1mm	Please see Nobel Branemark WP
Hiossen Regular 4.1mm	Please see Nobel Branemark RP
Hiossen SS Regular 4.8mm	Please see Straumann Tissue Level RN

IMPLANT PROSTHETICS



Compatible with Implant
Direct Legacy™



Implant Direct Legacy™ 3.0mm

Item #	Product
9005345	Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9005300	Analog
9005300-10	Analog (10-Pack)
9005300-25	Analog (25-Pack)
9005303	Digital Analog
9005353	ASC 3.0mm Titanium Screw
9005365	3.0mm ASC Legacy™-compatible Titanium Blank
9005367	3.0mm Legacy™-compatible Titanium Blank

Implant Direct Legacy™ 3.5mm

Item #	Product
9005645	Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9005600	Analog
9005600-10	Analog (10-Pack)
9005600-25	Analog (25-Pack)
9005603	Digital Analog
9005653	ASC 3.5mm/4.5mm/5.7mm Titanium Screw
9005665	3.5mm ASC Legacy™-compatible Titanium Blank
9005667	3.5mm Legacy™-compatible Titanium Blank

See Zimmer TSV 3.5mm for additional options

Implant Direct Legacy Compatible



Implant Direct Legacy™ 4.5mm

Item #	Product
9007545	Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9007500	Analog
9007500-10	Analog (10-Pack)
9007500-25	Analog (25-Pack)
9007503	Digital Analog
9007565	4.5mm ASC Legacy™-compatible Titanium Blank
9007567	4.5mm Legacy™-compatible Titanium Blank

Implant Direct Legacy™ 5.7mm


Item #	Product
9008545	Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9008500	Analog
9008500-10	Analog (10-Pack)
9008500-25	Analog (25-Pack)
9008503	Digital Analog
9008565	5.7mm ASC Legacy™-compatible Titanium Blank
9008567	5.7mm Legacy™-compatible Titanium Blank

See Zimmer TSV for additional options


Implant Direct Implant Systems

Implant Direct InterActive™ /Simply Iconic™ NP	Please see NobelActive / Conical NP
Implant Direct InterActive™ /Simply Iconic™ RP	Please see NobelActive / Conical RP
Implant Direct RePlant™ / RePlus™ / ReActive NP™	Please see Nobel Tri-Lobe NP
Implant Direct RePlant™ / RePlus™ / ReActive RP™	Please see Nobel Tri-Lobe RP
Implant Direct RePlant™ / RePlus™ / ReActive WP™	Please see Nobel Tri-Lobe WP
Implant Direct SwishPlant™ / SwishPlus™ (4.8mmD)	Please see Straumann Tissue Level RN

IMPLANT PROSTHETICS



Compatible with
Keystone® TiLobe®



Keystone TiLobe 3.5mm

Item #	Product
9005745	3.5mm Preat Choice Titanium Scan Body
9005700	3.5mm Analog
9005700-10	3.5mm Analog (10-Pack)
9005700-25	3.5mm Analog (25-Pack)
9005701	3.5mm Open-Tray Transfer
9005704	3.5mm Closed-Tray Transfer
9005752	3.5mm/4.1mm/5.0mm Closed-Tray Pin
9005703	3.5mm Digital Analog
9005703-10	3.5mm Digital Analog (10-Pack)
9005721	3.5mm X 2mm Healing Abutment
9005722	3.5mm X 2mm Healing Abutment (Wide)
9005723	3.5mm X 4mm Healing Abutment
9005724	3.5mm X 4mm Healing Abutment (Wide)
9005730	3.5mm Engaging Verification Cylinder
9005730-10	3.5mm Engaging Verification Cylinder (10-Pack)
9005735	3.5mm Non-Engaging Verification Cylinder
9005735-10	3.5mm Non-Engaging Verification Cylinder (10-Pack)
9005790	ASC 3.5mm x 9mm Engaging Ti Base
9005795	ASC 3.5mm x 9mm Non-Engaging Ti Base
9005753	ASC 3.5mm/4.1mm/5.0mm Titanium Screw
9005759	3.5mm x 4.5mm Non-Engaging Ti Base
9005761	3.5mm x 4.5mm Engaging Ti Base
9005762	3.5mm x 6mm Engaging Ti Base
9005750	3.5mm / 4.1mm / 5.0mm Titanium Screw
9005751	3.5mm / 4.1mm / 5.0mm Guide Pin
9005765	3.5mm ASC Keystone TiLobe®-compatible Titanium Blank
9005767	3.5mm Keystone TiLobe®-compatible Titanium Blank

Keystone TiLobe 3.5mm Compatible



Keystone TiLobe 4.1mm

Item #	Product
9005845	4.1mm Preat Choice Titanium Scan Body
9005800	4.1mm Analog
9005800-10	4.1mm Analog (10-Pack)
9005800-25	4.1mm Analog (25-Pack)
9005801	4.1mm Open-Tray Transfer
9005803	4.1mm Digital Analog
9005803-10	4.1mm Digital Analog (10-Pack)
9005804	4.1mm Closed-Tray Transfer
9005821	4.1mm X 2mm Healing Abutment
9005822	4.1mm X 2mm Healing Abutment (Wide)
9005823	4.1mm X 4mm Healing Abutment
9005824	4.1mm X 4mm Healing Abutment (Wide)
9005830	4.1mm Engaging Verification Cylinder
9005830-10	4.1mm Engaging Verification Cylinder (10-Pack)
9005835	4.1mm Non-Engaging Verification Cylinder
9005835-10	4.1mm Non-Engaging Verification Cylinder (10-Pack)
9005890	ASC 4.1mm x 9mm Engaging Ti Base
9005895	ASC 4.1mm x 9mm Non-Engaging Ti Base
9005859	4.1mm x 4.5mm Non-Engaging Ti Base
9005861	4.1mm x 4.5mm Engaging Ti Base
9005862	4.1mm x 6mm Engaging Ti Base
9005865	4.1mm ASC Keystone TiLobe®-compatible Titanium Blank
9005867	4.1mm Keystone TiLobe®-compatible Titanium Blank

IMPLANT PROSTHETICS



Keystone TiLobe 5.0mm

Item #	Product
9005945	5.0mm Preat Choice Titanium Scan Body
9005900	5.0mm Analog
9005900-10	5.0mm Analog (10-Pack)
9005900-25	5.0mm Analog (25-Pack)
9005901	5.0mm Open-Tray Transfer
9005903	5.0mm Digital Analog
9005903-10	5.0mm Digital Analog (10-Pack)
9005904	5.0mm Closed-Tray Transfer
9005921	5.0mm X 2mm Healing Abutment
9005922	5.0mm X 2mm Healing Abutment (Wide)
9005923	5.0mm X 4mm Healing Abutment
9005924	5.0mm X 4mm Healing Abutment (Wide)
9005930	5.0mm Engaging Verification Cylinder
9005930-10	5.0mm Engaging Verification Cylinder (10-Pack)
9005935	5.0mm Non-Engaging Verification Cylinder
9005935-10	5.0mm Non-Engaging Verification Cylinder (10-Pack)
9005990	ASC 5.0mm x 9mm Engaging Ti Base
9005995	ASC 5.0mm x 9mm Non-Engaging Ti Base
9005959	5.0mm Non-Engaging Ti Base x 4.5mm
9005961	5.0mm Engaging Ti Base x 4.5mm
9005962	5.0mm Engaging Ti Base x 6mm
9005965	5.0mm ASC Keystone TiLobe®-compatible Titanium Blank
9005967	5.0mm Keystone TiLobe®-compatible Titanium Blank

Keystone TiLobe 5.0mm Compatible

Hi-Tec® Implant Systems

Hi-Tec® Implant Systems	
Hi-Tec Tapered 3.5mm	<i>Please see Zimmer TSV 3.5mm</i>
Hi-Tec Tapered 4.5mm	<i>Please see Zimmer TSV 4.5mm</i>
Hi-Tec Vision 3.5mm	<i>Please see Nobel Tri-Lobe NP</i>
Hi-Tec Vision 4.3mm	<i>Please see Nobel Tri-Lobe RP</i>
Hi-Tec Vision 5.0mm	<i>Please see Nobel Tri-Lobe WP</i>
Hi-Tec Logic 3.5mm	<i>Please see Nobel Active / Conical NP</i>
Hi-Tec Logic 4.3 / 5.0mm	<i>Please see Nobel Active / Conical RP</i>
Hi-Tec Tapered Tite Fit 4.8mm	<i>Please see Straumann Tissue Level RN</i>
Hi-Tec Implex / Summitt 4.1mm	<i>Please see Nobel Branemark RP</i>

Imtec® Implant Systems

Imtec® Implant Systems	
Imtec 3.40 / 3.75 / 4.0 External Hex	<i>Please see Nobel Branemark RP</i>

IMZ Implant Systems

IMZ Implant Systems	
IMZ 3.3 / 4.0 External Hex	<i>Please see Nobel Branemark RP</i>
Innova Endopore 4.1mm	<i>Please see Nobel Branemark RP</i>

Lifecore Implant Systems

Lifecore Implant Systems	
Lifecore Stage One™ 3.3 / 4.1 / 4.8mm (RDS)	<i>Please see Straumann Tissue Level RN</i>

YOUR KEY TO RESTORATIVE CHOICE

The Preat Choice Library was built on the concept of giving you restorative choice from a single **scan body**. Now it's even easier to access our eighteen different options for fixture-level cases and nine restorative solutions for multi-unit level cases with the new titanium scan body. Whether you are scanning intraorally or in the laboratory, the scan body fits within your workflow and **allows you to design the best possible restoration for the case.**

TITANIUM SCAN BODY

MATTE FINISH

Easy scanning with no reflection

TITANIUM

Radiopaque to verify fit intraorally


LASER ETCHED

Clearly identify the system and size plus easily match with the product description in the CAD software



DESIGN RESTORATIONS WITH

3shape  exocad



**Compatible with
MegaGen® AnyRidge®**



8004418



8004421



8004405



8004404



9004445



9004400



9004403



9004435



9004495



9004490

MegaGen AnyRidge

Item #	Product
8004418	Healing Abutment 5 x 3mm
8004419	Healing Abutment 5 x 5mm
8004421	Healing Abutment 6 x 3mm
8004423	Healing Abutment 6 x 5mm
8004405	Closed Tray Impression Coping 4.5 x 12mm
8004404	Closed Tray Impression Coping 4.5 x 9mm
9004445	Preat Choice Titanium Scan Body
9004400	Analog
9004400-10	Analog (10-Pack)
9004400-25	Analog (25-Pack)
9004403	Digital Analog
9004403-10	Digital Analog (10-Pack)
9004435	Non-Engaging Verification Cylinder
9004435-10	Non-Engaging Verification Cylinder (10-Pack)
9004430	Engaging Verification Cylinder
9004430-10	Engaging Verification Cylinder (10-Pack)
9004490	ASC x 9mm Non-Engaging Ti Base
9004495	ASC x 9mm Engaging Ti Base
9004453	ASC MegaGen AnyRidge®-compatible Titanium Screw
9004459	Non-Engaging Ti Base x 4.5mm

MegaGen AnyRidge Compatible

IMPLANT PROSTHETICS




9004461




9004450




9004451




8004481



8004485



8004486



8004488




8004475



8004470



8004472



8004473



8004410



8004415



9004465

MegaGen AnyRidge (cont'd)

Item #	Product
9004461	Engaging Ti Base x 4.5mm
9004462	Engaging Ti Base x 6mm
9004450	Titanium Screw
9004450-10	Titanium Screw (10-Pack)
9004451	Guide Pin
8004489	Multi Unit Titanium Screw
8004480	17° Multi-Unit Abutment x 3.5mm
8004481	17° Multi-Unit Abutment x 4.5mm
8004482	17° Multi-Unit Abutment x 5.5mm
8004485	30° Multi-Unit Abutment x 3.5mm
8004486	30° Multi-Unit Abutment x 4.5mm
8004487	30° Multi-Unit Abutment x 5.5mm
8004488	30° Multi-Unit Abutment x 6.5mm
8004475	Straight Multi-Unit Abutment x 2mm
8004476	Straight Multi-Unit Abutment x 3mm
8004477	Straight Multi-Unit Abutment x 4mm
8004478	Straight Multi-Unit Abutment x 5mm
8004470	Esthetic Abutment 15° Angle, Anterior
8004472	Esthetic Abutment 15° Angle, Posterior
8004471	Esthetic Abutment Straight, Anterior
8004473	Esthetic Abutment Straight, Posterior
8004410	Engaging Plastic UCLA
8004415	Non-Engaging Plastic UCLA
9004465	ASC MegaGen AnyRidge®-compatible Titanium Blank
9004467	MegaGen AnyRidge®-compatible Titanium Blank

MegaGen AnyRidge Compatible

MegaGen Implant Systems


MegaGen Implant Systems	
MegaGen Anyone RP (Internal	<i>Please see Hiossen Regular</i>
MegaGen ExFeel External Hex	<i>Please see Nobel Branemark RP</i>
MegaGen ExFeel 4.8	<i>Please see Straumann Tissue Level RN</i>

MIS Implant Systems

MIS Implant Systems	
MIS C4 3.75 / 4.20mm	<i>Please see Zimmer TSV 3.5mm</i>
MIS SEVEN 3.75 / 4.20mm	<i>Please see Zimmer TSV 3.5mm</i>
MIS Lance External Hex 4.1mm	<i>Please see Nobel Branemark RP</i>
MIS Mistral Internal Morse Taper One Stage RN	<i>Please see Straumann Tissue Level RN</i>

IMPLANT PROSTHETICS

PREAT Multi-Unit Components are compatible with industry-standard Multi-Unit prosthetic connection and abutments.



Multi-Unit



8003706



9003701



9003702



9003743



9003700





9003703

Multi-Unit Components

Item #	Product
8003706	Multi-Unit Temporary Healing Cap
9003701	Multi-Unit Open-Tray Impression Coping
9003702	Multi-Unit Closed-Tray Impression Coping
9003743	Multi-Unit Preat Choice Titanium Scan Body w/.048 Screw
9003742	Multi-Unit Preat Choice Titanium Scan Body w/SRA Screw
9003741	Multi-Unit Preat Choice Titanium Scan Body w/.050 Screw
9003740	Multi-Unit Preat Choice Titanium Scan Body w/Unigrip Screw
9003700	Multi-Unit Abutment Analog
9003700-10	Multi-Unit Abutment Analog (10-Pack)
9003700-25	Multi-Unit Abutment Analog (25-Pack)
9003700-100	Multi-Unit Abutment Analog (100-Pack)
9003703	Multi-Unit Abutment Digital Analog
9003703-10	Multi-Unit Abutment Digital Analog (10-Pack)
9003730	Multi-Unit Verification Cylinder w/.050 Screw
9003730-10	Multi-Unit Verification Cylinder w/.050 Screw (10-Pack)
9003730-50	Multi-Unit Verification Cylinder w/.050 Screw (50-Pack)
9003730-100	Multi-Unit Verification Cylinder w/.050 Screw (100-Pack)

Multi-Unit Compatible



Multi-Unit Components (cont'd)

Item #	Product
9003731	Multi-Unit Verification Cylinder w/.048 Screw
9003731-10	Multi-Unit Verification Cylinder w/.048 Screw (10-Pack)
9003731-50	Multi-Unit Verification Cylinder w/.048 Screw (50 Pack)
9003731-100	Multi-Unit Verification Cylinder w/.048 Screw (100 Pack)
9003732	Multi-Unit Verification Cylinder w/Unigrip Screw
9003732-10	Multi-Unit Verification Cylinder w/Unigrip Screw (10-Pack)
9003732-50	Multi-Unit Verification Cylinder w/Unigrip Screw (50 Pack)
9003732-100	Multi-Unit Verification Cylinder w/Unigrip Screw (100 Pack)
9003733	Multi-Unit Verification Cylinder w/Neo Screw
9003733-10	Multi-Unit Verification Cylinder w/Neo Screw (10-Pack)
9003733-50	Multi-Unit Verification Cylinder w/Neo Screw (50-Pack)
9003733-100	Multi-Unit Verification Cylinder w/Neo Screw (100-Pack)
9003736	Multi-Unit Verification Cylinder w/SRA Screw
9003736-10	Multi-Unit Verification Cylinder w/SRA Screw (10-Pack)
9003736-50	Multi-Unit Verification Cylinder w/SRA Screw (50 Pack)
9003736-100	Multi-Unit Verification Cylinder w/SRA Screw (100 Pack)
9003795	ASC Multi-Unit x 9mm Non-Engaging Ti Base
9003752	ASC Multi-Unit Titanium Screw
9003762	Multi-Unit Ti Base w/.048 Screw
9003760	Multi-Unit Ti Base w/.050 Screw
9003761	Multi-Unit Ti Base w/SRA Screw
9003763	Multi-Unit Ti Base w/Unigrip Screw
9003768	Multi-Unit Ti Base w/Neo Screw
9003739	Multi-Unit Slim Tall Ti Base w/SRA Screw

IMPLANT PROSTHETICS



9003738



9003750



9003753



9000037



8003725



8003720



8003710



9003705

Multi-Unit Components (cont'd)

Item #	Product
9003735	Multi-Unit Slim Tall Ti Base w/.048 Screw
9003738	Multi-Unit Slim Tall Ti Base w/.050 Hex Screw
9003737	Multi-Unit Slim Tall Ti Base w/Unigrip Screw
9003757	Multi-Unit Titanium Screw (.048 Hex)
9003757-10	Multi-Unit Titanium Screw (.048 Hex)(10-Pack)
9003751	Multi-Unit Titanium Screw (.050 Hex)
9003751-10	Multi-Unit Titanium Screw (.050 Hex)(10-Pack)
9003750	Multi-Unit Titanium Screw (Unigrip)
9003750-10	Multi-Unit Titanium Screw (Unigrip)(10-Pack)
9003758	Multi-Unit Screw (Neo)
9003758-10	Multi-Unit Screw (Neo)(10-Pack)
9003753	Multi-Unit Guide Pin
9000037	Multi-Unit Driver
8003725	Multi-Unit Gold Bar UCLA Coping
8003720	Multi-Unit Gold / Plastic UCLA Coping
8003710	Multi-Unit Abutment Plastic UCLA Coping
9003704	Bidra Conversion Smart Polishing Cap
9003704-4	Bidra Conversion Smart Polishing Cap (4-pack)
9003705	Bidra Smart Polishing Cap Short
9003705-4	Bidra Smart Polishing Cap Short (4-pack)


Multi-Unit Compatible

THE BIDRA CONVERSION SMART POLISHING CAP

With an extended metal rim, **the Bidra Conversion Smart Polishing Cap fully protects both the inside and outside of the metal coping**, preventing damage to the coping when setting the acrylic, and while trimming, polishing and finishing the prosthesis.

External rough areas attract bacteria and may lead to soft tissue inflammation and bone loss around the implants.

Existing polishing caps only protect the inside of the coping, not the outside.




129

Protect Your Next Full Arch Implant Conversion



IMPLANT PROSTHETICS



Compatible with
Neodent® Helix™ GM



Neodent® Helix™ GM

Item #	Product
9007121	2mm Healing Abutment
9007122	2mm Healing Abutment (Wide)
9007123	4mm Healing Abutment
9007124	4mm Healing Abutment (Wide)
9007101	Open-Tray Transfer
9007104	Closed-Tray Transfer
9007152	Closed-Tray Pin
9007145	Preat Choice Titanium Scan Body
9007100	Analog
9007100-10	Analog (10-Pack)
9007100-25	Analog (25-Pack)
9007103	Digital Analog
9007103-10	Digital Analog (10-Pack)
9007130	Engaging Verification Cylinder


Neodent® Helix™ GM Compatible




Neodent® Helix™ GM (cont'd)

Item #	Product
9007130-10	Engaging Verification Cylinder (10-Pack)
9007135	Non-Engaging Verification Cylinder
9007135-10	Non-Engaging Verification Cylinder (10-Pack)
9007190	ASC 9mm Engaging Ti Base
9007195	ASC 9mm Non-Engaging Ti Base
9007153	ASC Titanium Screw
9007159	4.5mm Non-Engaging Ti Base
9007161	4.5mm Engaging Ti Base
9007162	6mm Engaging Ti Base
9007150	Titanium Screw
9007150-10	Titanium Screw (10-Pack)
9007151	Guide Pin
9007165	ASC Neodent® GM™-compatible Titanium Blank
9007167	Neodent® GM™-compatible Titanium Blank

IMPLANT PROSTHETICS



Compatible with
Neoss®



Compatible with
NobelActive™



9004701



9004704



9004745



9004700

Neoss® Regular Platform 3.5 / 4.0 / 4.5 / 5.0 / 5.5 mm



Item #	Product
9004145	3.5-6.0mm Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9004100	Analog
9004100-10	10 Analog (10-Pack)
9004100-25	25 Analog (25-Pack)
9004103	Digital Analog
9004153	ASC 3.5mm-6.0mm Titanium Screw
9004165	ASC Neoss®-compatible 3.5-6.0mm Titanium Blank
9004167	Neoss®-compatible 3.5-6.0mm Titanium Blank

NobelActive / Conical Connection 3.0mm

Item #	Product
9004701	3.0mm Open-Tray Transfer Nobel Biocare®
9004704	3.0mm Closed Tray Impression Coping Nobel Biocare®
9004745	3.0mm Preat Choice Titanium Scan Body
9004700	3.0mm Analog
9004700-10	3.0mm Analog (10-Pack)
9004700-25	3.0mm Analog (25-Pack)
9004703	3.0mm Digital Analog
9004721	3.0mm X 2mm Healing Abutment
9004722	3.0mm X 2mm Healing Abutment (Wide)
9004723	3.0mm X 4mm Healing Abutment
9004724	3.0mm X 4mm Healing Abutment (Wide)
8004730	3.0mm Engaging Verification Cylinder
8004730-10	3.0mm Engaging Verification Cylinder (10-Pack)
8004735	3.0mm Non-Engaging Verification Cylinder
8004715	3.0mm Non-Engaging Plastic UCLA

Neoss® Regular Platform 3.5 / 4.0 / 4.5 / 5.0 / 5.5 mm Compatible

NobelActive / Conical Connection 3.0mm Compatible



NobelActive / Conical Connection 3.0mm (cont'd)

Item #	Product
9004735-10	3.0mm Non-Engaging Verification Cylinder (10-Pack)
9004790	ASC 9.0mm Engaging Ti Base
9004795	ASC 9.0mm Non-Engaging Ti Base
004753	ASC 3.0mm Screw
9004759	3.0mm Non-Engaging Ti Base x 4.5mm
9004761	3.0mm Engaging Ti Base x 4.5mm
9004762	3.0mm Engaging Ti Base x 6.0mm
9004750	3.0mm Titanium Screw
9004750-10	3.0mm Titanium Screw (10-Pack)
9004751	3.0mm Guide Pin 20mm
9004752	3.0mm Closed-Tray Pin
8004786	3.0mm Straight Multi-Unit Abutment x 2.5mm
8004787	3.0mm Straight Multi-Unit Abutment x 3.5mm
8004770	3.0mm Esthetic Abutment 15° Angle, Anterior
8004771	3.0mm Esthetic Abutment Straight, Anterior
8004720	3.0mm Engaging Gold / Plastic UCLA Abutment
8004710	3.0mm Engaging Plastic UCLA
8004725	3.0mm Non-Engaging Gold / Plastic UCLA Abutment
8004715	3.0mm Non-Engaging Plastic UCLA
9004765	3.0mm ASC NobelActive™/Conical-compatible Titanium Blank
9004767	3.0mm NobelActive™/Conical-compatible Titanium Blank

IMPLANT PROSTHETICS



9000845



8000804



8000801




9000890



9000800



9000803



8000819



9000830



9000835



9000895



9000853



9000895





9000853

NobelActive / Conical Connection NP 3.5mm

Item #	Product
8000818	Healing Abutment 4 x 3mm
8000819	Healing Abutment 4 x 5mm
9000801	Open Tray Transfer
8000804	Closed Tray Transfer
9000845	Preat Choice Titanium Scan Body
9000800	Analog
9000800-10	Analog (10-Pack)
9000800-25	Analog (25-Pack)
9000803	Digital Analog
9000803-10	Digital Analog (10-Pack)
9000830	Engaging Verification Cylinder
9000830-10	Engaging Verification Cylinder (10-Pack)
9000835	Non-Engaging Verification Cylinder
9000835-10	Non-Engaging Verification Cylinder (10-Pack)
9000890	ASC 9mm Engaging Ti Base
9000895	ASC 9mm Non-Engaging Ti Base
9000853	ASC Titanium Screw
9000859	Non-Engaging Ti Base x 4.5mm
9000861	Engaging Ti Base x 4.5mm

NobelActive / Conical Connection NP 3.5mm Compatible



NobelActive / Conical Connection NP 3.5mm (cont'd)

Item #	Product
9000862	Engaging Ti Base x 6mm
9000850	Titanium Screw
9000850-10	Titanium Screw (10-Pack)
9000851	Guide Pin 20mm
8000889	Multi-Unit Titanium Screw
8000880	17° Multi-Unit Abutment x 2.5mm
8000881	17° Multi-Unit Abutment x 3.5mm
8000885	30° Multi-Unit Abutment x 3.5mm
8000886	30° Multi-Unit Abutment x 4.5mm
8000875	Straight Multi-Unit Abutment x 1.5mm
8000876	Straight Multi-Unit Abutment x 2.5mm
8000877	Straight Multi-Unit Abutment x 3.5mm
8000870	Esthetic Abutment 15° Angle, Anterior
8000871	Esthetic Abutment Straight, Anterior
8000825	Engaging Gold / Plastic UCLA Abutment
8000810	Engaging Plastic UCLA
8000820	Non-Engaging Gold / Plastic UCLA Abutment
8000815	Non-Engaging Plastic UCLA
9000865	ASC NobelActive™/Conical-compatible NP Titanium Blank
9000867	NobelActive™/Conical-compatible NP Titanium Blank

IMPLANT PROSTHETICS



9000945



8000901



8000822



9000900



9000903



9000930



9000935



9000990



9000995



9000962





9000950

NobelActive / Conical Connection RP 4.3 / 5.0mm

Item #	Product
8000822	Healing Abutment 5 x 3mm
8000823	Healing Abutment 5 x 5mm
8000901	Open Tray Transfer
8000904	Closed Tray Transfer
9000945	Preat Choice Titanium Scan Body
9000900	Analog
9000900-10	Analog (10-Pack)
9000900-25	Analog (25-Pack)
9000903	Digital Analog
9000903-10	Digital Analog (10-Pack)
9000930	Engaging Verification Cylinder
9000930-10	Engaging Verification Cylinder (10-Pack)
9000935	Non-Engaging Verification Cylinder
9000935-10	Non-Engaging Verification Cylinder (10-Pack)
9000990	ASC x 9mm Engaging Ti Base
9000995	ASC x 9mm Non-Engaging Ti Base
9000953	ASC Titanium Screw
9000959	Non-Engaging Ti Base x 4.5mm
9000961	Engaging Ti Base x 4.5mm
9000962	Engaging Ti Base x 6mm

NobelActive / Conical Connection RP 4.3 / 5.0mm Compatible



NobelActive / Conical Connection RP 4.3 / 5.0mm (cont'd)

Item #	Product
9000950	Titanium Screw
9000950-10	Titanium Screw (10-Pack)
9000951	Guide Pin 20mm
8000989	Multi-Unit Titanium Screw
8000980	17° Multi-Unit Abutment x 2.5mm
8000981	17° Multi-Unit Abutment x 3.5mm
8000985	30° Multi-Unit Abutment x 3.5mm
8000986	30° Multi-Unit Abutment x 4.5mm
8000975	Straight Multi-Unit Abutment x 1.5mm
8000976	Straight Multi-Unit Abutment x 2.5mm
8000977	Straight Multi-Unit Abutment x 3.5mm
8000978	Straight Multi-Unit Abutment x 4.5mm
8000970	Esthetic Abutment 15° Angle, Anterior
8000972	Esthetic Abutment 15° Angle, Posterior
8000971	Esthetic Abutment Straight, Anterior
8000973	Esthetic Abutment Straight, Posterior
8000920	Engaging Gold / Plastic UCLA Abutment
9000910	Engaging Plastic UCLA
8000925	Non-Engaging Gold / Plastic UCLA Abutment
8000915	Non-Engaging Plastic UCLA
9000965	ASC NobelActive™/Conical-compatible RP Titanium Blank
9000967	NobelActive™/Conical-compatible RP Titanium Blank

IMPLANT PROSTHETICS



8004904



9004903



9004900



8004930



8004935



8004910



8004962



8004915




8004925

NobelActive / Conical Connection WP

Item #	Product
8004901	Open Tray Transfer
8004904	Closed Tray Impression Coping
8004941	Scan Body <i>(for use with Ti Bases)</i>
9004900	Analog
9004900-10	Analog (10-Pack)
9004900-25	Analog (25-Pack)
9004903	Digital Analog
9004903-10	Digital Analog (10-Pack)
8004935	Non-Engaging Verification Cylinder
8004935-10	Non-Engaging Verification Cylinder (10-Pack)
8004930	Engaging Verification Cylinder
8004930-10	Engaging Verification Cylinder (10-Pack)
8004959	Non-Engaging Ti Base x 4.5mm
8004961	Engaging Ti Base x 4.5mm
8004962	Engaging Ti Base x 6mm
8004972	Esthetic Abutment Straight, Posterior
8004920	Engaging Gold / Plastic UCLA Abutment
8004910	Engaging Plastic UCLA
8004925	Non-Engaging Gold / Plastic UCLA Abutment
8004915	Non-Engaging Plastic UCLA

NobelActive / Conical Connection WP Compatible



Compatible with
NobelBiocare™ Tri-Lobe



NobelBiocare Tri-Lobe NP

Item #	Product
8000401	Open-Tray Transfer 3.5mm
8000404	Closed-Tray Transfer 3.5mm
9000445	Preat Choice Titanium Scan Body
9000400	Analog
9000400-10	Analog (10-Pack)
9000400-25	Analog (25-Pack)
9000403	Digital Analog
9000403-10	Digital Analog (10-Pack)
9000430	Engaging Verification Cylinder
9000430-10	Engaging Verification Cylinder (10-Pack)
9000435	Non-Engaging Verification Cylinder
9000435-10	Non-Engaging Verification Cylinder (10-Pack)
9000453	ASC Titanium Screw
9000459	Non-Engaging Ti Base x 4.5mm
9000461	Engaging Ti Base x 4.5mm

NobelBiocare Tri-Lobe NP Compatible

IMPLANT PROSTHETICS



9000462




9000453



9000459



9000451



8000481



8000476



9000467



8000410-10



8000415



9000465

NobelBiocare Tri-Lobe NP (cont'd)

Item #	Product
9000462	Engaging Ti Base x 6mm
9000450	Titanium Screw
9000450-10	Titanium Screw (10-Pack)
9000451	Guide Pin 20mm
8000480	17° Multi-Unit Abutment x 2mm
8000481	17° Multi-Unit Abutment x 3mm
8000475	Straight Multi-Unit Abutment x 1mm
8000476	Straight Multi-Unit Abutment x 2mm
8000477	Straight Multi-Unit Abutment x 3mm
8000470	Esthetic Abutment 15° Angle, Anterior
8000471	Esthetic Abutment Straight, Anterior
8000420	Engaging Gold / Plastic UCLA Abutment
8000410	Engaging Plastic UCLA
8000415	Non-Engaging Plastic UCLA
9000765	ASC NobelBiocare™ Tri-Lobe-compatible NP Titanium Blank
9000767	NobelBiocare™ Tri-Lobe-compatible NP Titanium Blank


NobelBiocare Tri-Lobe NP Compatible



NobelBiocare Tri-Lobe RP

Item #	Product
8000501	Open-Tray Transfer 4.3mm
8000504	Closed-Tray Transfer 4.3mm
9000545	Preat Choice Titanium Scan Body
9000500	Analog
9000500-10	Analog (10-Pack)
9000500-25	Analog (25-Pack)
9000503	Digital Analog
9000503-10	Digital Analog (10-Pack)
9000530	Engaging Verification Cylinder
9000530-10	Engaging Verification Cylinder (10-Pack)
9000535	Non-Engaging Verification Cylinder
9000535-10	Non-Engaging Verification Cylinder (10-Pack)
9000590	ASC x 9mm Engaging Ti Base <i>(not currently for sale)</i>
9000595	ASC x 9mm Non-Engaging Ti Base <i>(not currently for sale)</i>
9000553	ASC Titanium Screw
9000559	Non-Engaging Ti Base x 4.5mm
9000561	Engaging Ti Base x 4.5mm
9000562	Engaging Ti Base x 6mm
9000550	Titanium Screw

IMPLANT PROSTHETICS




8000576



9000551



9000550




8000587



8000572



8000571



8000573



8000520



8000510



8000515




9000567

NobelBiocare Tri-Lobe RP (cont'd)

Item #	Product
9000550-10	Titanium Screw (10-pack)
9000551	Guide Pin 20mm
8000589	Multi-Unit Titanium Screw
8000575	Straight Multi-Unit Abutment x 1mm
8000576	Straight Multi-Unit Abutment x 2mm
8000577	Straight Multi-Unit Abutment x 3mm
8000579	Straight Multi-Unit Abutment x 5mm
8000580	17° Multi-Unit Abutment x 2mm
8000581	17° Multi-Unit Abutment x 3mm
8000582	17° Multi-Unit Abutment x 4mm
8000586	30° Multi-Unit Abutment x 4mm
8000587	30° Multi-Unit Abutment x 5mm
8000570	Esthetic Abutment 15° Angle, Anterior
8000572	Esthetic Abutment 15° Angle, Posterior
8000571	Esthetic Abutment Straight, Anterior
8000573	Esthetic Abutment Straight, Posterior
8000520	Engaging Gold / Plastic UCLA Abutment
8000510	Engaging Plastic UCLA
8000515	Non-Engaging Plastic UCLA
9000565	ASC NobelBiocare™ Tri-Lobe-compatible RP Titanium Blank
9000567	NobelBiocare™ Tri-Lobe-compatible RP Titanium Blank


NobelBiocare Tri-Lobe RP Compatible




NobelBiocare Tri-Lobe WP

Item #	Product
8000601	Open-Tray Transfer 5.0mm
8000604	Closed-Tray Transfer 5.0mm
9000645	Preat Choice Titanium Scan Body
9000600	Analog
9000600-10	Analog (10-Pack)
9000600-25	Analog (25-Pack)
9000603	Digital Analog
9000603-10	Digital Analog (10-Pack)
9000630	Engaging Verification Cylinder
9000630-10	Engaging Verification Cylinder (10-Pack)
9000635	Non-Engaging Verification Cylinder
9000635-10	Non-Engaging Verification Cylinder (10-Pack)
9000659	Non-Engaging Ti Base x 4.5mm
9000553	ASC Titanium Screw
9000690	ASC x 9mm Engaging Ti Base <i>(not currently for sale)</i>
9000695	ASC x 9mm Non-Engaging Ti Base <i>(not currently for sale)</i>
9000661	Engaging Ti Base x 4.5mm
9000662	Engaging Ti Base x 6mm
	Screws and Guide Pins - See page 132 & 133

IMPLANT PROSTHETICS



8000675



8000677



9000665

NobelBiocare Tri-Lobe WP (cont'd)

Item #	Product
8000675	Straight Multi-Unit Abutment x 1mm
8000676	Straight Multi-Unit Abutment x 2mm
8000677	Straight Multi-Unit Abutment x 3mm
8000670	Esthetic Abutment 15° Angle, Anterior
8000672	Esthetic Abutment 15° Angle, Posterior
8000671	Esthetic Abutment Straight, Anterior
8000673	Esthetic Abutment Straight, Posterior
8000620	Engaging Gold / Plastic UCLA Abutment
8000610	Engaging Plastic UCLA
8000615	Non-Engaging Plastic UCLA
9000665	ASC NobelBiocare™ Tri-Lobe-compatible WP Titanium Blank
9000667	NobelBiocare™ Tri-Lobe-compatible WP Titanium Blank

NobelBiocare Tri-Lobe WP Compatible



NobelBiocare Tri-Lobe 6.0mm

Item #	Product
8000741	6.0mm Titanium Scan Body (for use with Ti Bases)
9000745	Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9000700	6.0mm Analog
9000700-10	6.0mm Analog (10-Pack)
9000700-25	6.0mm Analog (25-Pack)
9000703	6.0mm Digital Analog
8000730	6.0mm Engaging Verification Cylinder
8000735	6.0mm Non-Engaging Verification Cylinder
8000759	6.0mm Non-Engaging Ti Base x 4.5mm
8000761	6.0mm Engaging Ti Base x 4.5mm
8000762	6.0mm Engaging Ti Base x 6mm
8000773	6.0mm Esthetic Abutment Straight, Posterior
8000720	6.0mm Engaging Gold / Plastic UCLA Abutment
8000710	6.0mm Engaging Plastic UCLA
8000715	6.0mm Non-Engaging Plastic UCLA
9000765	ASC NobelBiocare™ Tri-Lobe-compatible 6.0mm Titanium Blank
9000767	NobelBiocare™ Tri-Lobe-compatible 6.0mm Titanium Blank

IMPLANT PROSTHETICS

NobelBiocare Implant Systems

Steri-Oss Hex Lock 3.8 / 4.5mm	Please see Nobel Branemark RP
--------------------------------	-------------------------------

Southern Implant Systems

Southern 12 Degrees 4.3 (3.5mmD)	Please see Nobel Tri-Lobe NP
Southern 12 Degrees 5.0 (4.3mmD)	Please see Nobel Tri-Lobe RP
Southern 12 Degrees 7.0 (5.0mmD)	Please see Nobel Tri-Lobe WP
Southern CoAxis / MAX-IT 4.0 / 5.0mm (4.8mmD)	Please see Straumann Tissue Level RN
Southern CoAxis / MAX-IT 6.0 / (6.5mmD)	Please see Straumann Tissue Level WN
Southern CoAxis TriNeck / TriMax 6.0	Please see Nobel Tri-Lobe 6.0
Southern CoAxis TriNeck / TriMax NP	Please see Nobel Tri-Lobe NP
Southern CoAxis TriNeck / TriMax RP	Please see Nobel Tri-Lobe RP
Southern CoAxis TriNeck / TriMax WP	Please see Nobel Tri-Lobe WP


Sterngold Implant Systems

Sterngold IC 3.3 / 4.1 / 4.8mm	Please see Straumann Tissue Level RN
Sterngold PUR NP	Please see Zimmer TSV 3.5
Sterngold PUR RP	Please see Zimmer TSV 4.5
Sterngold TRU NP	Please see NobelActive / Conical NP
Sterngold TRU RP	Please see NobelActive / Conical RP
Hex Cylinder 3.3 / 4.0mm	Please see Nobel Branemark RP
Hex Screws 3.75 / 4.0 / 5.0mm	Please see Nobel Branemark RP

NobelBiocare Implant System

Southern Implant System

Sterngold Implant System




**Compatible with
Straumann® BLX**

Straumann® BLX RB/WB

Item #	Product
9007021	2mm Healing Abutment
9007022	2mm Healing Abutment (Wide)
9007023	4mm Healing Abutment
9007024	4mm Healing Abutment (Wide)
9007001	Open-Tray Transfer
9007004	Closed-Tray Transfer
9007045	Preat Choice Titanium Scan Body
9007000	Analog
9007000-10	Analog (10-Pack)
9007000-25	Analog (25-Pack)
9007003	Digital Analog
9007003-10	Digital Analog (10-Pack)
9007030	Engaging Verification Cylinder
9007030-10	Engaging Verification Cylinder (10-Pack)
9007035	Non-Engaging Verification Cylinder
9007035-10	Non-Engaging Verification Cylinder (10-Pack)
9007090	ASC 9mm Engaging Ti Base
9007095	ASC 9mm Non-Engaging Ti Base


Straumann® BLX RB/WB Compatible





Straumann® BLX RB/WB (cont'd)	
Item #	Product
9007053	ASC Titanium Screw
9007059	Non-Engaging Ti Base x 4.5mm
9007061	Engaging Ti Base x 4.5mm
9007062	Engaging Ti Base x 6mm
9007050	Titanium Screw
9007050-10	Titanium Screw (10-Pack)
9007051	Guide Pin
9007065	ASC Straumann®™ BLX-compatible RB/WB Titanium Blank
9007067	Straumann® BLX-compatible RB/WB Titanium Blank


Straumann® BLX RB/WB Compatible




Straumann® BLX WB

Item #	Product
9007221	Healing Abutment x 2mm
9007223	Healing Abutment x 4mm
9007201	Open-Tray Transfer
9007204	Closed-Tray Transfer
9007245	Preat Choice Titanium Scan Body
9007200	Analog
9007200-10	Analog (10-Pack)
9007200-25	Analog (25-Pack)
9007203	Digital Analog
9007203-10	Digital Analog (10-Pack)
9007230	Engaging Verification Cylinder
9007235	Non-Engaging Verification Cylinder
9007259	Non-Engaging Ti Base x 4.5mm
9007261	Engaging Ti Base x 4.5mm
9007262	Engaging Ti Base x 6mm
9007290	ASC x 9mm Engaging Ti Base
9007295	ASC x 9mm Non-Engaging Ti Base
9007230-10	Eng. Verification Cylinder (10-Pack)
9007235-10	Non-Eng. Verification Cylinder (10-Pack)
9007265	ASC Straumann® BLX-compatible WB Titanium Blank
9007267	Straumann® BLX-compatible WB Titanium Blank

IMPLANT PROSTHETICS



**Compatible with
Straumann® Bone Level**




Straumann® Bone Level NC (Narrow Crossfit, 3.3mm)

Item #	Product
8001521	Healing Abutment 4.8 x 3mm
8001523	Healing Abutment 4.8 x 5mm
8001504	Closed Tray Impression Tray Coping Slim
9001545	Preat Choice Titanium Scan Body
9001500	Analog
9001500-10	Analog (10-Pack)
9001500-25	Analog (25-Pack)
9001503	Digital Analog
9001503-10	Digital Analog (10-Pack)
9001530	Engaging Verification Cylinder
9001530-10	Engaging Verification Cylinder (10-Pack)
9001535	Non-Engaging Verification Cylinder
9001535-10	Non-Engaging Verification Cylinder (10-Pack)
9001590	ASC 9mm Engaging Ti Base
9001595	ASC 9mm Non-Engaging Ti Base

Straumann® Bone Level NC (Narrow Crossfit, 3.3mm) Compatible

150



Straumann® Bone Level NC (Narrow Crossfit, 3.3mm) (cont'd)

Item #	Product
9001553	ASC NC Titanium Screw
9001559	Non-Engaging Ti Base x 4.5mm
9001561	Engaging Ti Base x 4.5mm
9001562	Engaging Ti Base x 6.0mm
9001550	Titanium Screw
9001550-10	Titanium Screw (10-Pack)
9001551	Guide Pin 20mm
8001575	Straight Multi-Unit Abutment x 1.5mm
8001576	Straight Multi-Unit Abutment x 2.5mm
8001577	Straight Multi-Unit Abutment x 3.5mm
8001580	17° Multi-Unit Abutment x 2.5mm
8001581	17° Multi-Unit Abutment x 3.5mm
8001570	Esthetic Abutment 15° Angle, Anterior
8001571	Esthetic Abutment Straight, Anterior
8001510	Engaging Plastic UCLA
8001515	Non-Engaging Plastic UCLA
9001565	ASC Straumann® BL-compatible NC Titanium Blank
9001567	Straumann® BL-compatible NC Titanium Blank




IMPLANT PROSTHETICS



Straumann® Bone Level RC (Regular Crossfit, 4.1mm)

Item #	Product
8001421	Healing Abutment 6 x 3mm
8001423	Healing Abutment 6 x 5mm
8001404	Closed Tray Impression Tray Coping 4.5 x 9mm
9001445	Preat Choice Titanium Scan Body
9001400	Analog
9001400-10	Analog (10-Pack)
9001400-25	Analog (25-Pack)
9001403	Digital Analog
9001403-10	Digital Analog (10-Pack)
9001430	Engaging Verification Cylinder
9001430-10	Engaging Verification Cylinder (10-Pack)
9001435	Non-Engaging Verification Cylinder
9001435-10	Non-Engaging Verification Cylinder (10-Pack)
9001490	ASC 9mm Engaging Ti Base
9001495	ASC 9mm Non-Engaging Ti Base
9001453	ASC RC Titanium Screw
9001459	Non-Engaging Ti Base x 4.5mm
9001461	Engaging Ti Base x 4.5mm
9001462	Engaging Ti Base x 6mm
9001450	Titanium Screw


Straumann® Bone Level RC (Regular Crossfit, 4.1mm) Compatible



Straumann Bone Level RC (Regular Crossfit, 4.1mm) (cont'd)

Item #	Product
9001450-10	Titanium Screw (10-Pack)
9001451	Guide Pin 20mm
9003756	Straumann Screw Receiving Abutment Screw 1.4mm
9003756-10	Screw Receiving Abutment Screw 1.4mm (10-Pack)
8001489	Multi-Unit Titanium Screw
8001480	17° Multi-Unit Abutment x 2.5mm
8001481	17° Multi-Unit Abutment x 3.5mm
8001485	30° Multi-Unit Abutment x 3.5mm
8001486	30° Multi-Unit Abutment x 4.5mm
8001482	30° Multi-Unit Abutment x 4.5mm
8001475	Straight Multi-Unit Abutment x 1.5mm
8001476	Straight Multi-Unit Abutment x 2.5mm
8001477	Straight Multi-Unit Abutment x 3.5mm
8001478	Straight Multi-Unit Abutment x 4.5mm
8001470	Esthetic Abutment 15° Angle, Anterior
8001472	Esthetic Abutment 15° Angle, Posterior
8001471	Esthetic Abutment Straight, Anterior
8001473	Esthetic Abutment Straight, Posterior
9001410	Engaging Plastic UCLA
8001415	Non-Engaging Plastic UCLA
9001465	ASC Straumann® BL-compatible RC Titanium Blank
9001467	Straumann® BL-compatible RC Titanium Blank

IMPLANT PROSTHETICS



**Compatible with
Straumann® Tissue Level**



8001100-10



8001150



9001245



9001200



9001203



9001230



9001235



9001290



9001295

Straumann® Tissue Level NN (3.5mm Platform Diameter)

Item #	Product
8001100-10	Analog (10-Pack)
8001100-25	Analog (25-Pack)
8001150	Titanium Screw
8001150-10	Titanium Screw (10-Pack)

Straumann® Tissue Level RN (4.8mm Platform Diameter)

Item #	Product
9001245	Preat Choice Titanium Scan Body
9001200	Analog
9001200-10	Analog (10-Pack)
9001200-25	Analog (25-Pack)
9001203	Digital Analog
9001203-10	Digital Analog (10-Pack)
9001230	Engaging Verification Cylinder
9001230-10	Engaging Verification Cylinder (10-Pack)
9001235	Non-Engaging Verification Cylinder
9001235-10	Non-Engaging Verification Cylinder (10-Pack)
9001290	ASC 9mm Engaging Ti Base
9001295	ASC 9mm Non-Engaging Ti Base

Straumann® Tissue Level NN (3.5mm Platform Diameter) Compatible

Straumann® Tissue Level RN
(4.8mm Platform Diameter) Compatible



Straumann® Tissue Level RN (4.8mm Platform Diameter) (cont'd)

Item #	Product
9001253	ASC RN/WN Titanium Screw
9001261	Engaging Ti Base x 4.5mm
9001259	Non-Engaging Ti Base x 4.5mm
9001262	Engaging Ti Base x 6mm
9001250	RN / WN Titanium Screw
9001250-10	RN / WN Titanium Screw (10-Pack)
9001251	Guide Pin
8001220	Engaging Gold / Plastic UCLA Abutment
8001225	Non-Engaging Gold / Plastic UCLA Abutment
9001265	ASC Straumann® TL-compatible RN Titanium Blank
9001267	Straumann® TL-compatible RN Titanium Blank



9001345



9001300




8001337



8001331



8001359



8001361



8001320




8001325


Straumann Tissue Level WN (6.5mm Platform Diameter)

Item #	Product
8001341	Scan Body (for use with Ti Bases)
9001345	Preat Choice Titanium Scan Body (for Preat Custom Solutions)
9001300	Analog
9001300-10	Analog (10-Pack)
9001300-25	Analog (25-Pack)
9001303	Digital Analog
8001337	Non-Engaging Verification Cylinder
8001337-10	Non-Engaging Verification Cylinder (10-Pack)
8001331	Engaging Verification Cylinder
8001331-10	Engaging Verification Cylinder (10-Pack)
8001359	Non-Engaging Ti Base x 4.5mm
8001361	Engaging Ti Base x 4.5mm
8001362	Engaging Ti Base x 6mm
8001320	Engaging Gold / Plastic UCLA Abutment
8001325	Non-Engaging Gold / Plastic UCLA Abutment

Straumann Tissue Level WN (6.5mm Platform Diameter) Compatible



**Compatible with
Zimmer® Tapered Screw-Vent®**



Zimmer Tapered Screw-Vent® 3.5mm

Item #	Product
8000109	3.5mm Healing Abutment 3.8 x 3mm
8000119	3.5mm Healing Abutment 3.8 x 5mm
8000121	3.5mm Healing Abutment 4.5 x 3mm
8000123	3.5mm Healing Abutment 4.5 x 5mm
8000101	3.5mm Open-Tray Transfer 3.5 x 10mm
8000102	3.5mm Open-Tray Transfer 4.7 x 10mm
8000103	3.5mm Open-Tray Transfer Screw
8000106	3.5mm Closed-Tray Transfer Screw
8000104	3.5mm Closed-Tray Transfer 3.5 x 8mm
8000105	3.5mm Closed-Tray Transfer 4.7 x 8mm
9000145	3.5mm Preat Choice Titanium Scan Body
9000100	3.5mm Analog
9000100-10	3.5mm Analog (10-Pack)
9000100-25	3.5mm Analog (25-Pack)
9000103	3.5mm Digital Analog
9000103-10	3.5mm Digital Analog (10-Pack)
9000130	3.5mm Engaging Verification Cylinder
9000130-10	3.5mm Engaging Verification Cylinder (10-Pack)
9000135	3.5mm Non-Engaging Verification Cylinder
9000135-10	3.5mm Non-Engaging Verification Cylinder (10-Pack)
9000190	ASC 3.5mm x 9mm Engaging Ti Base
9000195	ASC 3.5mm x 9mm Non-Engaging Ti Base

IMPLANT PROSTHETICS




9000190




9000161




9000151



8000181



8000176



8000170



8000171



9000110



8000115



9000165

Zimmer Tapered Screw-Vent® 3.5mm (cont'd)

Item #	Product
9000153	ASC 3.5mm/ 4.5mm/ 5.7mm Titanium Screw
9000159	3.5mm Non-Engaging Ti Base x 4.5mm
9000161	3.5mm Engaging Ti Base x 4.5mm
9000162	3.5mm Engaging Ti Base x 6mm
9000150	<u>3.5mm / 4.5mm / 5.7mm</u> Titanium Screw
9000150-10	<u>3.5mm / 4.5mm / 5.7mm</u> Titanium Screw (10-Pack)
9000151	<u>3.5mm / 4.5mm / 5.7mm</u> Guide Pin 20mm
8000189	<u>3.5mm / 4.5mm / 5.7mm</u> Multi Unit Titanium Screw
8000180	3.5mm 17° Multi-Unit Abutment x 2mm
8000181	3.5mm 17° Multi-Unit Abutment x 3mm
8000185	3.5mm 30° Multi-Unit Abutment x 3mm
8000175	3.5mm Straight Multi-Unit Abutment x 1mm
8000176	3.5mm Straight Multi-Unit Abutment x 2mm
8000170	3.5mm Esthetic Abutment 15° Angle, Anterior
8000171	3.5mm Esthetic Abutment Straight, Anterior
8000120	3.5mm Engaging Gold / Plastic UCLA Abutment
9000110	3.5mm Engaging Plastic UCLA
8000115	3.5mm Non-Engaging Plastic UCLA
9000165	ASC Zimmer® TSV-compatible 3.5mm Titanium Blank
9000167	Zimmer® TSV-compatible 3.5mm Titanium Blank

Zimmer Tapered Screw-Vent® 3.5mm Compatible



Zimmer Tapered Screw-Vent® 4.5mm

Item #	Product
8000221	4.5mm Healing Abutment 5.5 x 3mm
8000223	4.5mm Healing Abutment 5.5 x 5mm
8000201	4.5mm Open-Tray Transfer 4.5 x 10mm
8000202	4.5mm Open-Tray Transfer 5.7 x 10mm
8000204	4.5mm Closed-Tray Transfer 4.5 x 8mm
8000205	4.5mm Closed-Tray Transfer 5.7 x 8mm
9000245	4.5mm Preat Choice Titanium Scan Body
9000200	4.5mm Analog
9000200-10	4.5mm Analog (10-Pack)
9000200-25	4.5mm Analog (25-Pack)
9000203	4.5mm Digital Analog
9000203-10	4.5mm Digital Analog (10-Pack)
9000230	4.5mm Engaging Verification Cylinder
9000230-10	4.5mm Engaging Verification Cylinder (10-Pack)
9000235	4.5mm Non-Engaging Verification Cylinder
9000235-10	4.5mm Non-Engaging Verification Cylinder (10-Pack)
9000290	ASC 4.5mm x 9mm Engaging Ti Base
9000295	ASC 4.5mm x 9mm Non-Engaging Ti Base


IMPLANT PROSTHETICS




9000295




9000261




8000276



8000280




8000285



8000270



8000271



8000273



8000210



8000215



9000267

Zimmer Tapered Screw-Vent® 4.5mm (cont'd)

Item #	Product
9000259	4.5mm Non-Engaging Ti Base x 4.5mm
9000261	4.5mm Engaging Ti Base x 4.5mm
9000262	4.5mm Engaging Ti Base x 6mm
	Screws and Guide Pins - See page 149
8000275	4.5mm Straight Multi-Unit Abutment x 1mm
8000276	4.5mm Straight Multi-Unit Abutment x 2mm
8000277	4.5mm Straight Multi-Unit Abutment x 3mm
8000278	4.5mm Straight Multi-Unit Abutment x 5mm
8000280	4.5mm 17° Multi-Unit Abutment x 2mm
8000281	4.5mm 17° Multi-Unit Abutment x 3mm
8000283	4.5mm 17° Multi-Unit Abutment x 5mm
8000285	4.5mm 30° Multi-Unit Abutment x 3mm
8000286	4.5mm 30° Multi-Unit Abutment x 5mm
8000270	4.5mm Esthetic Abutment 15° Angle, Anterior
8000272	4.5mm Esthetic Abutment 15° Angle, Posterior
8000271	4.5mm Esthetic Abutment Straight, Anterior
8000273	4.5mm Esthetic Abutment Straight, Posterior
8000220	4.5mm Engaging Gold / Plastic UCLA Abutment
8000210	4.5mm Engaging Plastic UCLA
8000215	4.5mm Non-Engaging Plastic UCLA
9000265	ASC Zimmer® TSV-compatible 4.5mm Titanium Blank
9000267	Zimmer® TSV-compatible 4.5mm Titanium Blank

Zimmer Tapered Screw-Vent® 4.5mm Compatible



Zimmer Tapered Screw-Vent® 5.7mm

Item #	Product
9000345	5.7mm Preat Choice Titanium Scan Body
9000300	5.7mm Analog
9000300-10	5.7mm Analog (10-Pack)
9000300-25	5.7mm Analog (25-Pack)
9000303	5.7mm Digital Analog
9000303-10	5.7mm Digital Analog (10-Pack)
9000335	5.7mm Non-Engaging Verification Cylinder
9000335-10	5.7mm Non-Engaging Verification Cylinder (10-Pack)
9000330	5.7mm Engaging Verification Cylinder
9000330-10	5.7mm Engaging Verification Cylinder (10-Pack)
9000390	ASC 5.7mm x 9mm Engaging Ti Base
9000395	ASC 5.7mm x 9mm Non-Engaging Ti Base
9000359	5.7mm Non-Engaging Ti Base x 4.5mm
9000361	5.7mm Engaging Ti Base x 4.5mm
9000362	5.7mm Engaging Ti Base x 6mm
	Screws and Guide Pins - See page 149
8000310	5.7mm Engaging Plastic UCLA
8000315	5.7mm Non-Engaging Plastic UCLA
9000365	ASC Zimmer® TSV-compatible 5.7mm Titanium Blank
9000367	Zimmer® TSV-compatible 5.7mm Titanium Blank

IMPLANT BUDDY

RESTORE 99% OF IMPLANT CASES

Prosthetic Connection Included:

- Straight
- SCS
- UniGrip
- Hex 1.3mm/.050"
- Hex 1.2mm/.048"
- Hex 0.9mm/.035"
- ASC Angled Driver

TAPERED DRIVERS

for a stronger connection

SMART THUMB KNOB

inserts directly onto
torque driver

SHORT & LONG DRIVERS

for each prosthetic
screw connection

COLOR CODED

for easy ID




AUTOCLAVABLE

MAGNETIC BASE


all components stay in
place mount on any
metallic surface

IMPLANT BUDDY LT (LATCH-TYPE)

designed to work with the
Omega Adjustable Torque Wrench



Driver / Prosthetic Connection	Systems
 Straight Driver	Straight Systems
 SCS Driver	Straumann®
 UniGrip Driver	NobelBiocare™
 Hex 1.3mm/.050 inch	Lifecore, BioHorizons®, Astra, Zimmer
 Hex 1.2mm/.048 inch	Biomet 3i, MegaGen®, Hiossen®
 Hex 0.9mm/.035 inch	Friadent, Ankylos 3i Cover Screw
 Angled Driver	ASC – Angulated Screw Channel
 Square 1.3mm/.050 inch	Gold Screw/NobelBiocare™



Implant Components and Drivers

Item #	Product
9000000	Implant Buddy Driver Set
9000020	Omega Adjustable Torque Wrench: Works hand in hand with the Implant Buddy. This High Precision Calibrated Torque can adjust from 10 to 70 Ncm! Implant Buddy & Omega Wrench Set: Save when you buy the Implant Buddy Driver Set and Omega Torque Wrench Set together
4009020	30N-cm Torque and Thumb Knob Wrench Kit
4034391	20N-cm LOCATOR Torque and Thumb Knob Wrench Kit
1600035	Latch to Square Driver Insert
1600036	Square To Straumann Driver Insert
8000085	Cuff Height & Angulation Tool Set
8000080	Cuff Height Measuring Tool
9000041	Hand Driver

IMPLANT DRIVERS & TOOLS



9000030



9000032



9000033



9000034



9000035



9000036



9000037

Latch Type Drivers

Item #	Product
9000030	Implant Buddy LT - Latch-Type Driver Set
9000031	Latch-Type Straight Driver
9000032	Latch-Type SCS Driver
9000033	Latch-Type UniGrip Driver
9000034	Latch-Type .050 Hex Driver
9000035	Latch-Type .048 Hex Driver
9000036	Latch-Type .035 Hex Driver
9000037	Multi-Unit Latch Driver

Square Type Drivers

Square Type Prosthetic Connection	Short Total Length 18mm Shaft Length 11mm	Medium Total Length 23mm Shaft Length 16mm	Long Total Length 39mm Shaft Length 32mm
 Straight Driver	8000051	9000061	9000071
 SCS Driver	8000052	9000062	9000072
 UniGrip Driver	8000053	9000063	9000073
 .050 Hex Driver	8000054	9000064	9000074
 .048 Hex Driver	8000055	9000065	9000075
 .035 Hex Driver	8000056	9000066	9000076
 .050 Square Driver	8000057	8000067	8000077


THE CUFF HEIGHT AND ANGULATION SET*

The Must Have Tool for All Implant Professionals


Stop Guessing Abutment Heights & Angulations

EASY TO USE

1. Place the tool on top of the implant platform to measure gingival height.
2. Rotate angle measurement arm to identify if a straight abutment, 17° or 30° abutment is necessary.



Item #8000085



Abutments Too High?



Abutments Too Low?

See What our Customers Have to Say

"This is a must have tool for anyone doing anything implant related. I use this tool daily!"

- Sevan Pulurian, Pacific Dental

"This tool really saves me time and confusion on large implant cases. It takes away all the guesswork and permits me to plant these cases accurately from the laboratory, allowing for the proper ordering of necessary components."

- Joshua Polansky, MDC, Niche Dental Studio

"When things get tough small tools can have a big impact. Using this one almost every day. Thank you PREAT, great development!"

- Alexander Wünsche, Cdt, Zahntechnik Inc.

*Set includes two diameters of the tool, allowing use on implants from 3mm to 6mm Ø.

Implant Compatibility Chart

System	Compatible With	Page
3.5 (red), 4.0 (yellow), and 5.0mm (green)	OsseoLink Global Implant® Solutions OctaCone	
3I® Osseotite XP 3.75 & 4.1mm Standard	Nobel® Branemark® RP 3.75 / 4.1mm	94
AB Dental® 3.75mm NP	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
AB Dental® 4.50mm WP	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
ACE® Infinity Internal Hex 3.7 / 4.1mm	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
ACE® Infinity Internal Hex 4.7 / 5.1mm	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
ACE® Infinity TRI-CAM 3.5mm	Nobel® Replace Tri-Lobe NP 3.5mm	130
ACE® Infinity TRI-CAM 4.3mm	Nobel® Replace Tri-Lobe RP 4.3mm	132
ACE® Infinity TRI-CAM 5.0mm	Nobel® Replace Tri-Lobe WP 5.0mm	134
ACE® Infinity X-hex Implant 3.3mm	Nobel® Branemark NP 3.3mm	
ACE® Infinity X-hex Implant 3.7 / 4.0mm	Nobel® Branemark RP 3.75 / 4.1mm	94
ACE® Infinity X-hex Implant 4.75mm	Nobel® Branemark WP 5.1mm	
Adin® Swell	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Adin® Touareg	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Alfa Gate® Bioactive	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Alfa Gate® DFCI 3.5mm	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Almitech Implant System®	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Alpha Bio® Internal Hex	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Alpha Dent® Active / Plus	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Alpha Dent® Bio	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Alpha Dent® Classic	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Alpha Dent® Conus	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
American Dental Implant® Int. Hex Skinny 2.4 (3.5)	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
American Dental Implant® Int. Hex Small (3.5)	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
American Dental Implant® Int. Hex Wide (4.5)	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
American Dental Implant® Int. Hex XL (5.7)	Zimmer® Screw-Vent / Tapered Screw-Vent 5.7mm	152
ARDS® Classic 3.75mm	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
AS Technology® c-fix Cone Morse Regular (4.8)	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
AS Technology® e-fix Ext Hexagon-Narrow Plat. (3.5)	Straumann® SynOcta (Tissue Level) Standard NN (3.5mmD)	

All product names, logos, and brands are property of their respective owners. All company, product and service names are for identification purposes only. Use of these third-party names, logos, and brands does not imply, nor should it be understood as, any indication that the respective third-party brand owners have endorsed or sponsored, or are otherwise affiliated with or connected to, Preat or Preat's goods or services.

Implant Compatibility Chart *(cont'd)*

System	Compatible With	Page
AS Technology® e-fix Ext Hexagon-Regular Platform (4.1)	Nobel® Branemark® RP 3.75 / 4.1mm	94
AS Technology® e-fix Turque Interno-Regular Platform (4.1)	Nobel® Branemark® RP 3.75 / 4.1mm	94
Avinent® Externa®	Nobel® Branemark® RP 3.75 / 4.1mm	94
Avinent® HI/IC	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
B&B Dental Conexa™	NobelBiocare™ NobelActive™/Conical NP	125
B&W® External Hexagonal (5.0)	Biomet 3i® Osseotite XP 5.0mm WP	
B&W® External Hexagonal Cylindrical (3.3)	Steri-Oss® Hex Lock 3.25mm	
B&W® Tapered Morse MEI (4.8)	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Bioconcept® Tissue Level	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
BioHorizon® Internal 3.5	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
BioHorizon® Internal 4.5	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
BioHorizon® Internal 5.7	Zimmer® Screw-Vent / Tapered Screw-Vent 5.7mm	152
Biogenesis® Iticon	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
BIO-LOK® 4.0mm IC (Silhouette)	Nobel® Branemark® RP 3.75 / 4.1mm	94
BIO-LOK® External Hex	Nobel® Branemark® RP 3.75 / 4.1mm	94
Biomain® Combine Narrow 3.5	Nobel® Replace Tri-Lobe NP 3.5mm	130
Biomain® Combine Regular 4.3	Nobel® Replace Tri-Lobe RP 4.3mm	132
Biomain® Combine Wide 5.0	Nobel® Replace Tri-Lobe WP 5.0mm	134
Biomain® External Hexagon Narrow (3.5)	Nobel® Branemark® NP 3.3mm	
Biomain® External Hexagon Regular (4.1)	Nobel® Branemark® RP 3.75 / 4.1mm	94
Biomain® External Hexagon Wide (5.1)	Nobel® Branemark® WP 5.1mm	
Biomain® I-Fix Regular (4.8)	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Bioner® Hikelt 3.8	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Bioner® Hikelt 4.5	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
Bioner® TopDM	Nobel® Active / Conical Connection NP 3.5mm	125
Bionnovation® Hexagono Ext HE Classic Cilindrico SP (3.3)	Steri-Oss® Hex Lock 3.25mm	
Bionnovation® Hexagono Externo HE Classic WP (5.0)	Intra-Lock® Threaded Conic 4.75/ 5.0mm Ex Hex	
Bionnovation® Hexagono Externo HE RP (4.1)	Nobel® Branemark® RP 3.75 / 4.1mm	94
Biotech Korum® External (3.3) / TOP (3.4) / TRE (3.4)	Biomet 3i® 3.25mm & 3.3mm MicroMiniplant®	

The compatibility information referenced in these charts is based on published, available data and design equivalence. This information is intended to assist dental providers with treatment and diagnosis. Preat Corporation does not provide any diagnosis or suggested medical treatment in the report.

Implant Compatibility Chart_(cont'd)

System	Compatible With	Page
Biotech® Korum External / TOP / TRE (4.0)	Nobel® Branemark® RP 3.75 / 4.1mm	94
Biotech® Korum External / TOP / TRE (5.0)	Biomet 3i® Osseotite XP 5.0mm WP	
Biotech® SPR/SPTT/CIM	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Blue Sky Bio® Conus 12 3.5/4.0	Astra® Tech Osseospeed 3.5/4.0s (AQUA)	79
Blue Sky Bio® Conus 12 4.5/5.0	Astra® Tech Osseospeed 4.5/5.0s (LILAC)	81
Blue Sky Bio® Forte	Nobel® Active / Conical Connection NP 3.5mm	125
Blue Sky Bio® Internal Hex 3.5mm Green	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Blue Sky Bio® Internal Hex 4.5mm Purple	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
Blue Sky Bio® One Stage RN (4.1mm)	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Blue Sky Bio® Three	Astra® Tech Osseospeed 3.0s	78
Blue Sky Bio® Tri-Lobe 6.0mm	Nobel® Replace Tri-Lobe 6.0mm	136
Blue Sky Bio® Tri-Lobe NP	Nobel® Replace Tri-Lobe NP 3.5mm	130
Blue Sky Bio® Tri-Lobe RP	Nobel® Replace Tri-Lobe RP 4.3mm	132
Blue Sky Bio® Tri-Lobe WP	Nobel® Replace Tri-Lobe WP 5.0mm	134
Blue Sky® MAX	Nobel® Active / Conical Connection NP 3.5mm	125
Branemark® Compatible		94
BTI® Externa and Externa Plus RP	Nobel® Branemark® RP 3.75 / 4.1mm	94
Calcitek® Threadlock 4.0mm	Nobel® Branemark® RP 3.75 / 4.1mm	94
CORTEK® Dynamix/Classix/Saturn	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Cowell® Atlas Internal	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
DENTIN® Classic, Prestige®, Rapid®	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Dentium®	Hiossen® Regular 4.0 / 4.5 / 5.0mm	105
Dentis® Cleanlant	Hiossen® Regular 4.0/4.5/5.0mm	105
DIO Implant UFII	Hiossen® Regular 4.0/4.5/5.0mm	105
Ditron Dental MPI™	Zimmer® Screw-Vent/Tapered Screw-Vent 3.3/3.7/4.1mm (3.5mm)	148
DSI® Internal Hex Implant System	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
DSI® Conical Orange NP	Nobel® Active / Conical Connection NP 3.5mm	125
DSI® Conical Orange RP	Nobel® Active / Conical Connection RP 4.3/5.0mm	127
Easy Implant Master® C 4.5/5.0	Hiossen® Regular 4.0 / 4.5 / 5.0mm	105
Easy Implant Master® S 3.75/4.25/4.75	Hiossen® Regular 4.0 / 4.5 / 5.0mm	105
Easy Implant Master® L 3.3/3.75/4.25/4.75	Hiossen® Regular 4.0 / 4.5 / 5.0mm	105
Euroteknika® Aesthetica 4.1	Straumann® (Tissue Level) Standard RN (4.8mmD)	145

All product names, logos, and brands are property of their respective owners. All company, product and service names are for identification purposes only. Use of these third-party names, logos, and brands does not imply, nor should it be understood as, any indication that the respective third-party brand owners have endorsed or sponsored, or are otherwise affiliated with or connected to, Preat or Preat's goods or services.

Implant Compatibility Chart *(cont'd)*

System	Compatible With	Page
Euroteknika® Aesthetica 4.8	Straumann® (Tissue Level) Wide Neck (6.5mmD)	147
Euroteknika Naturall	Astra 3.5-4.0 (Aqua)	79
Glidewell Laboratories® Inclusive 3.5mm	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Glidewell Laboratories® Inclusive 4.5mm	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
GT Medical®	Straumann® (Tissue Level) Standard RN (4.8mmD)	145
Hahn Tapered® 3.0mm	Nobel® Active / Conical Connection 3.0mm	123
Hahn Tapered® 3.5, 4.3mm	Nobel® Active / Conical Connection NP 3.5mm	125
Hahn Tapered® 5.0mm	Nobel® Active / Conical Connection RP 4.3/5.0mm	127
Hahn Tapered® 7.0mm	Nobel® Active / Conical Connection WP 5.0mm	129
Hiossen® GS III Regular 4.0 / 4.5 / 5.0mm	Hiossen® Regular 4.0 / 4.5 / 5.0mm	105
Hiossen® GS Ultra Wide 6.0 / 7.0mm	Hiossen® Regular 4.0 / 4.5 / 5.0mm	105
Hiossen® SS I / II / III Regular 4.8mm	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Hiossen® SS III Wide 6.0mm	Hiossen® SS II Wide 6.0mm	
Hiossen® SS Ultra Wide 6.0mm	Hiossen® SS II Wide 6.0mm	
Hiossen® US II / III Plus Regular 4.1mm	Nobel® Branemark® RP 3.75 / 4.1mm	94
Hiossen® US II / III Plus Wide 5.1mm	Nobel® Branemark® WP 5.1mm	
Hiossen® US II / III Regular 4.1mm	Nobel® Branemark® RP 3.75 / 4.1mm	94
Hiossen® US II / III Wide 5.1mm	Nobel® Branemark® WP 5.1mm	
Hiossen® US II Mini 3.5mm	Nobel® Branemark® NP 3.3mm	
Hiossen® US II Plus Mini 3.5mm	Nobel® Branemark® NP 3.3mm	
Hiossen® US Mini 3.5mm	Nobel® Branemark® NP 3.3mm	
Hiossen® US Ultra Wide 5.1mm	Nobel® Branemark® WP 5.1mm	
Hi-Tec® Implex / Summitt 4.1mm	Nobel® Branemark® RP 3.75 / 4.1mm	94
Hi-Tec® Logic 3.5mm	Nobel® Active / Conical Connection NP 3.5mm	125
Hi-Tec® Logic 4.3 / 5.0mm	Nobel® Active / Conical Connection RP 4.3/5.0mm	127
Hi-Tec® Logic Plus™	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Hi-Tec® Smooth & Slim Fit	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Hi-Tec® Tapered Self Thread 3.5mm	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Hi-Tec® Tapered Self Thread 4.5mm	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
Hi-Tec® Tapered Tite Fit 4.8	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Hi-Tec® Tite Fit™	Straumann® SynOcta (Tissue Level) Standard RN (4.8mm Ø)	145
Hi-Tec® Vision 3.5mm	Nobel® Replace Tri-Lobe NP 3.5mm	130

The compatibility information referenced in these charts is based on published, available data and design equivalence. This information is intended to assist dental providers with treatment and diagnosis. Preat Corporation does not provide any diagnosis or suggested medical treatment in the report.

Implant Compatibility Chart_(cont'd)

System	Compatible With	Page
Hi-Tec® Vision 4.3mm	Nobel® Replace Tri-Lobe RP 4.3mm	132
Hi-Tec® Vision 5.0mm	Nobel® Replace Tri-Lobe WP 5.0mm	134
IBS® Magic FC/N.R Fix	Hiossen® Regular 4.0 / 4.5 / 5.0mm	105
IDO Implant®	Hiossen® Regular 4.0 / 4.5 / 5.0mm	105
IDS Reflect® Aspire 3.5/4.0mm	Astra® Osseospeed 3.5/4.0s (AQUA)	79
IDS Reflect® Aspire 4.5/5.0mm	Astra® Osseospeed 4.5/5.0s (LILAC)	81
IDS Reflect® Rapid 3.5mm	Nobel® Active / Conical Connection NP 3.5mm	125
IDS Reflect® Rapid 4.3mm	Nobel® Active / Conical Connection RP 4.3/5.0mm	127
IDS Reflect® Rapid 5.0mm	Nobel® Active / Conical Connection WP 5.0mm	129
IDS Reflect® Recover 3.5mm	Nobel® Active / Conical Connection NP 3.5mm	125
IDS Reflect® Recover 4.3mm	Nobel® Active / Conical Connection RP 4.3/5.0mm	127
IDS Reflect® Recover 5.0mm	Nobel® Active / Conical Connection WP 5.0mm	129
IDS Reflect® Tapered 3.7-4.1	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
IDS Reflect® Tapered 4.7	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
Implanet®	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Implant Direct® InterActive™/Simply Iconic™ NP	Nobel® Active / Conical Connection NP 3.5mm	125
Implant Direct® InterActive™/Simply Iconic™ RP	Nobel® Active / Conical Connection RP 4.3/5.0mm	127
Implant Direct® Legacy™ 1/2/3/4 3.5mmD	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Implant Direct® Legacy™ 1/2/3/4 4.5mmD	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
Implant Direct® Legacy™ 1/2/3/4 5.7mmD	Zimmer® Screw-Vent / Tapered Screw-Vent 5.7mm	152
Implant Direct® RePlant® / RePlus® / ReActive® NP	Nobel® Replace Tri-Lobe NP 3.5mm	130
Implant Direct® RePlant® / RePlus® / ReActive® RP	Nobel® Replace Tri-Lobe RP 4.3mm	132
Implant Direct® RePlant® / RePlus® / ReActive® WP	Nobel® Replace Tri-Lobe WP 5.0mm	134
Implant Direct® SwishPlant / SwishPlus (4.8mmD)	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Implant Genesis® Aktiv System	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Imtec® 3.40/ 3.75/ 4.0mm external hex	Nobel® Branemark RP 3.75 / 4.1mm	94
IMZ® 3.3 / 4.0mm External Hex	Nobel® Branemark RP 3.75 / 4.1mm	94
Innova® Endopore 4.1 / 5.0mm (One-Stage)	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Innova® Endopore 4.1mm	Nobel® Branemark RP 3.75 / 4.1mm	94

All product names, logos, and brands are property of their respective owners. All company, product and service names are for identification purposes only. Use of these third-party names, logos, and brands does not imply, nor should it be understood as, any indication that the respective third-party brand owners have endorsed or sponsored, or are otherwise affiliated with or connected to, Preat or Preat's goods or services.

Implant Compatibility Chart *(cont'd)*

System	Compatible With	Page
Innova Endopore 5.0mm / Entegra 4.7mm (0.7mm) Ext Hex	Biomet 3i® Osseotite XP 5.0mm WP	
Innova Entegra 3.25 - 4.1mm (0.7mm) External Hex	Nobel® Branemark® RP 3.75 / 4.1mm	94
Innova Hex Plus 4.1mm	Nobel® Branemark® RP 3.75 / 4.1mm	94
Innova Pro Series (TL) Tissue Level 4.1 / 4.8mm	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Innova Pro Series XRT Bone Level 4.8mm	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Intra-Lock 3.75 / 4.0mm External Hex	Nobel® Branemark RP 3.75 / 4.1mm	94
Intra-Lock 4.3mm Standard Platform Conic Body	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
J Dental Care JDEvolution	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Lifecore Stage One 3.3 / 4.1 / 4.8mm (RDS)	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Leader Tixos Internal Hex	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Luna RP 2.5mm	Hiossen® Regular 4.0 / 4.5 / 5.0mm	105
MegaGen Anyone® RP (Internal)	Hiossen® Regular 4.0 / 4.5 / 5.0mm	105
MegaGen ExFeel External Hex	Nobel® Branemark RP 3.75 / 4.1mm	94
MegaGen ExFeel External Hex Small (3.5)	Nobel® Branemark NP 3.3mm	
MegaGen ExFeel External Hex Wide (5.0)	Biomet 3i® Osseotite XP 5.0mm WP	
MegaGen ExFeel Internal Straight/Taper (4.8)	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
MIS C4 3.75 / 4.20mm	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
MIS Lance External hex Standard Platform 4.1mm	Nobel® Branemark® RP 3.75 / 4.1mm	94
MIS Lance External hex Wide Platform (5.0)	Biomet 3i® Osseotite XP 5.0mm WP	
MIS Mistral Internal Morse Taper One Stage RN	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
MIS Seven 3.75 / 4.20mm	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
MIS Seven 4.5	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
Neobiotech IS Implant: System	Hiossen® Regular 4.0 / 4.5 / 5.0mm	105
Neodent Titamax TI Medular (4.1)	Nobel® Branemark RP 3.75 / 4.1mm	94
Neodent Titamax TI Medular (5.0)	Intra-Lock® Threaded Conic 4.75/ 5.0mm Ex Hex	
Noris Medical	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
OsteoReady	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
OSTEO Implant System External Hex	Nobel® Branemark® RP 3.75 / 4.1mm	94
Osteoplus Shi	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148

The compatibility information referenced in these charts is based on published, available data and design equivalence. This information is intended to assist dental providers with treatment and diagnosis. Preat Corporation does not provide any diagnosis or suggested medical treatment in the report.

Implant Compatibility Chart^(cont'd)

System	Compatible With	Page
PALTOP® Advanced Classic	Zimmer Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
PALTOP® Advanced +	Zimmer Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
PALTOP® Dynamic	Zimmer Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
PALTOP® DIVA	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
PALTOP® Conical Active	Hiossen® GS II & GS III Mini 3.5mm	103
Paragon® Bio-Vent 3.5mm Diameter (3.5mm D)	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Paragon® Bio-Vent 4.5mm Diameter (4.5mm D)	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
Paragon® Core-Vent 4.3mm Diameter (3.5mm D)	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Paragon® Core-Vent 5.3mm Diameter (4.5mm D)	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
Paragon® Micro-Vent 3.25mm Diameter (3.5mm D)	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Paragon® Micro-Vent 4.25mm Diameter (4.5mm D)	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
Paragon® Swede-Vent (TL) 4.0mm	Nobel® Branemark® RP 3.75 / 4.1mm	94
Paragon® Swede-Vent 4.0mm	Nobel® Branemark® RP 3.75 / 4.1mm	94
Paragon® Swiss Plus 4.8mm	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Paragon® Taperlock	Nobel® Branemark® RP 3.75 / 4.1mm	94
Park Dental Research® Star Vent External Hex Screw (3.3)	Steri-Oss® Hex Lock 3.25mm	
Park Dental Research® Star Vent External Hex Screw (4.1)	Nobel® Branemark® RP 3.75 / 4.1mm	94
Park Dental Research® Star Vent External Hex Screw (5.0)	Intra-Lock® Threaded Conic 4.75/ 5.0mm Ex Hex	
Park Dental Research® Star Vent External Hex Screw (5.0)	Intra-Lock® Threaded Conic 4.75/ 5.0mm Ex Hex	
ProClinic® Cylindrical/Conical Internal	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
ProClinic® SP Octa 4.8	Straumann® (Tissue Level) Standard RN (4.8mmD)	145
ProClinic® SP Octa 6.5	Straumann® (Tissue Level) Wide Neck (6.5mmD)	147
Renew Biocare® Re-Root SD 3.5mm	Straumann® SynOcta (Tissue Level) Standard NN (3.5mmD)	
Restore® / Sustain® SD 3.3 mm External Hex	Steri-Oss® Hex Lock 3.25mm	
Restore® / Sustain® SD 3.75 / 4.0 / 4.2mm External Hex	Nobel® Branemark® RP 3.75 / 4.1mm	94
Restore® / Sustain® SD 5.0mm External Hex	Intra-Lock® Threaded Conic 4.75/ 5.0mm Ex Hex	
Ritter™ Implants (3.5-6.0mm)	Zimmer® Screw-Vent / Tapered Screw-Vent®	148
SGS Dental® P1/P7	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
SGS Dental P1, P7	Zimmer® Screw-Vent / Tapered Screw-Vent® 3.3 / 3.7 / 4.1mm (3.5mm)	148
Sigma Active Classic	Zimmer® Screw-Vent / Tapered Screw-Vent®	148
Sigma Mega X	Hiossen® Regular 4.0/4.5/5.0mm	105
Simply Implants®	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Southern® 12 degrees 4.3 (3.5mmD)	Nobel® Replace Tri-Lobe NP 3.5mm	130

All product names, logos, and brands are property of their respective owners. All company, product and service names are for identification purposes only. Use of these third-party names, logos, and brands does not imply, nor should it be understood as, any indication that the respective third-party brand owners have endorsed or sponsored, or are otherwise affiliated with or connected to, Preat or Preat's goods or services.

Implant Compatibility Chart *(cont'd)*

System	Compatible With	Page
Southern® 12 degrees 5.0 (4.3mmD)	Nobel® Replace Tri-Lobe RP 4.3mm	132
Southern® 12 degrees 7.0 (5.0mmD)	Nobel® Replace Tri-Lobe WP 5.0mm	134
Southern® CoAxis OCT-XT / MAX-IT 4.0 / 5.0mm (4.8mmD)	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Southern® CoAxis OCT-XT / MAX-IT 6.0mm (6.5mmD)	Straumann® SynOcta (Tissue Level) Standard WN (6.5mmD)	147
Southern® CoAxis TriNeck / TriMax 6.0 / 8.0 / 9.0 (6.0mmD)	Nobel® Replace Tri-Lobe 6.0mm	136
Southern® CoAxis TriNeck / TriMax NP	Nobel® Replace Tri-Lobe NP 3.5mm	130
Southern® CoAxis TriNeck / TriMax RP	Nobel® Replace Tri-Lobe RP 4.3mm	132
Southern® CoAxis TriNeck / TriMax WP	Nobel® Replace Tri-Lobe WP 5.0mm	134
Southern® Parallel / Tapered External Hex NP	Nobel® Branemark NP 3.3mm	
Southern® Parallel / Tapered External Hex WP	Biomet 3i® Osseotite XP 5.0mm WP	
Spiral Tech Premium	NobelBiocare™ NobelReplace™ Tri-Lobe	130
Spiral Tech Ultimate	NobelBiocare™ NobelActive™ Conical	123
Steri-Oss® (no-hex) Flat Top 3.8mm	Calcitek® Integral 4.0mm	
Steri-Oss® Hex Lock 3.8 / 4.5mm	Nobel® Branemark® RP 3.75 / 4.1mm	94
Sterngold® Hex Cylinder 3.3 / 4.0mm	Nobel® Branemark® RP 3.75 / 4.1mm	94
Sterngold® Hex Screw 3.75 / 4.0 / 5.0mm	Nobel® Branemark® RP 3.75 / 4.1mm	94
Sterngold® IC 3.3 / 4.1 / 4.8mm	Straumann® SynOcta (Tissue Level) Standard RN (4.8mmD)	145
Sterngold® PUR NP	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Sterngold® PUR RP	Zimmer® Screw-Vent / Tapered Screw-Vent 4.7mm Diameter (4.5mm Ø)	150
Sterngold® TRU NP	Nobel® Active / Conical Connection NP 3.5mm	125
Sterngold® TRU RP	Nobel® Active / Conical Connection RP 4.3/5.0mm	127
Surgikor® Fixation™ Conical	NobelBiocare™ NobelActive™ Conical RP	127
Surgikor® Immediate™ Hex, Versatile™ Hex	Zimmer® Screw-Vent®/Tapered Screw-Vent 3.3/3.7/4.1mm (3.5mm)	148
Surgikor® Hex Regular and Wide	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Sustain® 4.0 / 4.7mm Internal Bevel (flat top)	Calcitek® Integral 4.0mm	
Sybron® Hex Plus 3.5mm	Nobel® Branemark® NP 3.3mm	
Sybron® Hex Plus 5.1mm	Nobel® Branemark® WP 5.1mm	
Sybron® Pitt Easy 3.8 (3.75mm D) (Yellow)	Sybron® Pitt Easy 3.3 (3.25mm D) (Green)	
Sybron® Pitt Easy 4.1 (4.0mm D) (Red)	Sybron® Pitt Easy 3.3 (3.25mm D) (Green)	
Sybron® Pitt Easy 4.9mm (Blue)	Sybron® Pitt Easy 3.3 (3.25mm D) (Green)	
The Little Implant Company®	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Tre-Oss® Simple	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Tri Dental Implants®	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148
Uni-Dental Direct® IN-SIX	Zimmer® Screw-Vent / Tapered Screw-Vent 3.3 / 3.7 / 4.1mm (3.5mm Ø)	148

The compatibility information referenced in these charts is based on published, available data and design equivalence. This information is intended to assist dental providers with treatment and diagnosis. Preat Corporation does not provide any diagnosis or suggested medical treatment in the report.

Torque Values

System	Prosthetic Screw Connection	Prosthetic Screw Thread	Manufacturer Recommended Torque
3i®			
3i Certain Internal	0.048" (1.22mm) Hex	M1.6 x 0.35	20 Ncm
3i Certain External	Square	M2 x 0.40	32-35 Ncm
Ankylos®			
C/X	1mm Hex	M1.8 x 0.35	15 Ncm
Astra®			
Astra 3.0	0.050" (1.27mm) Hex	M1.4 x 0.30	15 Ncm
Astra 3.5-4.0 (Aqua)	0.050" (1.27mm) Hex	M1.6 x 0.35	20 Ncm
Astra 4.5-5.0 (Lilac)	0.050" (1.27mm) Hex	M2 x 0.40	25 Ncm
Astra EV 3.0, 3.6, 4.2, 4.8, 5.4	0.050" (1.27mm) Hex		25 Ncm
BioHorizons®			
External, Internal, Internal Tapered	0.050" (1.27mm) Hex		30 Ncm
BlueSky Bio®			
One Stage	Star & 0.048" (1.20mm) Hex		30 Ncm
Trilobe	Unigrip & 0.048" (1.20mm) Hex		30 Ncm
Internal Hex & Molar	0.050" (1.27mm) Hex		30 Ncm
Conus 12 & Three	0.050" (1.27mm) Hex		30 Ncm
Max	0.048" (1.20mm) Hex		30 Ncm
Quatro	0.048" (1.20mm) Hex		30 Ncm
Camlog® ScrewLine			
3.3, 3.8, 4.3mm	0.050" (1.27mm) Hex	M1.6 x 0.35	20 Ncm
5.0, 6.0mm	0.050" (1.27mm) Hex	M2 x 0.40	20 Ncm
Endopore			
Endopore (Sybron)	0.050" (1.27mm) Hex		30 Ncm
Dentsply®			
Fridadent/Frialit & Xive	0.048" (1.20mm) Hex		24 Ncm
Hiossen®			
Mini	0.048" (1.20mm) Hex	M1.6 x 0.35	20 Ncm
Regular	0.048" (1.20mm) Hex	M2 x 0.40	30 Ncm
Implant Direct®			
Legacy™ 3.0mm	0.050" (1.27mm) Hex	M1.6 x 0.35	30 Ncm

Torque Values *(cont'd)*

System	Prosthetic Screw Connection	Prosthetic Screw Thread	Manufacturer Recommended Torque
Keystone®			
Keystone Prima Connex	Square	M1.8 x 0.35	30 Ncm
SD Restore	0.048" (1.20mm) Hex, Square	Square 1-72 UNF	30 Ncm
RD Restore	0.048" (1.20mm) Hex, Square	M2 x 0.40	30 Ncm
WD Restore 0.048" (1.22mm) Hex, Square	0.048" (1.20mm) Hex, Square	M2.5 x 0.45	30 Ncm
Keystone Renova	0.048" (1.20mm) Hex	M1.8 x 0.35	30 Ncm
MegaGen®			
Anyridge®	0.048" (1.20mm) Hex	M1.8 x 0.35	25 Ncm
Neodent®			
Neodent® CM	0.048" (1.20mm) Hex	M1.8 x 0.35	32 Ncm
Neodent® GM	Neo Screwdriver		20 Ncm
Neoss®			
Neoss®	Neoss Proprietary	M2 x 0.40	32 Ncm
Nobel Biocare®			
NobelActive™ 3.0	Unigrip	M1.4 x 0.30	15 Ncm
NobelActive™ NP	Unigrip	M1.6 x 0.35	35 Ncm
NobelActive™ RP/WP	Unigrip	M2 x 0.40	35 Ncm
Branemark NP	Unigrip	M1.6 x 0.35	35 Ncm
Branemark RP	Unigrip	M2 x 0.40	35 Ncm
Branemark WP	Unigrip	M2.5 x 0.45	35 Ncm
NobelReplace™ NP	Unigrip	M1.8 x 0.35	35 Ncm
NobelReplace™ RP	Unigrip	M2 x 0.40	35 Ncm
Straumann® Tissue-level			
Narrow Neck SCS	SCS	M1.8 x 0.35	35 Ncm
Regular Neck SCS	SCS	M2 x 0.40	35 Ncm
Wide Neck	SCS	M2 x 0.40	35 Ncm
Straumann® Bone-level			
Narrow Cross-fit SCS	SCS	M1.6 x 0.35	35 Ncm
Regular Cross-fit SCS	SCS	M1.6 x 0.35	35 Ncm
Straumann® BLX			
RB/WB	SCS	M1.6 x 0.66	35 Ncm
Zimmer®			
Zimmer® ScrewVent	0.050" (1.27mm) Hex	1-72 UNF	30 Ncm

Preat

PATIENT-SPECIFIC ABUTMENTS

Elevating Your Restorations TO THE HIGHEST DEGREE

Preat is setting the standard with advanced Patient-Specific Abutment solutions. With our innovative angled screw channel abutment, offering up to 30° of angulation, we're not just supporting your workflow – we're empowering your ability to provide the best restorations possible!

From fully validated in-house milling to expert overflow and outsourcing support, you have the freedom to own your own strategy, meeting your needs today and accelerating your growth tomorrow.



COMPATIBILITY

All Major Implant Systems



QUALITY

Fit Meets or Exceeds OEM





SPEED

Two-Day Turnaround Time



PROTECTION

All OEM Implant Warranties Matched



LEARN MORE

visit us at www.preat.com/preat-patient-specific-abutments
to discover Preat Patient-Specific Abutments

Your Honorary Team Member

Let us be your implant restoration expert. We've spent decades developing knowledge of the many restorative solutions available on the market and instilling that expertise within our team. We believe in helping patients access the best care possible and approach each question with that in mind.

Connect on Your Terms



phone



email



mobile



live chat

**Within The USA
Customer Service**

6am - 5pm PST

t: 1.800.232.7732 | f: 805.202.3076 | www.preat.com

**Outside The USA
Customer Service**

6am - 5pm PST

t: +1.805.202.3070 | f: +1.805.202.3076

Technical Support

technical@preat.com

Customer Support

info@preat.com

Orders

orders@preat.com



PREAT.COM • 800.232.7732

All product names, logos, and brands are property of their respective owners. All company, product and service names are for identification purposes only. Use of these third-party names, logos, and brands does not imply, nor should it be understood as, any indication that the respective third-party brand owners have endorsed or sponsored, or are otherwise affiliated with or connected to, Preat or Preat's goods or services.

©2025 Copyright Preat Corporation • All rights reserved • Specifications are subject to change without notice.

Rev. 1/26